



-  Surgery & Operation Products
-  Intensive Care Products
-  New Born Products
-  Urology Products
-  Wound Care Support Products
-  Patient Safety Products
-  General Products




**General Catalogue**



Medikokim established 1988 in Istanbul and started production process concerning its first product in 1990. In this period developed its product range and specifications.

Today, with its modern technology and experienced staff Medikokim is manufacturing disposable products of its own according to the international standards, all of them are first manufactured in Turkey, becomes a contributor to the economy and the health sector of the country.

Products of the Medikokim are used with success in the region of Europe, Asia, Africa and Middle America in 28 countries, in addition to Turkey.




# Index

<b>Surgery &amp; Operation Products</b>	06	<b>SAUGERMED</b> Surgical Brain Pads
	10	<b>SAUGERMED</b> Surgical Brain Strips
	12	<b>SAUGERMED</b> Surgical Nasal Pads
	13	<b>CAPAD</b> Drainage Cannula Pads
<b>Intensive Care Products</b>	16	<b>NASOMED</b> Nasal Cannula Fixation Bands
	17	<b>TRAPAD</b> Tracheal Cannula Pads
	18	<b>TRAEGERMED</b> Tracheal Cannula Holders
	20	<b>STAENDERMED</b> Endotracheal Cannula Holders
	22	<b>DECKELMED</b> Ear Support Band
	23	<b>STOPPERMED</b> Nasal Tampon Holder
	24	<b>PORTMED</b> Multiple Catheter Holder
	24	<b>PULSOMED</b> Pulse Oximeter Probe Wrap
	25	<b>ASMED</b> Angio Support Band
	25	<b>SCHÜTZERMED</b> Eye Protective Band
	26	<b>FIXELMED</b> Urine Bag Holder
	27	<b>FIXERMED</b> Male Urine Catheter Holder
	28	<b>MITTMED</b> Finger Control Mitts
	29	<b>RESTMED</b> Freedom Splint
	30	<b>FESTERMED</b> Hand & Foot Fixation Bands
	32	<b>Fall Management Socks</b>
	33	<b>Fall Precaution Bracelets</b>
	33	<b>Essential Fall Management Kit</b>
<b>New Born Products</b>	36	<b>BINDERMED</b> Baby Limb Holder
	36	<b>NASOMED</b> Baby Nasal Cannula Fixation Band
	37	<b>FIXERBOARD</b> Baby Limb Boards
	37	<b>MAMASAFE</b> Breast Protection Pad
	38	<b>SCHÜTZERMED</b> Phototherapy Eye Bands
<b>Urology Products</b>	42	<b>CLAMPMED</b> Urine Control Clamp
	43	<b>BAREMED</b> Urine Preventer Band
	43	<b>CONDOFIX</b> Condom Fixer
<b>Wound Care Support Products</b>	46	<b>FERSEMED</b> Heel Protector
	47	<b>FERSEMED Rg</b> Foot Elevator
	48	<b>WUNPAD</b> Wound Care and Medical Dressing Pad
<b>Patient Safety Products</b>	51	<b>SAFEBELT</b> Bed & Stretchers Safety Belts
	52	<b>SAFEBELT</b> Wheelchair Safety Belts
	54	<b>SAFEBELT</b> Patient Transfer Safety Belts
	54	<b>SAFEBELT</b> Optional Belts
	55	<b>SAFEBELT</b> Accessories
	57	<b>SAFEBELT</b> Wrist - Ankle Fixation Belts
	58	<b>SAFEBELT</b> Wheelchair Safety Belts
	59	<b>SAFEBELT</b> Patient Support Belts
<b>General Product</b>	62	<b>UNI-GLOVE</b> Operating Room Light Handle Cover
	63	<b>PASSPOP</b> Patient Identification Band



## Surgery & Operation Products



06

**SAUGERMED**

Surgical Brain Pads



10

**SAUGERMED**

Surgical Brain Strips



12

**SAUGERMED**

Surgical Nasal Pads



13


**CAPAD**

Drainage Cannula Pads

# SAUGERMED

## Surgical Brain Pads

For neurosurgical operations  
special produced haemostatic  
featured pads.



### Enables safe use

Also, all the pads have an atraumatic featured Trace Rope with smooth touch. Pad locations can easily be traced and identified by these trace ropes.  
Pls. see figure 1

To avoid any distention Trace rope looped with a special technics\*. Pls. see figure 2

All the pads consists of radio-pack blue strap which is X-ray detectable. Pls. see figure 3

To recognize the number of used pads during the operation it has a Special Counter Card 1 to 10. Pls. see figure 4

To prevent the pads from deformation it consists of a protective belt.  
Pls. see figure 5

### Atraumatic

Produced from hydrofiber structured organic and inorganic materials in nonwoven form Atraumatic and soft.

Unfuzzy and unfeathered, does not leave loofah.

### Sterile and dirtiness protected

Gamma sterilized.

Tyek made of sterile pouch prevents unbounded tear and can only be open from its opening part.

It has an anticontamination bag in the sterile pouch.

There are 16 different dimensions of each pad made of 11 various raw material.

\* Patent rights of the special loop protected by the certification of the Turkish Patent Institute as Utility Model.



Figure 1



Figure 2



Figure 4





Figure 3



Figure 5



# SAUGERMED

Surgical Brain Pads

Various options for different use of purpose 2 main featured and structured pad produced from 11 various raw material.

## Viscose Pads



### SAUGERMED Thinpad

- Absorption capacity 8 times of its weight.
- Very thin, smooth and silky structured.
- Made of viscose and synthetic additive material.




### SAUGERMED Plus

- Absorption capacity 11 times of its weight .
- Thin, closely structured.
- Made of viscose and synthetic additive material.



### SAUGERMED Thincot

- Absorption capacity 9 times of its weight.
- Thin, smooth structured, laminated surface.
- Made of viscose and synthetic additive material.




### SAUGERMED Laminat

- Absorption capacity 10 times of its weight.
- Standard thickness, smooth structured, laminated surface.
- Made of viscose and synthetic additive material.



### SAUGERMED Soft

- Absorption capacity 10 times of its weight
- Standard thickness, closely structured, laminated surface.
- Made of viscose and synthetic additive material.




### SAUGERMED Ultra Soft

- Absorption capacity 15 times of its weight.
- Smooth , delicate structured.
- Made of viscose and synthetic additive material.



### SAUGERMED Cotton Soft

- Absorption capacity 16 times of its weight.
- Very smooth, delicate structured.
- Made of viscose and synthetic additive material.




## Cotton Pads




### SAUGERMED Pure Cotton

- Absorption capacity 15 times of its weight.
- Very thin, very light structured.
- Made of 100% cotton raw materials




### SAUGERMED Pure Cotton Thick

- Absorption capacity 16 times of its weight
- Smooth, delicate structured.
- Made of 100% cotton raw materials



### SAUGERMED Pure Cotton Medium

- Absorption capacity 16 times of its weight
- Smooth, delicate structured.
- Made of 100% cotton raw materials



### SAUGERMED Pure Cotton Soft


- Absorption capacity 17 times of its weight
- Very soft, very absorbative, closely structured.
- Made of 100% cotton raw materials.

\*Please contact with our sales department for order informations of Saugarmed.


# SAUGERMED

## Surgical Brain Strips

7 different strips for different purposes.



- Produced by nonwoven form with materials of hydrofiber speciality.
- Atraumatic and soft.
- Unfeathered, unfuzzy, does not leave loofah.
- They can haemostatize smallest compressions.
- All strips contains x-ray detectable Radio-pack speciality blue strap
- Strips are attached to the counter card, which consists a protective. band made of carton to prevent from damage.



1:1 scale product table (mm.)

\*Please contact with our sales department for order informations of Saugermesdv strips.

CE 1783




### **SAUGERMED Thinpad**

- Absorption capacity 8 times of its weight.
- Very thin, smooth and silky structured.
- Made of viscose and synthetic additive material.



### **SAUGERMED Plus**

- Absorption capacity 11 times of its weight .
- Thin, closely structured.
- Made of viscose and synthetic additive material.



### **SAUGERMED Thincot**

- Absorption capacity 9 times of its weight.
- Thin, smooth structured, laminated surface.
- Made of viscose and synthetic additive material.




### **SAUGERMED Laminat**

- Absorption capacity 10 times of its weight.
- Standard thickness, smooth structured, laminated surface.
- Made of viscose and synthetic additive material.




### **SAUGERMED Soft**

- Absorption capacity 10 times of its weight
- Standard thickness, closely structured, laminated surface.
- Made of viscose and synthetic additive material.



### **SAUGERMED Ultra Soft**

- Absorption capacity 15 times of its weight.
- Smooth , delicate structured.
- Made of viscose and synthetic additive material.





### **SAUGERMED Cotton Soft**

- Absorption capacity 16 times of its weight.
- Very smooth, delicate structured.
- Made of viscose and synthetic additive material.

# SAUGERMED

## Surgical Nasal Pads

Special produced hemostatic featured surgical pad for the use of nose and ear surgery.



### Enables safe use

It consists of a trace rope which is x-ray detectable and does not cause any atraumatic effect by touching the tissue. Trace rope also x-ray detectable. With a special method knotted trace rope, to avoid any fabric pot, can not detached from the pad.

By its trace rope, easy to identify the place of the pad and to take off to recognize the number of the used pads during the operation it has a Special Counter Card 1 to 10.

### Atraumatic

Both sides micro laminated to be smooth.  
Absorption capacity 8 times of its weight.

### Advantage

Both sides micro laminated to be smooth.  
Absorption capacity 8 times of its weight.

### Sterile and dirtiness protected

Gamma sterilized Tyvek made of sterile pouch provides unbounded tear and can only be open from its opening part

It has an anticontamination bag in sterilized pouch.


### Saugermed Ordering Information

Model	Size (mm)	Piece / Pouch	Pouch / Box	Box Size (cm)	Box Weight	Order Code
Saugermed	10 x 20	10	24	15 x 20 x 30	330 gr.	NP 1020
	10 x 30	10	24	15 x 20 x 30	330 gr.	NP 1030
	10 x 40	10	24	15 x 20 x 30	333 gr.	NP 1040
	10 x 60	10	24	15 x 20 x 30	338 gr.	NP 1060
	10 x 75	10	24	15 x 20 x 30	352 gr.	NP 1075
	12 x 20	10	24	15 x 20 x 30	312 gr.	NP 1220
	12 x 30	10	24	15 x 20 x 30	319 gr.	NP 1230
	12 x 40	10	24	15 x 20 x 30	340 gr.	NP 1240
	12 x 50	10	24	15 x 20 x 30	347 gr.	NP 1250
	12 x 60	10	24	15 x 20 x 30	350 gr.	NP 1260
	12 x 75	10	24	15 x 20 x 30	357 gr.	NP 1275
	20 x 30	10	24	15 x 20 x 30	345 gr.	NP 2030
	20 x 40	10	24	15 x 20 x 30	362 gr.	NP 2040
	20 x 50	10	24	15 x 20 x 30	376 gr.	NP 2050
	20 x 60	10	24	15 x 20 x 30	369 gr.	NP 2060
	20 x 75	10	24	15 x 20 x 30	393 gr.	NP 2075



# CAPAD

## Drainage Cannula Pads

Designed to use together with the drainage cannula.



Each 3 raw material has separately 2 different model big and small.




### Capad Standart

Both sides absorbent pad.



### Capad Ultra

Double side pad, absorbent on both sides.



### Capad Blocker

One side absorbent, the other side nonpermeable pad.

### Practical

Can fit the every size and model of cannula by its elastic cannula holes.

### Offers alternative options

Different optional pads, each of them produced from various raw material.

### Atraumatic

Produced from smooth, hydrofiber structured nonwoven or hydrophil foam material.

### Sterile

Gamma sterilized.

\*Please contact with our sales department for order informations of Capad.

CE 1783



16

## **NASOMED**

Nasal Cannula Fixation Bands



17

## **TRAPAD**

Tracheal Cannula Pads



22

## **DECKELMED**

Ear Support Band



23

## **STOPPERMED**

Nasal Tampon Holder



25

## **ASMED**

Angio Support Band



25

## **SCHÜTZERMED**

Eye Protective Band



28

## **MITTMED**

Finger Control Mitts




29

## **RESTMED**

Freedom Splint

## Intensive Care Products



18

### TRAEGERMED

Tracheal Cannula Holders



20

### STAENDERMED

Endotracheal Cannula Holders



24

### PORTMED

Multiple Catheter Holder



24

### PULSOMED

Pulse Oximeter Probe Wrap



26

### FIXELMED

Urine Bag Holder



27

### FIXERMED

Hand & Foot Fixation Bands



30

### FESTERMED

Hand & Foot Fixation Bands



32

### Fall Management Socks



33

### Fall Precaution Bracelets



33

### Essential Fall Management Kit

# NASOMED


Nasal Cannula Fixation Bands

3 different model for adults,  
child and babies



## NASOMED Adult

- Covers the head smoothly.
- Safe and easy to use with its velcro bands.
- Skin touch cotton surface absorbs the sweat.
- In case of long term use does not cause allergy and irritation.
- Gamma sterilized.



## NASOMED Child

- Consists of two parts supported both right and left ears.
- Safe and easy to use with its velcro bands.
- Can be used together with every size and model of cannula.
- Skin touch cotton surface absorbs the sweat.
- In case of long term use does not cause allergy and irritation.
- Gamma sterilized.



## NASOMED Baby

- Practical and easy to use with its velcro bands.
- Keeps stable the cannula at the nose of the baby getting through the hoops of the nasal cannula.
- Can be used together with every size and model of cannula.
- Covers the head smoothly, skin touch cotton surface absorbs the sweat.
- In case of long term use does not cause allergy and irritation.
- Gamma sterilized.

### Nasomed Ordering Information

Model	Size (cm)	Piece / Pouch	Pouch / Box	Box Size (cm)	Box Weight	Order Code
Adult	65	1	25	11 x 20 x 25	355 gr.	NTH 5002
Child	32	1	25	11 x 20 x 25	220 gr.	NTH 5001 TP
Baby	43	1	25	11 x 20 x 25	245 gr.	NTH 5001

CE 1783

# TRAPAD

## Tracheal Cannula Pads



3 different raw materials,  
4 models from each raw material,  
Different sizes for each model for  
adults and children.

### Practical

By its elastic cannula holes can be used with every  
type and model of cannula.

### Atraumatic

Produced from smooth, hydrofiber structured of  
nonwoven or hydrophil foam material.

### Offers alternative options

4 different models, each of them produced from  
various structured raw material. Each model has  
various dimensions for adult and children.


### Sterile

Gamma sterilized.



#### **Trapad Standart**

Both sides absorbent.




#### **Trapad Ultra**

Double side pad,  
absorbent on both sides.




#### **Trapad Blocker**


One side absorbent  
the other side  
nonpermeable pad.




Straight slit Cannula  
hole at the top.



Without slit Cannula  
hole at the top.



Zigzag slit Cannula  
centered



Without slit  
Cannula centered.

\*Please contact with our sales department for order informations of Trapad.

CE 1783

# TRAEGERMED

## Tracheal Cannula Holders

Two different structured of holders: Sponge supported and rigid holders.

- Skin colour does not take attention.
- Skin touch cotton surface absorbs the sweat.
- In case of long term use does not cause allergy and irritation.
- Gamma sterilized.



### Sponge Supported (Flexible) Fixation Holders

The sponge layer in the middle provides flexibility and softness to patient 2 different models.



### Single piece model

Suitable for use in self.




### Two pieces model

Better to use for the inpatient. Can be intervene to the cannula without moving.




## Velor Supported Rigid Fixation Holders

Non flexible and fixes canula in the first position thin structured fixation tapes 2 different models.



### Single piece model

Proper for the use of the patients except inpatient.



### Two pieces model

Better to use for the inpatient. Can be intervene to the cannula without moving.

Traegermed Sponge Supported Fixation Holders Ordering Information						
Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
Single Piece Baby	23	1	25	11 x 20 x 25	200 gr.	TCHS 01
Single Piece Child	44	1	25	11 x 20 x 25	230 gr.	TCHS 02
Single Piece Adult	56	1	25	11 x 20 x 25	260 gr.	TCHS 03
Two Pieces	36 + 6	1	25	11 x 20 x 25	290 gr.	TCHS TP


Traegermed Velor Supported Fixation Holders Ordering Information						
Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
Single Piece Baby	23	1	25	11 x 20 x 25	200 gr.	TCHR 01
Single Piece Child	44	1	25	11 x 20 x 25	230 gr.	TCHR 02
Single Piece Adult	56	1	25	11 x 20 x 25	260 gr.	TCHR 03
Two Pieces	36 + 6	1	25	11 x 20 x 25	290 gr.	TCHR TP

# STAENDERMED


## Endotracheal Cannula Holders



CLP



CLP TP



### STAENDERMED - CLP

Clemm supported holder bands


Holder band fixes the cannula by gripping it with its locked clemm (CLP). There are 2 different models.

**CLP - Single piece, velcro connection is on the back of the head model.**

It grips and fixes the endotracheal tube with its clemm, it wraps around the neck and velcro s connects at the back of the head. There are 2 different sizes.

**CLP TP - Two pieces, velcro connection is on the side of the head model.**

It grips and fixes the endotracheal tube with its clemm, it wraps around the neck and velcro s connects at the side of the head. There is only one standard single size.




V1



V2



V3 TP



### STAENDERMED - V

Velcro Supported Holders

Holds the canula with the agency of velcro (V) lock, 3 different models

**V1 Single piece mid-supported**

The middle portion of the tube wraps, 2 sizes.

**V2 Single piece side-supported**


The part close to the tube-to-end wraps 3 sizes.

**V3 Two pieces model**

One size.

- Cotton surfaces touching the skin absorbs the sweat.
- Does not cause allergy and irritation in case of long term use.
- Cannula can be placed and removed without moving the patient's neck

- by the usage of single piece side-supported and two pieces models.
- Three different types of Endotracheal Cannula Holders: Velcro supported, silicon supported and clomp supported.
- Gamma sterilized.



## STAENDERMED - SC

Silicone Supported Holder Bands

Holder band fixes the cannula by gripping it with its silicone surface (S). There are 3 different models.

### **S1** Single piece, supported at middle part model.

The middle part wraps and grips the cannula. There are 2 sizes.

### **S2** Single piece, supported by side part model.

The part close to one end wraps and grips the cannula. There are 3 sizes.

### **S3 TP** Two pieces model

There is only one standard single size.

Staendermed Ordering Information						
Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
V1	55	1	25	11 x 20 x 25	330 gr.	ETH 101
	65	1	25	11 x 20 x 25	345 gr.	ETH 102
V2	45	1	25	11 x 20 x 25	330 gr.	ETH 201
	55	1	25	11 x 20 x 25	345 gr.	ETH 202
	65	1	25	11 x 20 x 25	330 gr.	ETH 203
V3 TP	55	1	25	11 x 20 x 25	330 gr.	ETH TP
S1	55	1	25	11 x 20 x 25	330 gr.	ETH-SC 101
	65	1	25	11 x 20 x 25	345 gr.	ETH-SC 102
S2	55	1	25	11 x 20 x 25	330 gr.	ETH-SC 201
	65	1	25	11 x 20 x 25	345 gr.	ETH-SC 202
	55	1	25	11 x 20 x 25	330 gr.	ETH-SC 203
S3 TP	50 + 28	1	25	11 x 20 x 25	330 gr.	ETH-SC TP
CLP	55	1	25	11 x 20 x 25	330 gr.	ETH CLP 101
	65	1	25	11 x 20 x 25	345 gr.	ETH CLP 102
CLP TP	50 + 25	1	25	11 x 20 x 25	330 gr.	ETH CLP TP

CE 1783

# DECKELMED

## Ear Support Band



### Practical

- D shape protects ear against hits with its high density foam structure.
- Design allows air and voice permeable.
- Two different dimensions for adult and children.
- Covers the head smoothly.
- Easy to use and safe with its velcro bands.
- Skin touch surface absorbs the sweat.

### Offers alternative options

- In case of long term use does not cause allergy and irritation.

### Sterilized

- Gamma sterilized.
- Single use.,

#### Deckelmed Ordering Information

Model	Size (cm)	Piece / Pouch	Pouch / Box	Box Size (cm)	Box Weight	Order Code
Deckelmed	64 cm.	1	25	20 x 30 x 20 cm.	825 gr.	EPH 7001
	54 cm.	1	25	20 x 30 x 20 cm.	775 gr.	EPH 7002

# STOPPERMED

## Nasal Tampon Holder



### Safety

- Supports safely and stabilizes the tampon inside the nose.
- Transparent silicon strip in the middle allows to see possible bleeding and supports tampon, which is not effected when it gets wet.



- It absorbs sweat on cotton surface texture that contacts with the skin.
- It does not cause allergy or irritation on the skin even if it stays longer times.

### Practical

- There are elastic loops that provides the fixation of the plaster below of the patient's nose by passing through both of the ears.
- Velcro locking pieces allow adjusting the length and fixing after placing then around ears, so fits all sizes of heads.
- It is easy and safe to use thanks to velcro tape.

### Sterilized

- Gamma sterilized.
- Single use

Stoppermed Ordering Information						
Model	Size (cm)	Piece / Pouch	Pouch / Box	Box Size (cm)	Box Weight	Order Code
Stoppermed	2.5 x 32	1	25	11 x 20 x 25	240 gr.	NTH 01

CE 1783

# PORTMED

## Multiple Catheter Holder



- By its sewed catheter support band enables to fix the 3 different catheter at the same time.
- Easy to use, pratique and safe.
- Latex free.
- Can be fit to every arm size, does not hinder blood circulation.
- With the catheter pad dimensions of 15 x 8 cm, can be fix 3 different catheter at the same time.
- To distinguish the catheters easily it has 3 different colour of velcro binding.
- Does not cause allergy and irritation.
- Washable at 30 °C, reusable.

### Portmed Ordering Information

Model	Size (cm)	Piece / Pouch	Pouch / Box	Box Size (cm)	Box Weight	Order Code
Portmed	45 x 8	1	25	20 x 30 x 20	820 gr.	HP MED



# PULSOMED

## Pulse Oximeter Probe Wrap



- Covers the probe safely and fix stable.
- Very simple to use.
- Can be used with every model and marked probe.
- Not sterile.

### Pulsomed Ordering Information

Model	Size (cm)	Piece / Pouch	Pouch / Box	Box Size (cm)	Box Weight	Order Code
Pulsomed	3 x 12	1	50	11 x 20 x 25	275 gr.	PPW 9001



# ASMED

## Angio Support Band

Supports the pelvic area after angio operations and prevents swelling.



- Supports the pelvic area after angio operations and prevents swelling.
- Elastic, does not hinder the blood circulation
- Does not cause allergy and irritation in case of long term use.
- Special texture enables the skin breathe.
- Easy to use and safe with its velcro bands.
- Smooth structured, skin touch surface absorbs the sweat.
- Two different type for left and right pubis
- Single size.
- Does not cause allergy and irritation.
- Washable at 30 °C, reusable.

### Asmed Ordering Information

Model	Size (cm)	Piece / Pouch	Pouch / Box	Box Size (cm)	Box Weight	Order Code
Asmed Right	110/60x14	1	12	20 x 30 x 30	2330 gr.	ASBR
Asmed Left	110/60x14	1	12	20 x 30 x 30	2330 gr.	ASBL



# SCHÜTZERMED

## Eye Protective Band



- Prevents from the UV rays up to 450 nanometer.
- Easy to use and practical, can be adjust to every head size.
- Fits the face perfectly and covers softly.
- Maximum safe with its connection elements.
- Gamma sterilized.

### Schützermmed Ordering Information

Model	Size (cm)	Piece / Pouch	Pouch / Box	Box Size (cm)	Box Weight	Order Code
Schützermmed	Standart	1	25	11 x 20 x 25	450 gr.	PEB 304



# FIXELMED

## Urine Bag Holder

Safely fix the urine bag holder to the leg and bears.



### Practical

- A pair of fixation bands fits every size and model of leg urine bag.
- Safe and easy to use, can be washed at 30°C reusable.

### Safe to use

- Elastic structure fits all leg size, covers the leg safely and keep stable.
- Wide part of the skin touch side which consists of silicone strap prevents the slide from the leg.

### Comfort

- Special texture enables the skin breathe, reduces the sweat.
- Skin touch smooth feathery surface absorbs the sweat.
- Does not cause allergy and irritation.

#### Fixelmed Ordering Information

Model	Size (cm)	Piece / Pouch	Pouch / Box	Box Size (cm)	Box Weight	Order Code
Fixelmed	60 x 2 / 60 x 2	1	10	11 x 20 x 25	320 gr.	UBH 1001



# FIXERMED

## Male Urine Catheter Holder

Stables the urine catheter and prevents traumatic effects.



### Practical

- Easy to use, pratique and safe.
- Fits all urine cannula types and models.
- Machine washable at 30 °C reusable.

### Safe to use

- Latex free.
- Elastic structure fits all leg size, covers the leg safely and keep stable.

### Comfort

- Special texture enables the skin breathe, reduces the sweat.
- Skin touch smooth feathery surface absorbs the sweat.
- Does not cause allergy and irritation.

### Fixermed Ordering Information

Model	Size (cm)	Piece / Pouch	Pouch / Box	Box Size (cm)	Box Weight	Order Code
Fixermed	60 x 5	1	25	11 x 20 x 25	680 gr.	CH 4001



# MITTMED

## Finger Control Mitts

It is specifically designed for patients who try to remove catheters, attached cannulas or tend to hurt themselves.



### MITTMED MITTFC 101 Fishnet

- It avoids the trauma with well rounded and soft wad part that placed hand palm.
- The upper part that consisting complete synthetic mesh restrict the movement of the fingers.
- It allows inserted use of IV cannula when get attached on the hand because of its netting tissue.
- Easy and safe with its velcro bands.
- In case of long term using does not cause allergy and irritation.
- One model for both right and left arms.
- Washable at 30 °C.



### MITTMED MITTOTW Open

- It avoids the trauma with well rounded and soft wad part that placed hand palm.
- If necessary , can be interfered to the hand and fingers by opening the upper portion.
- Easy and safe with its velcro bands.
- In case of long term using does not cause allergy and irritation.
- One model for both right and left arms.
- Washable at 30 °C.



### MITTMED MITTCTW Closed

- It avoids the trauma with well rounded and soft wad part that placed hand palm.
- Easy and safe with its velcro bands.
- In case of long term using does not cause allergy and irritation.
- One model for both right and left arms.
- Washable at 30 °C.



### MITTMED MITTSRVTW Therapy

- It helps to decrease patient stress as grabbing balls in the closed part of the mitt.
- Easy and safe with its velcro bands.
- In case of long term using does not cause allergy and irritation.
- One model for both right and left arms.
- Washable at 30 °C.

#### Mittmed Ordering Information

Model	Size (cm)	Piece / Pouch	Pouch / Box	Box Size (cm)	Box Weight	Order Code
Fishnet	26 x 18	1	12	42 x 30 x 30	1100 gr.	MITTFC101
Open	28 x 21	1	12	42 x 30 x 30	1600 gr.	MITTOTW
Closed	29 x 21	1	12	42 x 30 x 30	1950 gr.	MITTCTW
Therapy	29 x 23	1	15	42 x 30 x 30	850 gr.	MITTSRVTW

# RESTMED

## Freedom Splint

It is specifically designed for patients who try to remove catheters, attached cannulas or tend to hurt themselves. By preventing of the patient's elbow movement prevents hands touching through the head and body.



### Standart Series

- 100% cotton lining; waterproof outer covering.
- In case of long term using does not cause allergy and irritation.
- It allows to using in place during X - Rays.

- Easy and safe with its velcro bands and plastic clips.
- Clipped to the patient's sleeve to help prevent sliding.
- One model for both right and left arms.
- Washable at 30 °C.



### Soft Series

- 100% cotton lining and soft, permeable to air, vinly mesh outer covering.
- In case of long term using does not cause allergy and irritation.

- It allows to using in place during X - Rays.
- Easy and safe with its velcro bands and plastic clips.
- Clipped to the patient's sleeve to help prevent sliding.
- One model for both right and left arms.
- Washable at 30 °C.

#### Restmed Ordering Information

Model	Size (cm)	Piece / Pouch	Pouch / Box	Box Dimension (cm)	Box Weight	Order Code
Standart XS	23 x 34	1	10	25 x 50 x35	1000 gr.	RMRB 01
Standart S	28 x 37	1	10	25 x 50 x35	1500 gr.	RMRB 02
Standart M	30 x 47	1	10	25 x 50 x35	2250 gr.	RMRB 03
Standart L	36 x 64	1	10	25 x 50 x35	4000 gr.	RMRB 04
Standart XL	49 x 69	1	10	25 x 50 x35	5000 gr.	RMRB 05
Soft XS	23 x 34	1	10	25 x 50 x35	1000 gr.	RMSW 01
Soft S	28 x 37	1	10	25 x 50 x35	1500 gr.	RMSW 02
Soft M	30 x 47	1	10	25 x 50 x35	2250 gr.	RMSW 03
Soft L	36 x 64	1	10	25 x 50 x35	4000 gr.	RMSW 04
Soft XL	49 x 69	1	10	25 x 50 x35	5000 gr.	RSFRSW 05



# FESTERMED

## Hand & Foot Fixation Bands

Bands to fix the wrists of the patient in order to prevent the movement  
5 different models for various purposes.



### FESTERMED Metal

- For use of patients in bed or litter.
- With its perdurable connection ribbon through the metal hoop can be fixed to the needed place.
- Supported with its all round overlog seam extremely durable.



- Skin touch cotton surface absorbs the sweat.
- In case of long term use does not cause allergy and irritation.
- Can be used after washed 30 °C and making sterilization at autoclave.

CE



### FESTERMED QRB

- Patients assessed to be at risk of disrupting life-saving treatment or whose picking, pulling, scratching or peeling exacerbates a skin condition, causes self-injury, or compromises wound site integrity.
- Soft, and supported material.



- Initial hook and loop closure for additional security.
- Secures to the bed with quick-release buckles or ties.
- Adjustable size fits pediatric or adult limb.
- Wrist circumference 27 cm.

CE



### FESTERMED Deluxe

- Especially for the patient who applies high strength.
- Skin touch cotton surface absorbs sweat.



- In case of long term use does not cause allergy and irritation.
- Washable at 30°C reusable.

CE



## FESTERMED Twain

- For the use of the patients carried by ambulance and first aid.
- Brings the two wrist together and enables the hands at the front of the body .
- Skin touch cotton surface absorbs the sweat.



- In case of need with the velcro lock at the bottom surface can be stable by fixing to the bad or litter.
- In case of long term use does not cause allergy and irritation.

CE



## FESTERMED Classic

- For the use of inpatient both bed or litter.
- Skin touch cotton surface absorbs the sweat.



- In case of long term use does not cause allergy and irritation.
- Gamma sterilized.

CE<sub>1783</sub>

### Festermmed Ordering Information

Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
Deluxe	34	1	25	20 x 30 x 20	1215 gr.	ALH 2006
Deluxe / Çift	34	2	12	20 x 30 x 20	1175 gr.	ALH 2007
Twain	32	1	25	20 x 30 x 20	1275 gr.	TC 2009
QRB	27	1	25	20 x 30 x 20	1790 gr.	QRB 2013
QRB / Çift	27	2	12	20 x 30 x 20	1726 gr.	QRB 2014
Klasik El	25	1	25	20 x 30 x 20	850 gr.	ALH 2001
Klasik El / Çift	25	2	12	20 x 30 x 20	815 gr.	ALH 2008
Klasik Ayak	40	1	25	20 x 30 x 20	910 gr.	ALH 2002
Klasik Ayak / Çift	40	2	12	20 x 30 x 20	885 gr.	ALH 2010
Klasik Bacak	65	1	25	20 x 30 x 20	1100 gr.	ALH 2003
Metal El	25	1	25	20 x 30 x 20	1215 gr.	ALH 2011
Metal El / Çift	25	2	12	20 x 30 x 20	1130 gr.	ALH 2012

# Fall Management Socks

Designed for fall risk patients.




- Patients requiring non-skid footwear to help prevent falls.
- Soft terry cloth.
- Skid resistant rubber tread on both sides.
- Available in 5 colors.
- 3 size: large, medium and small.  
Small: up to 38  
Medium: up to 45  
Large: up to 47
- Latex free.
- 30 °C washable.

## Fall Management Socks Ordering Information

Color	Size	Ordering Information
Red	Small Socks Fits Up too Youth 38	MS35SR
	Medium Sooks Fits Up too Adult 45	MS45MR
	Large Socks Fits Up too Adult 47	MS47LR
Orange	Small Socks Fits Up too Youth 38	MS35SO
	Medium Socks Fits Up too Adult 45	MS45MO
	Large Socks Fits Up too Adult 47	MS47LO
Purple	Small Socks Fits Up too Youth 38	MS35SP
	Medium Socks Fits Up too Adult 45	MS45MP
	Large Socks Fits Up too Adult 47	MS47LP
Green	Small Socks Fits Up too Youth 38	MS35SG
	Medium Socks Fits Up too Adult 45	MS45MG
	Large Socks Fits Up too Adult 47	MS47LG
Yellow	Small Socks Fits Up too Youth 38	MS35SY
	Medium Socks Fits Up too Adult 45	MS45MY
	Large Socks Fits Up too Adult 47	MS47LY

# Fall Precaution Bracelets

Use for identify patients at risk of easily fall.



- Produced from the Polymer material compatible to the medical devices material terms and standards\*.
  - Antiallergic, antibacterial and antifungus hygienic structure that doesn't irritate the skin, not composed bacteria and fungal on it.
  - Easy and rapid application.
  - Safety lock with its poppered system.
  - It has various sizes, length and width for different kind of purposes.
- \*Not used cancerogenic and PVC/DOP featured materials harmful to human health.

Fall Precaution Bracelets Ordering Information							
Model	Dimensidon	Color	Size	Piece/Box	Box Size	Box Weight	Ordering Information
DRHTB 1	Child	Yellow	2.5 x 24.2 cm	250	11 x 27 x 7 cm.	410 gr.	PSP - F01
DRHTB 2	Adult	Yellow	3.2 x 30 cm	250	11 x 33 x 7 cm.	500 gr.	PSP - F02



## Essential Fall Management Kit\*

- Easily identify patients at risk for falls.
- All components for high fall risk patient identification in one prepackaged kit.
- Takes the guesswork out of what products to use.

\*Please contact with our sales department for order informations of products.



## New Born Products



36

### **BİNDERMED**

Baby Limb Holder



36

### **NASOMED**

Baby Nasal Cannula Fixation Band



37

### **FIXERBOARD**

Baby Limb Boards



37

### **MAMASAFE**

Göğüs Koruma Padi



38

### **SCHÜTZERMED**

Phototherapy Eye Bands

# BINDERMED

## Baby Limb Holder



- Easy to use and pratique.
- Both use for hand and leg wrist.
- Skin touch surface absorbs the sweat.



- Covers the wrist safely and stable with its connection ribbons.
- Does not irritate the skin and cause allergy.
- Gamma sterilized.

Bindermed Ordering Information						
Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
Bindermed	3 x 11	1	25	11 x 20 x 25	263 gr.	BLH 3001

CE 1783

# NASOMED

## Baby Nasal Cannula Fixation Band



- Pratique and easy to use with its velcro bands.
- Can be used with every dimension and model of cannula.
- Covers the head smoothly, skin touch surface absorbs the sweat.
- Keeps the cannula safely at the babies nose getting through the hoops of the nasal cannula.
- Does not irritate the skin and cause allergy.
- Gamma sterilized.

Nasomed Ordering Information						
Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
Nasomed	43	1	25	11 x 20 x 25	245 gr.	NTH 5001

CE 1783

# FIXERBOARD

## Baby Limb Boards



- During the IV cannula application to the babies prevents the sprain of the wrist and support.
- Easy to use and pratique.
- Can be shaped to 3 dimension (downwards, upwards, right, left) by hand and can be stable at that shape.
- Both use for hand and leg wrist.
- Supports the wrists smoothly, cotton structured skin touch surface absorbs the sweat.
- Does not cause allergy and irritation.
- Prevents the sprain of the wrist.
- With its connection ribbons and velcro lock enables to fix the wrist safely.
- Gamma sterilized.

Fixerboard Ordering Information						
Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
Small	3 x 8,5	1	25	11 x 20 x 25	288 gr.	LHB 01
Medium	4 x 10	1	25	11 x 20 x 25	355 gr.	LHB 02
Large	4,3 x 12,3	1	25	11 x 20 x 25	415 gr.	LHB 03

CE 1783

# MAMASAFE

## Breast Protection Pad

The product is placed onto breasts in order to protect nipples from trauma or any physical contact.



- The product consists of two different materials: Annular polyurethane foam and overlaying napped cotton woven fabric.
- The product is absorbent and soft.
- It absorbs milk leakage and sweat owing to its absorbent quality.
- It does not cause irritation on the skin and is anti-allergic.

Mamasafe Ordering Information						
Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
4 Set	Inner Dia: 5 h:3	1 Pair + Strapless	4 Pairs	26 x 10 x 11	172 gr.	BPP 104
20 Set	Inner Dia:5 h:3	20 Pairs	20 Pairs	32 x 38 x 10	476 gr.	BPP 120

CE

# SCHÜTZERMED

## Phototherapy Eye Bands

- Prevents from the UV rays up to 450 Nm.
- Does not irritate the skin and cause allergy.
- Eye protective bands developed for the use of phototherapy treatment 4 different model.
- Gamma sterilized.



### SCHÜTZERMED Classic

- Over head banded and supported model.



- Absorbative cotton inner surface.
- 3 size: large, medium and small.



### SCHÜTZERMED Sticker

- Sticker model adhesive to the temples.



- 3 size: large, medium and small.



## SCHÜTZERMED Mask

- Model consists of head piece and protective mask.
- Adjustable for every size.



## SCHÜTZERMED Yota

- Model consists of elastic band and protective eye cover.
- Net structure allows the UV rays to the head.
- 3 size: large, medium and small.



### Festermed Ordering Information

Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
Klasik Small	5	1	25	11 x 20 x 25	265 gr.	PEB 301
Klasik Medium	31	1	25	11 x 20 x 25	298 gr.	PEB 302
Klasik Large	33	1	25	11 x 20 x 25	320 gr.	PEB 303
Sticker Small	12	1	25	11 x 20 x 25	240 gr.	PEB STC 301
Sticker Medium	17	1	25	11 x 20 x 25	248 gr.	PEB STC 302
Sticker Small	20	1	25	11 x 20 x 25	258 gr.	PEB STC 303
Mask	18 + 25	1	25	11 x 20 x 25	325 gr.	PEB TPD
Yota Small	24	1	25	11 x 20 x 25	218 gr.	PEB YT 301
Yota Medium	28	1	25	11 x 20 x 25	240 gr.	PEB YT 302
Yota Large	33	1	25	11 x 20 x 25	255 gr.	PEB YT 303



## Urology Products



42

**CLAMPMED**

Urine Control Clamp



43

**BAREMED**

Urine Preventer Band



43

**CONDOFIX**

Condom Fixer

# CLAMPMED

## Urine Control Clamp

The ideal solution is applied for 23 years in the incontinence problem!



### Comfort

Body structure of Clampedmed consists of semi-hard polymer of soft foam provides easiness and comfort during use.

With the balanced elastic bandage spreads its radial strength to the cellwall equally and reduces the clamp pressure to minimum level.

\*Antiallergic featured of Clampedmed is certified by the university report which is accredited by the Notify Body.

### Safe

When applied to the penis special designed cushion at the bottom part pressured to the urine system to reduce the incontinence and does not effects the blood circulation.

### Practical to use

In case of need easy to open without support and fix again.

### Healthy and Hygienic

Sterilized for use, produced from antiallergic special materials so no irritation will be caused to the skin.\* Also antibacterial and antifungal featured no bacteria and fungus will be composed\*\*.

### Confidant

Not occupy to much place unforeseen easy to hide. Without thinking its presence very easy to use in daily life provides confidence to the patient.

### Economic

Taking account of its using time, It is a more economical choice of product than absorbent ones.

\*\*Antibacterial and antifungus featured of Clampedmed is certified by the laboratuar report which is accredited by the Notify Body.

### Clampedmed Ordering Information

Model	Size	Piece / Pouch	Pouch/ Box	Box Size	Box Weight	Order Code
Clampedmed	Normal	1	1	4 x 10 x 16 cm	31 gr.	UCC 01
	Large	1	1	4 x 10 x 16 cm	32 gr.	UCC 02



# BAREMED

## Urine Preventer Band

To prevent the urine drip involuntary by pressing to the urethra.



### Practical

Easy and pratique to use.

Special pressure unit at the bottom of the band enables the user easiness.

Elastic wrap band enables extra pressure on to the urethra.

Connection ribbons attached to the sides enables to fix to the body safely.

Fits all size of a body with its elastic structure.



### Comfort

Skin touch surface is made of airpermeable structured antibacterial material.

### Safe to use

Does not hinder the blood circulation.

Does not cause any allergy and irritation.

Gamma sterilized.


### Baremed Ordering Information

Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
Baremed	3 x 12	1	50	11 x 20 x 25	235 gr.	UCB 6001



# CONDOFIX

## Condom Fixer



### Practical

- Fixes the male external urine catheter safely and securely in place.
- Single size for adult and children.
- Non sterile and single use.

### CondoFix Ordering Information

Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
CondoFix	3 x 12	5	10	15 x 25 x 21	325 gr.	CF 8001





## Wound Care Support Products



46

**FERSEMED**

Heel Protector



47

**FERSEMED Rg**

Foot Elevator



48

**WUNPAD**

Wound Care and Medical Dressing Pad

# FERSEMED

Heel Protector

It is used for heel protection decubitus ulcers (pressure sores) while the patient is reclining for a long term.



- It is used for heel protection decubitus ulcers (pressure sores) while the patient is reclining for a long term.
- Prevents the heel contact with sheets and prevents erosion of involuntary motion.
- Easy to use and simple.
- Wraps the heel comfortably and fixed by Velcro straps.
- Provides effective relief in patients.
- Can be washed at 30 °C and reuseable.

## Fersemmed Ordering Information

Model	Size* (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
Fersemmed	17 x 20	1	25	51 x 35 x 50	1750 gr	HPS 101



# FERSEMED Rg

Foot Elevator

This foot elevation device is designed to reduce the risk of pressure ulcers (pressure sores) while the patient is in bed or reclining for a long term.



- It fits around the ankle to raise the heel off of any surface.
- It has a spiral-cut foam ring which easily fits over liner that provides foot elevation and comfort and can be fixed to the wrist with Velcro locks.
- Cylindrical body allows the patient's feet to the right to the left and back and forth.
- With its sponge structure provides softness and comfort.
- Cotton finished cover will not harm the skin.
- The product can be cleaned by hand or machine washed at 30 °C.

## Fersemmed Rg Ordering Information

Model	Size* (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
Fersemmed Rg	R: 24 x h: 10	1	1	25 x 11 x 25 cm	304 gr.	HPC 101
	R: 24 x h:15	1	1	25 x 16 x 25 cm	344 gr.	HPC 102



# WUNPAD

## Wound Care and Dressing Pad

Wound care and dressing pad for emergency situation for bleeding cuts and wounds.



### Atraumatic

Produced from hydro fiber structured organic and inorganic materials in nonwoven form, atraumatic and soft.

### Pratique and advantageous

Absorption capacity 8 times of its weight.

Do not adhere to the sore, do not leave fibre, feather, particle.

Enables the user options to choose with its different dimensions.

### Sterile

Gamma sterilized.

#### Wunpad Wound Care and Dressing Pad Ordering Information

Model	Size* (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
Wunpad	5 x 10	1	25	10 x 11 x 16	145 gr.	WP 5x10
	8 x 12	1	25	10 x 16 x 17	175 gr.	WP 8x12
	10 x 10	1	25	10 x 16 x 17	175 gr.	WP 10x10
	10 x 20	1	25	10 x 16 x 23	240 gr.	WP 10x20
	15 x 15	1	25	10 x 21 x 27	280 gr.	WP 15x15

\*Upon request special sizes can be produced.

CE 1783

## Patient Safety Products



51

### **SAFEBELT**

Bed & Stretchers Safety Belts



52

### **SAFEBELT**

Chair Safety Belts



53

### **SAFEBELT**

Patient Transfer Safety Belts



54

### **SAFEBELT**

Support Belts



55

### **SAFEBELT**

Accessories



57

### **SAFEBELT**

Bed Safety Belts Group



58

### **SAFEBELT**

Wheelchair Safety Belts



59

### **SAFEBELT**

Patient Support Belts

# SAFEBELT

Patient Safety Belts

## Steel Lock Serie



Safebelt patient safety belts used to control the movements of the patients with high safety risk (mental patients, alcohol addicts and patients at risk for escape) that might harm to their environments.

It has a wide product range to prevent and limited all undesired movements to respond all cases.

### General Specifications

- By providing a secure fixation offers a maximum safety both for patients and around.
- Provides maximum freedom sensation to the patients where the patient movement is restricted as much as necessary.
- Application is easy and practical through effective magnetic switches.
- They are quite reliable with special steel pin locks.
- According to their usage areas, are divided into 3 main groups;
  - Used bad and stretcher
  - Used in chair
  - Used in patient transfer
- Products are made of 100% terylene and polyester material.
- 2 years warranty.
- Washable at 95 °C.



# SAFEBELT

## Bed & Stretchers Safety Belts

This group belts are used in order to fix patients securely to the bed and stretcher.



### Waist Roll Belt

- Provides secure bedding for the patients.
- Provides freedom to turn to sides besides lying back and sitting positions.
- Meanwhile since patients are free to use their hands and legs, they can move comfortably for their social or personal needs (reading, eating etc.).




### Shoulder Restraints

- Prevents shoulder and body movement and fixes the patients laid back. Used together with Waist Roll Belt or Multi - Function Extension Belt.



### Wrist Restraints

- Securely fixes the wrist of the patient to the bed frame and restricts the movement certainly.
- Provides comfort for the patients due to its aesthetic and safe design.



### Ankle Restraints

- Securely fixes the ankle of the patient to the bed frame and restricts the movement certainly.
- By using together with the extension, provides limited restriction for the patient's leg movement.
- Provides comfort for the patients due to its aesthetic and safe design.



### Full Body Crisis Set

- Used to provide movement restriction for the patients during crisis.
- Consist of the items below;  
2 pcs. wrist restraints, 2 pcs. ankle restraints.  
Waist roll belt  
Shoulder restraints  
Multifunctional extension belts

# SAFEBELT

## Chair Safety Belts

This group belts are used in order to fix patients securely on the chair.



### Wheelchair Waist Belts

- Secures the patient sitting on the wheelchair from slipping.
- If necessary can be used together with Pelvic Belt.



### Wheelchair Shoulder Restraints


- Secures the patient sitting on the wheelchair from bending forward or backward also slipping.



# SAFEBELT

## Patient Transfer Safety Belts

Belt is used to restrict the patient's arm and leg movements while being taken from one place to another place. Belt is composed of three parts; including hand, foot and link belts.



### Wrist Movement Restriction Belt

It is formed from two wrist belts and a belt wrapped around the waist and sewn onto the belt conveyor. Wrist belts wrap securely and will not harm the skin ,conveyor belt ensure the hands at waist height by fixing the patient and prevents damage to itself and environment .




### Ankle Restriction Belts

It is formed from conveyor belt and two wrist belts sewn onto the belt .

Wrist belts wrap the ankle comfortably without any pressure and in confidence that will not harm the skin.

The conveyor belt controls the patient's walking speed and stride length and prevents damage to the environment by restricting leg movement.



### Binding Band for Wristbands and Ankle Bands

Binding band is used to join wristband set with ankle band set.

# SAFEBELT

## Support Belts

Supporting belts provides extra safety and comfort and used together with main safety belts.



### Multi - Function Belt

Used together with wrist and ankle restraints to provide limited movement restriction.



### Standard Extender

Used together with other main belts where longer ends needed.



### Waist Extender

Provides extension for the waist belts for XXL size patients.



### Pelvic Extender

Provides extra safety to the patients, used together with the waist roll belt and provides additional support in sitting position as well.

#### Safebelt Support Belts Ordering Information

Model	Size*	Order Code
Waist Roll Belt	S-M (60-110 cm.)	SBB101
	L (85-135 cm.)	SBB102
Wrist Restraints	Standart (14-26 cm.)	SBB201
Wrist Restraints Set	Standart (14-26 cm.)	SBB202
Ankle Restraints	Standart (16-32 cm.)	SBB301
Ankle Restraints Set	Standart (16-32 cm.)	SBB302
Ankle & Waist Restraints Set	Standart (14/16-24/32 cm.)	SBB605
Shoulder Restraints	Standart (216 cm.)	SBB401
Wheelchair Shoulder Restraints	Standart (216 cm.)	SBC 102
Multi Function Belt	Standart (180 cm.)	SBS501
Standard Extender	Standart (60 cm.)	SBS502
Waist Extender	Standart (25 cm.)	SBS503
Pelvic Belt	Standart (95 cm.)	SBS504
Standart Crisis Set	S-M	SBB603
	L	SBB604
Full Body Crisis Set	S-M	SBB601
	L	SBB602
Wheelchair Waist Belt	Standart (150 cm.)	SBC101
Wrist Movement Restriction Belts	Standart	SBT101
Ankle Restriction Belts	Standart (140 cm.)	SBT201
Binding Band for Wrist Bands and Ankle Bands		SBT501

# SAFEBELT

## Accessories



### Steel Lock

- Lock with magnetic control, special design, made of steel, can be unlocked by a special key.



### Magnetic Key

- Easy to use, magnetically controlled key.



### Bag

- Designed in 2 sizes can be able to carry complete set of crisis or small enough to carry single sets as separated.
- Existing of special cases which the locks and keys easily controlled.

#### Safebelt Accessoires Ordering Information

Model	Size* (cm)	Order Code
Steel Lock	Standart	SBA701
Magnetic Key	Standart	SBA702
Bag	Standart	SBA703

# SAFEBELT

Patient Fixation Belts

## Buckled Soft Serie



Safebelt patient fixation belts used to control the movements in sensitive and precisely of the patients with secondary safety risk that might harm themselves. It has a wide product range to prevent and limited all undesired movements to respond all cases.

### General Specifications

- They provide maximum patient safety to itself and the environment by fixing in safe.
- They provide comfort and sense of freedom to the user even is restricted or controlled at the time of movement.
- According to their usage areas, are divided into 2 groups used in the chair and used in the bed.
- They are applied easily and practical because of their front fastening straps made of Velcro and shatter - resistant material made of POM buckled locks.
- Products are made of %100 terylene and polyester material.
- Can be washed at 95 °C.
- 2 years warranty.


# SAFEBELT

## Bed Safety Belts Group

This group belts are used in order to fix elderly patients in a safe and comfortably way to the bed.



### Wrist Safety Belt



- Constants patient's wrist to the bed frame and restrict movement completely or partially.
- Provide comfort to the patients through its convenient model and aesthetic.
- Doesn't harm the skin cause to be made of soft material and provides comfort to the patient.



### Ankle Safety Belt



- Constants patient's ankle to the bed frame and restrict movement completely or partially.
- Provide comfort to the patients through its convenient model and aesthetic.
- Doesn't harm the skin cause to be made of soft material and provides comfort to the patient.

# SAFEBELT

## Wheelchair Safety Belts

This group belts are used in order to fix elderly patients securely and comfortably to the wheelchair.



### Wheelchair Safety Belt

- Prevents sliding and enables leaning forward for the patient sitting in the wheelchair.
- Provides support between a patient's leg with soft filled portion in the front.

#### Safebelt Bed Safety Belts Ordering Information

Model	Size* (cm)	Order Code
Wrist Fixation Belt	Buckled (45 cm)	YBM 201
Ankle Fixation Belt	Buckled (65 cm)	YBM 301
Wheelchair Fixation Belt	Buckled (164 -146 cm)	YBS 102



# SAFEBELT

## Patient Support Belts

Safebelt Patient Support Belts are designed for patients who need support and at risk of falling.



### Walking Belts

- Walking belt helps to support the patients who have difficulty to walk alone.
- Soft and durable belt has got six handles as being sewn on the vertical and horizontal ergonomic grip straps are available.
- By assistance of those handles on the belt wrapped above the waist helps prevent caregiver back injuries and assists in the safe transfer, stumble and fall risk is avoided.
- The belt is easily to grasp the patient's waist by its enough width.
- By courtesy of plastic buckle can be adjustable according to each body, opening and closing is simple, fast and secure.
- It can be cleaned by hand or machine washed at 95 °C.



### Lift Assist Belt

- The belt used to lift up the patient from bed or wheelchair.
- Ergonomic hand grips are mounted at both ends of the belt.
- The belt is wrapped around the patient's waist and given support to try to pull up, go back (fall risk) is blocked.
- It was made of durable soft, special flexible neoprene material.
- It can be easily carried by each officer because of lightness and takes up little space.
- If necessary, hand or machine washable at 30 °C.

#### Safebelt Patient Support Belts Ordering Information

Model	Size* (cm)	Order Code
Walking Support Belt	Standart (140 cm)	SBS 101
Patient Lift Assist Belt	Standart (46 cm)	SBL 201



## General Products



62

### **UNI-GLOVE**

Operating Room Light Handle Cover



63

### **PASSPOP**

Identification Band

# UNI-GLOVE

Operating Room Light Handle Cover

Practical and economical solution for operating room light system\* handles.



## Functional

Eliminates the need to remove, clean, sterilize and reconnect the universal adaptors of light systems in the operating room after every operation.

Thank to its circle shaped guard, prevents the light system against any stains caused by accidental touch of the hand.

## Useful

Fits on all currently used universal adapters in the market.

Fits on handles of light systems easily and securely.  
Easily removed and disposed after the operations.

## Ekonomical

Saves cleaning and sterilization time after operations.

Saves energy required for sterilization process.

Eliminates the risk of wear on handles and light systems due to constant removal and reconnection.

Eliminates the ever increasing and costly need to purchase and stock adaptors and handles.

## Environmental friendly

Made of enviromental friendly plastic.

\* Patented


### Uni-Glove Ordering Information

Model	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
Uni-Glove	Total h: 16 Envelope ø: 6.5 Cover ø: 12.5	1	25	20 x 31.5 x 17	660 gr.	UG 01

# PASSPOP

## Identification Band

## New model PASSPOP with popper Identification Band.



### Practical

Easily applicable like similar standard bands with popper.

### Hygienic

Antibacterial and antifungus featured. Prevents to compose bacteria and fungal\*.

### Antiallergic

In case of long term use on patients arm does not cause allergy and irritation\*\*.

### Useful

Barcode label is indelible, cannot scratched, cannot wet, ables to patient even take a bath.

### Safe

Cannot take out of the arm by stretching.  
Unavailable for the use of another patient after take out.

### Advantage of dimension

Provides a space to stick the barcode label easily.  
More long and wide from the similar ones.

### Special Marking

Due to request logo and name of the user firm can be written on the band.



\*Antibacterial and antifungal featured of Passpop is certified by the laboratuar report which is accredited by the Notify Body.

\*\* Antiallergic featured and nonirritation of Passpop is certified by the university report which is accredited by the Notify Body.

### Passpop Ordering Information

Boy	Renk	Size (cm)	Piece / Pouch	Pouch/ Box	Box Size (cm)	Box Weight	Order Code
Child	White	2.5 x 24.2	250	-	11 x 27 x 7 .	410 gr.	PSP - B01
Adult		3.2 x 30	250	-	11 x 33 x 7 .	500 gr.	PSP - B02
Child	Red	2.5 x 24.2	250	-	11 x 27 x 7 .	410 gr.	PSP - K01
Adult		3.2 x 30	250	-	11 x 33 x 7 .	500 gr.	PSP - K02
Child	Yellow	2.5 x 24.2	250	-	11 x 27 x 7 .	410 gr.	PSP - S01
Adult		3.2 x 30	250	-	11 x 33 x 7 .	500 gr.	PSP - S02
Child	Navy	2.5 x 24.2	250	-	11 x 27 x 7 .	410 gr.	PSP - L01
Adult		3.2 x 30	250	-	11 x 33 x 7 .	500 gr.	PSP - L02
Mother-Baby	Pink	3x27.5 - 2.2x18	125	-	11 x 33 x 7 .	400 gr.	PSP - P.AB
Mother-Baby	Blue	3x27.5 - 2.2x18	125	-	11 x 33 x 7 .	400 gr.	PSP - M.AB





Medikokim Tıbbi Kimyevi Malzeme San. Tic. Ltd. Şti  
İkitelli O.S.B. Demirciler Sanayi Sitesi B/7 Blok No:157  
Başakşehir PK 34306 İstanbul / TÜRKİYE  
Phone: +90 212 549 57 78 • Fax: +90 212 549 57 79  
E-mail: info@medikokim.com

[medikokim.com](http://medikokim.com)

Rev.: 03