


## ORTHOPEDIC DRILL BITS




BJ11021

### Standard straight shank

Length(mm)	Connect/Shank Length(mm)	Effective Cutting Length(mm)	Head Diameter(mm)
115	80	35	Φ1.5
115	80	35	Φ2.0
115	60	55	Φ2.5
115	60	55	Φ2.7
115	60	55	Φ3.2
115	60	55	Φ3.5
115	60	55	Φ4.0
115	60	55	Φ4.2
115	60	55	Φ4.5
115	60	55	Φ5.0

lot 47  
lot 50




BJ41021SAO

### AO Drill bits (Synthes quick connect coupling)

Length(mm)	Connect/Shank Length(mm)	Effective Cutting Length(mm)	Head Diameter(mm)
125	100	25	Φ1.1
125	100	25	Φ1.5
125	100	25	Φ1.8
125	95	30	Φ2.0
125	95	30	Φ2.5
125	95	30	Φ2.7
125	85	40	Φ3.2
125	85	40	Φ3.5
125	85	40	Φ4.0
125	85	40	Φ4.5

## ORTHOPEDIC DRILL BITS



### Standard straight shank cannulated drill bits

Length (mm)	Connect/Shank Length (mm)	Effective Cutting Length (mm)	Head Diameter (mm)	Inner cannulation (mm)
150	120	30	Φ2.7	Φ1.4
BJZCZ-2.7-30-150				
150	115	35	Φ3.2	Φ1.7
BJZCZ-3.2-35-150				
150	100	50	Φ3.8	Φ1.7
BJZCZ-3.8-50-150				
150	130	20	Φ2.0	Φ1.1
BJZCZ-2.0-20-150				
150	130	20	Φ2.2	Φ0.9
BJZCZ-2.2-20-150				
250	210	40	Φ4.5	Φ2.1
BJZCZ-4.5-40-250				
300	230	70	Φ6.0	Φ2.7
BJZCZ-6.0-70-300				