



常州金龙医用塑料器械有限公司
Changzhou Jinlong Medical Plastic Appliance Co.,Ltd

ADDRESS: NO.6 ZHENNAN ROAD, ZHENGLU TOWN,
CHANGZHOU CITY
TEL: +86-519-85785001 **FAX:** +86-519-88737581
EMAIL: JINLONGMED@JINLONGMED.COM.CN/ ZLY@
JINLONGMED.COM.CN
WEBSITE: WWW.JINLONGMED.COM

地址：中国江苏省常州市郑陆镇南路6号
电话：+ 86-519-85785001 传真：+86-519-88737581
电子邮件：jinlongmed@jinlongmed.com.cn/ zly@
jinlongmed.com.cn
网站：www.jinlongmed.com



**JINLONG
MEDICAL**
PLASTIC APPLIANCE CO.,LTD

常州金龙医用塑料器械有限公司
Changzhou Jinlong Medical Plastic Appliance Co.,Ltd



简介 公司 About us

Established in 1992, CHANGZHOU JINLONG MEDICAL PLASTIC APPLIANCE CO., LTD. is located in Zhenglu Town, Changzhou, Jiangsu Province, China. Our company is a professional medical device enterprise approved by National Medical Products Administration. The company covers an area of 15000 m², including 3813 m² purification workshops. We have 200 staff members, 20 administrators, 10 engineers and technicians. The annual sales amount reaches to fifty million CNY. Our main products are as follows: Disposable sterile syringe; Disposable infusion set; Disposable blood transfusion; Disposable hypodermic needle; Disposable scalp vein set and Disposable endoscopic biopsy forceps. Our company roundly actualizes the quality management system of ISO13485:2016. We waited through TUV ISO13485:2016, CE Quality Authentication. Nowadays, our products have been exported to various countries and regions in Mid-Asia, Mid-East, South east Asia, Africa, South America, Europe and so on, and have won a good reputation.

常州金龙医用塑料器械有限公司成立于一九九二年，位于中国江苏省常州市郑陆镇，是经国家药品监督管理局批准设立的医疗器械专业企业。公司占地面积为15000平方米，其中十万级净化面积3813平方米。现有员工200人，管理人员30人，工程技术人员10人。拥有固定资产2000万，年销售额达到5000万人民币。公司全面实施和运行ISO13485:2016质量管理体系，并通过德国莱茵公司EN ISO13485:2016和CE产品质量认证。

公司主导产品有“畅怡”牌一次性使用无菌注射器带针、一次性使用配药用注射器、一次性使用无菌注射针、一次性使用输液器带针、一次性使用输血器带针、一次性使用静脉输液针、一次性使用内窥镜活检取样钳。目前公司产品远销中亚、中东、东南亚、欧洲、非洲、南美等多个国家和地区，并获得了良好的声誉。公司全体员工以“客户至上、信誉第一”为宗旨，本着诚信、务实、开拓、敬业的经营作风，充分利用广泛的国外客户资源，互惠互利，为实现企业的繁荣和发展，为回报支持我们的合作伙伴，不懈努力，共同发展。



Changzhou Jinlong Medical Plastic Appliance Co., Ltd



DISPOSABLE HYPODERMIC NEEDLE

一次性使用无菌注射针



DISPOSABLE INTRAVENOUS NEEDLE (SCALP VEIN NEEDLE)

一次性使用静脉输液针—头皮针



规格		色 标		常用针管长度 Regular Needle Length				
				15mm 5/8"	20mm 3/4"	25mm 1"	32mm 1 1/4"	38mm 1 1/2"
0.D	Gauge	Colour Code						
0.40	27G	中灰色	grey					
0.45	26G	褐色	brown					
0.50	25G	橙色	orange					
0.60	23G	深蓝色	blue					
0.70	22G	黑色	black					
0.80	21G	深绿色	green					
0.90	20G	黄色	yellow					
1.20	18G	粉红色	pink					
1.60	16G	白色	white					

ADDRESS: NO.6 ZHENNAN ROAD, ZHENGLU TOWN, CHANGZHOU CITY.

TEL: +86-519-85785001

FAX: +86-5189-88737581

EMAIL: JINLONGMED@JINLONGMED.COM.CN/ ZLY@JINLONGMED.COM.CN WEBSITE: WWW.JINLONGMED.COM



JINLONG MEDICAL

3-PART DISPOSABLE SYRINGE WITH LUER LOUK
三件套一次性使用无菌注射器—螺口式



3-PART DISPOSABLE SYRINGE WITH LUER SLIP
三件套一次性使用无菌注射器—直口式



IRRIGATION SYRINGE 冲洗器



ADDRESS: NO.6 ZHENNAN ROAD, ZHENGLU TOWN, CHANGZHOU CITY.

TEL: +86-519-85785001

FAX: +86-5189-88737581

EMAIL: JINLONGMED@JINLONGMED.COM.CN/ ZLY@JINLONGMED.COM.CN WEBSITE: WWW.JINLONGMED.COM



JINLONG MEDICAL

Disposable infusion set 一次性使用输液器




Urine Bag 尿袋



HUB 6003 No outlet
单通



HUB 6005 Pull push valve
推拉阀



HUB 6006 Screw valve
螺口阀



HUB 6007 T cross valve
十字阀

Disposable Face Mask 一次性使用口罩



ADDRESS: NO.6 ZHENNAN ROAD, ZHENGLU TOWN, CHANGZHOU CITY.

TEL: +86-519-85785001


FAX: +86-519-88737581

EMAIL: JINLONGMED@JINLONGMED.COM.CN/ ZLY@JINLONGMED.COM.CN WEBSITE: WWW.JINLONGMED.COM



JINLONG MEDICAL

Disposable Cytology Brush 一次性使用细胞刷

	Model 型号	Balloon O.D.(mm) 外管直径	Brush head O.D.(mm) 刷头直径	Length(mm) 有效长度	Working channel (mm) 适用钳道孔
	XBS-1.8 x (600-2500)	1.8	2/ 2.5/ 3/ 4/ 5	1300-2800	$\geq \phi 2.0$
	YWQ-2.3 x (600-2500)	2.3	2/ 2.5/ 3/ 4/ 5	1300-2800	$\geq \phi 2.8$

Disposable Injection Needles 一次性使用内窥镜注射针

	Model 型号	diameter of the outer tube(mm) 外管直径	length(mm) 有效长度	The diameter of needle (mm) 针管外径	The length of needle(mm) 针管长度	Working channel(mm) 适用钳道孔
	ZSZ —1.8 x (300-2500) x 25 x 4	1.8	300-2500	25G	4	$\geq \phi 2.0$
	ZSZ—2.3 x (300-2500) x 23x 5	2.3	300-2500	23G	5	$\geq \phi 2.8$

ADDRESS: NO.6 ZHENNAN ROAD, ZHENGLU TOWN, CHANGZHOU CITY.

TEL: +86-519-85785001

FAX: +86-519-88737581

EMAIL: JINLONGMED@JINLONGMED.COM.CN/ ZLY@JINLONGMED.COM.CN WEBSITE: WWW.JINLONGMED.COM



JINLONG MEDICAL

Disposable grasping forceps 一次性使用内窥镜异物钳

Classification 类别	Model 型号	Outer Diameter(mm) 外管直径	Length(mm) 长度	Working channel (mm) 适用钳道孔
	YWQ-3SZ-1.8 x (300-2500)	1.8	300-2500	$\geq \phi 2.0$
	YWQ-3SZ-2.3 x (300-2500)	2.3	300-2500	$\geq \phi 2.8$
	YWQ-4SZ-1.8 x (300-2500)	1.8	300-2500	$\geq \phi 2.0$
	YWQ-4SZ-2.3 x (300-2500)	2.3	300-2500	$\geq \phi 2.8$
	YWQ-5WZ-1.8 x (300-2500)	1.8	300-2500	$\geq \phi 2.0$
	YWQ-5WZ-2.3 x (300-2500)	2.3	300-2500	$\geq \phi 2.8$
	YWQ-EK-1.8 x (300-2500)	1.8	300-2500	$\geq \phi 2.0$
	YWQ-EK-2.3 x (300-2500)	2.3	300-2500	$\geq \phi 2.8$
	YWQ-SC-1.8 x (300-2500)	1.8	300-2500	$\geq \phi 2.0$
	YWQ-SC-2.3 x (300-2500)	2.3	300-2500	$\geq \phi 2.8$

ADDRESS: NO.6 ZHENNAN ROAD, ZHENGLU TOWN, CHANGZHOU CITY.

TEL: +86-519-85785001

FAX: +86-519-88737581

EMAIL: JINLONGMED@JINLONGMED.COM.CN/ ZLY@JINLONGMED.COM.CN WEBSITE: WWW.JINLONGMED.COM



JINLONG MEDICAL

Disposable endoscopic biopsy forceps

一次性使用内窥镜活体取样钳

Jaws 钳头	Model 型号	Specification 规格详情	Diameter of jaws (mm) 钳头直径	Length (mm) 长度	For working channel 适用内窥镜钳道孔
	QYQ-FA 1.8*(300-2500)	Disposable biopsy forceps without PE coated, without spike, oval/alligator cup 一次性使用内窥镜活体取样钳 不包塑, 不带针, 平口/齿口	1.8	300-2500	$\geq \Phi 2.0\text{mm}$
	QYQ-FA 2.3*(300-2500)		2.3	300-2500	$\geq \Phi 2.8\text{mm}$
	QYQ-FB 2.3*(300-2500)	Disposable biopsy forceps without PE coated, with spike, oval/alligator cup 一次性使用内窥镜活体取样钳 不包塑, 带针, 平口/齿口	2.3	300-2500	$\geq \Phi 2.8\text{mm}$
	QYQ-FC 1.8*(300-2500)	Disposable biopsy forceps with PE coated, without spike, oval/alligator cup 一次性使用内窥镜活体取样钳 包塑, 不带针, 平口/齿口	1.8	300-2500	$\geq \Phi 2.0\text{mm}$
	QYQ-FC 2.3*(300-2500)		2.3	300-2500	$\geq \Phi 2.8\text{mm}$
	QYQ-FD 2.3*(300-2500)	Disposable biopsy forceps with PE coated, with spike, oval/alligator cup 一次性使用内窥镜活体取样钳 包塑, 带针, 平口/齿口	2.3	300-2500	$\geq \Phi 2.8\text{mm}$


ADDRESS: NO.6 ZHENNAN ROAD, ZHENGLU TOWN, CHANGZHOU CITY.

TEL: +86-519-85785001

FAX: +86-519-88737581

EMAIL: JINLONGMED@JINLONGMED.COM.CN/ ZLY@JINLONGMED.COM.CN WEBSITE: WWW.JINLONGMED.COM

Disposable grasping forceps 一次性使用内窥镜异物钳

Classification 类别	Model 型号	Outer Diameter(mm) 外管直径	Length(mm) 长度	Working channel (mm) 适用钳道孔	Opening width (mm) 开幅
	YWQ-BY-1.8 x (300-2500)	1.8	300-2500	≥φ2.0	15, 20, 25, 30
	YWQ-BY-2.3 x (300-2500)	2.3	300-2500	≥φ2.8	15, 20, 25, 30
	YWQ-TY-1.8 x (300-2500)	1.8	300-2500	≥φ2.0	15, 20, 25, 30
	YWQ-TY-2.3 x (300-2500)	2.3	300-2500	≥φ2.8	15, 20, 25, 30
	YWQ-LB-1.8 x (300-2500)	1.8	300-2500	≥φ2.0	15, 20, 25, 30
	YWQ-LB-2.3 x (300-2500)	2.3	300-2500	≥φ2.8	15, 20, 25, 30
	YWQ-ZJ-1.8 x (300-2500)	1.8	300-2500	≥φ2.0	15, 20, 25, 30
	YWQ-ZJ-2.3 x (300-2500)	2.3	300-2500	≥φ2.8	15, 20, 25, 30
	YWQ-TY-1.8 x (300-2500)	1.8	300-2500	≥φ2.0	15, 20, 25, 30
	YWQ-TY-2.3 x (300-2500)	2.3	300-2500	≥φ2.8	15, 20, 25, 30
	YWQ-WD-1.8 x (300-2500)	1.8	300-2500	≥φ2.0	20 x 50
	YWQ-WD-2.3 x (300-2500)	2.3	300-2500	≥φ2.8	30 x 60

ADDRESS: NO.6 ZHENNAN ROAD, ZHENGLU TOWN, CHANGZHOU CITY.




TEL: +86-519-85785001 FAX: +86-519-88737581

EMAIL: JINLONGMED@JINLONGMED.COM.CN/ ZLY@JINLONGMED.COM.CN WEBSITE: WWW.JINLONGMED.COM



JINLONG MEDICAL

Endoscopy series 内镜微创系列



ADDRESS: NO.6 ZHENNAN ROAD, ZHENGLU TOWN, CHANGZHOU CITY.

TEL: +86-519-85785001

FAX: +86-519-88737581

EMAIL: JINLONGMED@JINLONGMED.COM.CN/ ZLY@JINLONGMED.COM.CN WEBSITE: WWW.JINLONGMED.COM