

# ACUSON Sequoia Ultrasound System

## Installation Checklist

Basic Applications – General System 3.5 (VB30)

[siemens-healthineers.com/ultrasound](http://siemens-healthineers.com/ultrasound)



This checklist represents features discussed during system training by the Application Specialist/Clinical Sales Specialist. It is recommended that the customer read and review all sections of the manual for further knowledge and explanation of the system and its usage.

✓	General System
	Power On/Off/Hibernate/Breaker/Power Cord
	Wheels/Braking System
	Monitor Tilt/Adjustable Arm/Lock Lever
	Touch Screen Tilt/Lighting Adjustment
	Control Panel Height/Rotate/Lighting Adjustment
	Foot Switch/Physio Module/Gel Warmer
	Transducers/Ports/Transducer Holders (Connect and Remove) and Cord Management

	Accessing the Retractable Keyboard (Optional)
	USB Ports (Left Side of Touch Screen, Back Location for Peripherals)
	Peripherals (B/W Printer, DVD/CD, Fusion)
	Input/Output/Panel Connections/Universal Video Converter (Off System)
	Storage Bins (Front, Back)
	Cleaning the Filter
	Safety and Care (System, Transducers, Touch Screen, Trackball)
	Recommended Cleaners and Disinfectants

# ACUSON Sequoia Ultrasound System

## Installation Checklist

### Basic Applications – General System 3.5 (VB30)

This checklist represents features discussed during system training by the Application Specialist/Clinical Sales Specialist. It is recommended that the customer read and review all sections of the manual for further knowledge and explanation of the system and its usage.

✓	<b>System Settings</b>
	Set Date/Time/Confirmation Beep
✓	<b>Transducers</b>
	Gesture Detection
	Orientation/Changing Frequencies
	Frequency Range/Depth
	Field of View (FOV)
	Wide Field of View (WFOV) – 5C1/11M2
	Applications/Exam Types
	Needle Guides (Option)
✓	<b>Image Screen &amp; Image Settings</b>
	Screen Saver
	Patient Banner
	Dual: Entering Full Screen – Configuration/Workflow
	Capture Settings
	Imaging Parameters (Customize under System Configuration > Imaging Settings > Image Parameter Display)
	Transducer/Exam/System Status

	Orientation Marker/Measured Results
	Status Indicators (Cine, Trackball)/Thumbnails
	Network Connection and System Status Display
	Printing and Deleting Thumbnails Controls
✓	<b>Touch Screen Layout</b>
	Customize Touch Screen User-defined
	Log-in Options (Person Icon)/ Help (Question Mark Icon)/ System Configuration Menu (Cog Icon)
	Workflow Sections: Fusion/Panoramic/ Protocols/Modality Compare/Physio
	Additional Pages of Selections
	Soft Key Functions and Corresponding Press and Rotate Controls
	Page Selection/Active Page
	Main Navigation Selections (Patient, Imaging, Review, Report, End Exam)
	Additional Selections Icon/Swipe Gesture
	White Lighting: Available Selection Dithered: Not Available
	Orange Lighting: Active Selection

# ACUSON Sequoia Ultrasound System

## Installation Checklist

### Basic Applications – General System 3.5 (VB30)

This checklist represents features discussed during system training by the Application Specialist/Clinical Sales Specialist. It is recommended that the customer read and review all sections of the manual for further knowledge and explanation of the system and its usage.

✓ <b>Control Panel/User Interface</b>	
Power/Hibernate	Activating an Operating Mode/CINE
DGC Controls	Capturing Images and Clips
Home Base Controls (Clip, Image, ABC, Caliper)	Completing a Patient Study
Trackball and Trackball Keys (Set, Update)	✓ <b>Workflow Protocols</b>
Press and Rotate Keys (Correspond with Soft Key Functions on Touch Screen)	Configure protocols under System Configuration > Features > Protocols
Operating Modes (2D, C, PW/CW, M, VT, CEUS, 3D/4D)	Protocol Selections/Workflow
Depth, Zoom, Dual, Freeze, Forward/Backward, Delete, Pointer	Importing and Exporting Protocols
User-defined (U1/U2) controls	Configuring Protocols on a Personal Computer
Retractable Keyboard (Optional)	✓ <b>2D Image Optimization</b>
Special Characters (Alt, Alt Gr)	2D/Overall Gain/DGC/Auto TEQ
Caps Lock/Text Start (Space Bar)	Speed of Sound/UltraArt/Harmonics/Compounding
Blank (Displays Keyboard Language)	U/D Flip/L/R Flip/Live Dual/Timer/Biopsy/Image Size
Beginning a Study	Clarify/Frequency/Line Density/DR/Steer
Preregister a Patient/Register a New Patient/ Register from Worklist	In Focus/Selectable Focus
Entering/Correct Patient Registration Information	Transmit Power/Persist/Maps/Tint
Selecting Transducer/Exam	Clip Length/Capture Type (Prospective, Retrospective)
Appending a Study/Selecting an Exam Package	2D Post-processing/Create or Modify User-defined Exams
Pause/Resume/Discontinue exam	

# ACUSON Sequoia Ultrasound System

## Installation Checklist

### Basic Applications – General System 3.5 (VB30)

This checklist represents features discussed during system training by the Application Specialist/Clinical Sales Specialist. It is recommended that the customer read and review all sections of the manual for further knowledge and explanation of the system and its usage.

✓	Color Doppler/Power Doppler/Slow Flow
	Color/Power/Auto Color/Directional Power/Slow Flow
	Resize/Reposition Color Region of Interest (ROI)
	2D + Color/2D Only/Color Only
	L/R Flip/U/D Flip/Biopsy/Invert/Live Dual
	Optimize for Flow States/Anatomy (Low, General, High, Kidney, Aorta, Slow Flow)
	Filter/Freq/Line Density/Tint/Scale/Steer
	Smooth/Baseline/Persist/Maps/Priority
	Color Post-processing
	Velocity Variance Color Map
✓	Pulsed-waved (PW) or Continuous Wave (CW) Doppler
	PW/CW/Doppler Gain/Format
	Gate Size/Sweep Speed/Baseline/Scale/Angle Correct/Steer/UltraArt
	Auto Spectral/Configure Auto Spectral > System Configuration > Imaging Settings > Imaging Mode
	Freq/Filter/DR/Maps/Tint
	Invert/Edge/HPRF/Audio Volume
	Display Formats – Live/Live (Triplex)
	Spectral Doppler Post-processing
	Auto Trace/Sensitivity

✓	M-mode
	M/M-mode Gain
	Anatomical M-mode/Color M-mode/Format
	Sweep Speed/DR/Steer/Maps/Tint/Edge/Harmonics/Display Formats
✓	Text/Arrows/Body Markers
	Press ABC to Activate Text
	Using Pre-defined Text on Touch Screen Text A–D, Additional Text
	Clear Screen/Delete Line/Previous Annotation/Home/Set Home/Delete Key on Control Panel
	Reposition/Edit/Overwrite Text
	Customizing Text > System Configuration > Annotations > Text
	Press ABC > Arrow on Touch Screen
	Using Arrow Widget/Adding/Deleting/Moving Arrows/Clear Screen/Delete/Previous Arrows
	Press ABC > Body Marker on Touch Screen
	Select/Reposition Body Marker/Transducer Marker Orientation and Position/Delete/Widget/Exam Package
	Customizing Body Marker > System Configuration > Annotations > Body Markers

# ACUSON Sequoia Ultrasound System

## Installation Checklist

### Basic Applications – General System 3.5 (VB30)

This checklist represents features discussed during system training by the Application Specialist/Clinical Sales Specialist. It is recommended that the customer read and review all sections of the manual for further knowledge and explanation of the system and its usage.

✓	Measurements and Reports
	Press Caliper to Activate the Measurement Function
	Positioning and Anchoring Measurement Markers
	Performing and Labeling a Measurement
	Editing a Measurement/Viewing Measurement Results
	Basic Measurement Tools/More/Shadow
	Calculating Ratios/Exam Specific Measurements
	Customizing Measurements: System Configuration > M & R > Measurement > General/Exam Package/Exam Specific/Per Label Configuration
	Customizing Measurements: System Configuration > M & R > Touch Screen Measurement Config > Package Configuration/Touch Screen Layout Editor
	Customizing Measurements: System Configuration > M & R > Custom Measurement Label > [SR Mapping]
	Customizing Measurements: System Configuration > M&R > Custom Calculations
	Customizing Measurements: System Configuration > M & R > Reports > Report Content/Assessments/Comments/Report Template
	Customizing Measurements: System Configuration > M & R > OB Tables > Default OB Author/Custom GA Tables/Custom Growth Curves

	Reports > Editing/Comments/Images/Clips and Protocol Views/Deleting
	Reports > Method of Measurement
	Store/Print/Transfer Report
	Summary Measurement Customization/ Drop-down Menu/Blank column option
✓	Modality Compare
	Compatibility (Mammography, CT, MRI, Ultrasound)
	Importing Data and/or Query Retrieve
	Activating Modality Compare (Under Workflow)
	Comparing Images and Clips
	Magnify/Adjust Position of the Reference Image
	Adjust Contrast and Brightness
✓	Panoramic Imaging
	Activate Panoramic Imaging (2D, 2D with Power, 2D with Color)
	Acquiring an Image/Erase/Review during Cine
	Rotate/Resize after Acquisition/Ruler
	Adjust Color ROI/Remove Color

# ACUSON Sequoia Ultrasound System

## Installation Checklist

### Basic Applications – General System 3.5 (VB30)

This checklist represents features discussed during system training by the Application Specialist/Clinical Sales Specialist. It is recommended that the customer read and review all sections of the manual for further knowledge and explanation of the system and its usage.

✓	Review and Patient Browser
	Review Current Study/Display Format/Thumbnails
	Pause/Discontinue/Active exam icons
	Single Select/Multiple Select
	Flag/Delete Images
	Reports/OB Tables
	Select All/Select None
	Compare/Page
	Exporting (with and without DICOM viewer)
	Exporting (with and without acoustic data and GSDF)
	Exporting Images from Thumbnails
	Query/Retrieve/Send to Network
✓	Printing and Recording
	Recording a Patient Study > System Configuration > Peripherals > DVR
	Printing Images/Report
	System Configuration > Peripheral Devices > Print/Store
	Automatic Disk Management > System Configuration > System Management > System Information
	Configuring a Report Printer
	Record Using DVR

✓	DICOM/Network Set-up
	DICOM Set-up
	DICOM Structured Reporting Set-up
✓	Documentation Devices
	Universal Video Converter
	Foot Switch, DVR, Print/Store > System Configuration > Peripheral Devices
	Printer (Configure under System Configuration/Peripheral Devices/Local Printer)
	DVR (Configure under System Configuration > Peripheral Devices > DVR – Enable DVR)
✓	System Management
	System Management
✓	System Information (SN/Version)
	System Information (SN/Version)
✓	Service
	Kinectus
	Back-up/Restore
	Sales/Service and Applications Support
✓	Security
	Security



I acknowledge that I have received training on items checked above on the ACUSON Sequoia 3.5 (VB30) ultrasound system and had the opportunity to ask questions.

Facility: \_\_\_\_\_

Date(s) of Training: \_\_\_\_\_

Customer Name: \_\_\_\_\_ Signature: \_\_\_\_\_

Clinical Specialist Name: \_\_\_\_\_ Signature: \_\_\_\_\_

---

ACUSON Sequoia, Clarify, TEQ, UltraArt, and Kinectus  
are trademarks of Siemens Medical Solutions, USA, Inc.

---

**Siemens Healthineers Headquarters**

Siemens Healthineers AG  
Siemensstr. 3  
91301 Forchheim, Germany  
Phone: +49 9191 18-0  
siemens-healthineers.com

**Manufacturer**

Siemens Medical Solutions USA, Inc.  
Ultrasound  
22010 S.E. 51st Street  
Issaquah, WA 98029, USA  
Phone: 1-888-826-9702  
siemens-healthineers.com/ultrasound