

CRDLIGHT **PROTECTIVE EQUIPMENT**

PRODUCT CATALOG

Sales Manager: **Wilson Overseas Manager**

Phone Number: **+86 139 2468 8922**

E-mail Address: **wilson@gdyanyang.com**

Jiangmen Yanyang Trading Co., Ltd.



Jiangmen Yanyang Trading Co., Ltd. was established in 2008. The company is headquartered in Jiangmen City, Guangdong Province. We are a modern hightech enterprise integrating sales, testing and service.

Committed to the benchmark of the Chinese nursing brand industry chain, in addition to a number of new technologies and perfect quality management systems leading the international industry standards, the company is also equipped with a technical management team and quality inspection team of 50 people, and has 10 automatic production lines, not only ensures the output of products, but also guarantees the high quality standards of products




营 业 执 照

(副 本) (副本号:1-1)

统一社会信用代码 91440704MA52768M32	扫描二维码登录“国家企业信用信息公示系统”了解更多登记、备案、许可、监管信息。
名 称 江门岩羊工贸有限公司	注 册 资 本 人民币伍拾万元
类 型 有限责任公司(自然人投资或控股)	成 立 日 期 2018年08月30日
法 定 代 表 人 胡杰	营 业 期 限 长期
经 营 范 围 研发、生产、加工、销售：LED、驱动电源、开关电源、电子整流器、智能化设备系列产品及配件、电器产品附件、低压电器产品、照明电器附件、电子产品、电气元器件、家用电器、五金、塑料制品、医疗器械、消毒产品（不含危险化学品）；医疗器械技术研发；生产、销售；机械设备、橡胶制品、日用品、口罩；货物或技术进出口（国家禁止或涉及行政审批的货物和技术进出口除外）。（依法须经批准的项目，经相关部门批准后方可开展经营活动。）	住 所 江门市江海区信义路18号2幢4层自编1号（信息申报制，一址多照）（经营场所：江门市江海区信义路18号5幢5层自编A04（信息申报制））（一照多址）
	登 记 机 关 
	2020 年 5 月 12 日

Production Factory

CRDLIGHT
Jiangmen Yanyang Trading Co., Ltd.



Daily output: 500,000PCS





**Gas Exchange
Pressure Tester**



Resistance Tester



**Ventilation
Resistance & Pressure
Difference Tester**



**Flame Retardant
Tester**



**Synthetic Blood
Penetration Tester**



**Surface Moisture
Resistance Tester**



**Particulate Filtration
Efficiency Tester**



PFE Tester



Electronic Tensile Tester

Filtering Half Mask


CRDLIGHT
Jiangmen Yanyang Trading Co., Ltd.

YY0525 FILTERING HALF MASK

Packing quantity: 25PCS/BOX; 1,000PCS/CARTON



CE 2163



EN 149:2001 + A1:2009
FFP2 NR
Standard



Breathable



Skin friendly



5 Layers

Filtering Half Mask

CRDLIGHT
Jiangmen Yanyang Trading Co., Ltd.

YY0525 FILTERING HALF MASK

Packing quantity: 25PCS/BOX; 1,000PCS/CARTON



YY0525

25pcs/color box; color box size: 12.5x9.5x12.5cm;
40boxes/carton; carton size: 51.5x49x26.5cm;
20GP: 352ctns;
40GP: 768ctns;
40HQ: 960ctns.

Filtering Half Mask

CRDLIGHT
Jiangmen Yanyang Trading Co., Ltd.

YY0525

FILTERING HALF MASK



UNIVERSAL



EU TYPE EXAMINATION CERTIFICATE

Certificate No: 2163-PPE-834

Respiratory protective devices, filtering half masks to protect against particles manufactured by
JIANGMEN YANYANG TRADING CO., LTD.
No.1, 4th Floor, Building 2, No. 18 Xinyi Road, Jianghai District, Jiangmen City, Guangdong Province, China
are tested and evaluated according to

**EN 149:2001 + A1:2009 Respiratory Protective Devices -
Filtering Half Masks to Protect Against Particles -
Requirements, Testing, Marking**

Based on the type examination conducted with the evaluation of test reports, technical file according to Personal Protective Equipment Regulation (EU) 2016/425 Annex 5, it is approved that the product meets the requirements of the regulation.

Product Definition

Brand Name: CRDLIGHT Model: YY0525
Filtering half mask
Classification: FFP2 NR

Here by the manufacturer is allowed to use notified body number (2163) and can fix CE mark, as shown below, on the Category III product models given above, with:

- Issuing an appropriate EU Declaration of Conformity according to Personal Protective Equipment Regulation (EU) 2016/425 Annex 9.
- Ongoing successful performance in fulfillment of the requirements set out in Personal Protective Equipment Regulation (EU) 2016/425 and harmonised standards, ensured by assessments based on Annex 7 (Module C2) or Annex 8 (Module D) of the regulation no later than 1 year from the beginning of serial production

This certificate is initially issued on 25/06/2020 and will be valid for 5 years, if there is no change in the relevant harmonised standard affecting the essential health and safety requirements.



Suat KACMAZ
Suat KACMAZ
UNIVERSAL CERTIFICATION
Director



id: Paraffin Oil Testing

No. of Sample	Paraffin Oil Testing 95 L/min max (%)	Requirements in accordance with EN 149:2001 + A1:2009	Result
-	3,0		Filtering half masks fulfill the requirements of the standard EN EN 149:2001 + A1:2009 given in 7.9.2 in range of the FFP1, FFP2 classes.
-	3,4	FFP1 ≤ 20 %	
-	3,3		
-	3,5	FFP2 ≤ 6 %	
-	3,5		
-	3,4	FFP3 ≤ 1 %	
-	4,7		
-	4,1		
-	4,5		

anical Strength
status Conditioning
ceived, original
ated wearing treatment

Practical Performance report, the likelihood of mask materials in contact with the skin causing irritation or other not reported.

Sample	Visual inspection	Requirements in accordance with EN 149:2001 + A1:2009	Result
	Burn for 1s	Filtering half mask shall not burn or not continue to burn for more than 5 s after removal from the flame	Passed
	Burn for 1s		
	Burn for 1s		
	Burn for 2s		

ceived, original
status Conditioning

be inhalation air:

CO ₂ content of the inhalation air [%] by volume	An average CO ₂ content of the inhalation air	Requirements in accordance with EN 149:2001 + A1:2009	Result
0,2	0,3 [%]	CO ₂ content of the inhalation air shall not exceed an average of 1,0% by volume	Passed
0,3			
0,3			

ceived, original

Performance and TTL test reports no adverse effects have been reported for donning and remove of the mask also the s that the ear loops / head harness are capable of holding the mask firmly enough.

Performance report, no adverse effects were reported for the field of vision availability when the mask is worn.

odel under inspection have no valves.

lation

e figures gathered for 9 different samples 3 as received, 3 with temperature conditioning and 3 simulated wearing tes with the limits given in the standard for FFP1, FFP2 and FFP3 classes. This is valid for inhalation results for 30 ize at 160 L/min.

Necip Fazıl Bulvarı Keyap Sitesi E2 Blok No:44/84 Yukarı Duðulu Ümranye - İSTANBUL - TÜRKİYE T: +90 216 455 80 80

UNIVERSALCERT.COM



Filtering Half Mask

CRDLIGHT
Jiangmen Yanyang Trading Co., Ltd.


YYC1028

COLORED FILTERING HALF MASK

Packing quantity: 25PCS/BOX; 1,000PCS/CARTON



CE 2163



EN 149:2001 + A1:2009

FFP2 NR

Standard



Breathable



Skin friendly



5 Layers

Filtering Half Mask

CRDLIGHT
Jiangmen Yanyang Trading Co., Ltd.

YYC1028

COLORED FILTERING HALF MASK

Packing quantity: 25PCS/BOX; 1,000PCS/CARTON



YYC1028

25pcs/color box; color box size: 12.5x9.5x12.5cm;

40boxes/carton; carton size: 51.5x49x26.5cm;

20GP: 352ctns;

40GP: 768ctns;

40HQ: 960ctns.

YYC1028

COLORED FILTERING HALF MASK



UNIVERSAL



EU TYPE EXAMINATION CERTIFICATE

Certificate No: 2163-PPE-1732

Respiratory protective devices, filtering half masks to protect against particles manufactured by
JIANGMEN YANYANG TRADING CO., LTD
No. 1, 4th Floor, Building 2, No.18 Xinyi Road, Jianghai District, Jiangmen City,
Guangdong Province, China
are tested and evaluated according to

**EN 149:2001 + A1:2009 Respiratory Protective Devices -
Filtering Half Masks to Protect Against Particles -
Requirements, Testing, Marking**

Based on the type examination conducted with the evaluation of test reports, technical file according to Personal Protective Equipment Regulation (EU) 2016/425 Annex 5, it is approved that the product meets the requirements of the regulation.

Product Definition

Single shift use particle filtering half mask for protection against solid and liquid aerosols, is a folding type, 5 layers cotton and polypropylene fabrics, without valve, fitted with ear loops, with internal nose clip and inside sponge strip.

Brand Name: CRDLIGHT **Model:** YYC1028 **Classification:** FFP2 NR
Model have white, black, grey, pink, and blue versions

Here by the manufacturer is allowed to use notified body number (2163) and can fix CE mark, as shown below, on the Category III product models given above, with;

- Issuing an appropriate EU Declaration of Conformity according to Personal Protective Equipment Regulation (EU) 2016/425 Annex 9.
- Ongoing successful performance in fulfilment of the requirements set out in Personal Protective Equipment Regulation (EU) 2016/425 and harmonised standards, ensured by assessments based on Annex 7 (Module C2) or Annex 8 (Module D) of the regulation no later than 1 year from the beginning of serial production

This certificate is initially issued on 04/12/2020 and will be valid for 5 years, if there is no change in the relevant harmonised standard affecting the essential health and safety requirements.



Suat KACMAZ
UNIVERSAL CERTIFICATION
Director

Necip Fazıl Bulvarı Keypay Sitesi E2 Blok No:44/84 Yukarı Dudullu Ümraniye - İSTANBUL - TURKEY T:+90 216 455 80 80

UNIVERSALCERT.COM



Table 1: Sodium Chloride Testing

No. of sample	Sodium Chloride Testing 95 L/min max (%)	Requirements in accordance with EN 149:2001 + A1:2009	Result
-	0,1	FFP1 ≤ 20 %	Filtering half masks fulfill the requirements of the standard EN EN 149:2001 + A1:2009 given in 7.9.2 in range of the first, and second protection classes. FFP1, FFP2, FFP3
-	0,2		
-	0,1		
-	0,2	FFP2 ≤ 6 %	
-	0,1	FFP3 ≤ 1 %	
-	0,2		
-	0,6		
-	0,5		
-	0,6		

Practical Strength 95 L/min = 1,6 dm³·an⁻¹

Table 2: Paraffin Oil Testing

No. of Sample	Paraffin Oil Testing 95 L/min max (%)	Requirements in accordance with EN 149:2001 + A1:2009	Result
-	0,2	FFP1 ≤ 20 %	Filtering half masks fulfill the requirements of the standard EN EN 149:2001 + A1:2009 given in 7.9.2 in range of the first and second protection classes. FFP1, FFP2, FFP3
-	0,3		
-	0,2		
-	0,3	FFP2 ≤ 6 %	
-	0,7	FFP3 ≤ 1 %	
-	0,3		
-	0,6		
-	0,8		

Practical Strength 95 L/min = 1,6 dm³·an⁻¹

Practical Performance report, the likelihood of mask materials in contact with the skin causing irritation or other not reported.

No. of sample	Visual inspection	Requirements in accordance with EN 149:2001 + A1:2009	Result
-	0,4s	Filtering half mask shall not burn or not continue to burn for more than 5 s after removal from the flame	Passed
-	0,4s		
-	0,4s		Filtering half masks fulfill requirements of the standard

Practical Performance report, the likelihood of mask materials in contact with the skin causing irritation or other not reported.

No. of sample	CO ₂ content of the inhalation air [%] by volume	As average CO ₂ content of the inhalation air	Requirements in accordance with EN 149:2001 + A1:2009	Result
-	0,6827	0,68	CO ₂ content of the inhalation air shall not exceed an average of 1,0% by volume	Passed
-	0,6813			
-	0,6810			

Practical Performance report, the likelihood of mask materials in contact with the skin causing irritation or other not reported.

Performance and TIL test reports no adverse effects have been reported for donning and remove of the mask also that the head harnesses are capable of holding the mask firmly enough.

7.14	FIELD OF VISION: in accordance with Performance report, no adverse effects were reported for the field of vision availability when the mask is worn.
Article 7.16	Breathing Resistance: Inhalation The overall evaluation in the figures gathered for 9 different samples 3 as received, 3 with temperature conditioning and 3 Simulated wearing conditioned, complies with the limits given in the standard for FFP1, FFP2 and FFP3 classes. This is valid for inhalation results for 30 L/min, 95 L/min and exhalation at 160 L/min. Passed.



Filtering Half Mask

CRDLIGHT
Jiangmen Yanyang Trading Co., Ltd.

YY1029

FFP3 FILTERING HALF MASK

Packing quantity: 25PCS/BOX; 1,000PCS/CARTON




FFP3

99% filtration efficiency



CE 2163



EN 149:2001 + A1:2009

FFP3 NR

Standard



Breathable



Skin friendly



5 Layers

Filtering Half Mask

CRDLIGHT
Jiangmen Yanyang Trading Co., Ltd.

YY1029

FFP3 FILTERING HALF MASK

Packing quantity: 25PCS/BOX; 1,000PCS/CARTON



YY1029

25pcs/color box; color box size: 16x11x12.5cm;
40boxes/carton; carton size: 56.5x51.5x33.5cm;
20GP: 240ctns;
40GP: 480ctns;
40HQ: 720ctns.

Glove Introductions

CRDLIGHT
Jiangmen Yanyang Trading Co., Ltd.

COLOR: CLEAR / BLACK / BLUE



SIZES: S / M / L / XL

Suggestion S Size

Palm of your hand is 7.1-8cm wide

Suggestion M Size

The palm of your hand is 8.1-9cm wide


Suggestion L Size

The palm of your hand is 9.1-10cm wide

Suggestion XL Size

The palm of your hand is 10.1-11cm wide

Most women use size S and those with larger hands use size M
Most men use size M and those with larger hands use size L



Manufacturing

CRDLIGHT
Jiangmen Yanyang Trading Co., Ltd.



Vinyl Gloves

CRDLIGHT
Jiangmen Yanyang Trading Co., Ltd.

YYPC01

PVC EXAMINATION GLOVE (Clear)

YYPC02

PVC EXAMINATION GLOVE (Blue)

SIZE: S / M / L / XL

- NON-STERILE
- POWDER-FREE
- AMBIDEXTROUS
- SINGLE USE
- LATEX-FREE

YYPC01/02

100pcs/color box (23x12.5x6.5cm)
10boxes/carton (35x27x25cm)
20GP: 1,024ctns;
40GP: 2,112ctns;
40HQ: 2,640ctns.



Nitrile Gloves

CRDLIGHT
Jiangmen Yanyang Trading Co., Ltd.

YYNG01

NITRILE EXAMINATION GLOVE (Black)



YYNG02

NITRILE EXAMINATION GLOVE (Blue)

SIZE: S / M / L / XL



- NON-STERILE
- POWDER-FREE
- AMBIDEXTROUS
- SINGLE USE
- LATEX-FREE

YYNG01/02

100pcs/color box (23x12.5x6.5cm)
10boxes/carton (35x27x25cm)
20GP: 1,024ctns;
40GP: 2,112ctns;
40HQ: 2,640ctns.



Wide usage
application

Blended Gloves

CRDLIGHT
Jiangmen Yanyang Trading Co., Ltd.

YYBG01

PVC/NITRILE BLENDED GLOVE (Black)



YYBG02

PVC/NITRILE BLENDED GLOVE (Black)

SIZE: S / M / L / XL



- NON-STERILE
- POWDER-FREE
- AMBIDEXTROUS
- SINGLE USE
- LATEX-FREE

YYBG01/02

100pcs/color box (23x12.5x6.5cm)
10boxes/carton (35x27x25cm)
20GP: 1,024ctns;
40GP: 2,112ctns;
40HQ: 2,640ctns.

