

SURGICAL CAP ANNIE



Ref. no	Description	RET/TRP
621500	Surgical Cap L (elastic), mixed colours	150/600
621715	Surgical Cap M (elastic), blue	150/600
621725	Surgical Cap M (elastic), green	150/600
621735	Surgical Cap M (elastic), white	150/600
621815	Surgical Cap L (elastic), blue	150/600
621825	Surgical Cap L (elastic), green	150/600
621835	Surgical Cap L (elastic), white	150/600
621925	Surgical Cap XL (elastic), green	150/600

SURGICAL CAP CHARLOTTE



Ref. no	Description	RET/TRP
620900	Surgical Cap M (elastic), white	-/200

SURGICAL CAP EUROCAP



Ref. no	Description	RET/TRP
620800	Surgical Cap one-size (elastic), blue	100/500
620810	Surgical Cap one-size (elastic), green	100/500

SURGICAL CAP KOSACK



Ref. no	Description	RET/TRP
621000	Surgical Cap one-size (elastic), green	100/500
621001	Surgical Cap one-size (elastic), purple	100/500
621005	Surgical Cap one-size (elastic), blue	100/500

SURGICAL CAP JACK



Ref. no	Description	RET/TRP
620300	Surgical Cap one-size (elastic), green	100/500

SURGICAL CAP PHILIP



Ref. no	Description	RET/TRP
621301	Surgical Cap one-size (tie-band), green	100/500

SURGICAL CAP FLORY



Ref. no	Description	RET/TRP
620500	Surgical Cap one-size (tie-band), green	100/500



Tray ID	97095457
Tray name	SET CHIRURGICAL ABORD CEREBRAL
Colour code	
Speciality filter	NEUROSURGERY
Intervention filter	Craniotomy
Hospital filter	Tehnomedica SRL

Description		Qty
	Crepe paper 60x60cm 60g/m2 White	1
	BNS Craniotomy Drape 230x300cm Inc. 30x20	1
	Adh. op towel 50x50cm	4
	Hand towel 18x25cm	2
	Op tape 9x49cm	1
	Suction And Diathermy Bag, 40x35 cm	1
	Neuro swab gauze 12,7x38,1mm 6P XRD dot String White 10/1	10



Microscope Drape Zeiss 117 x 267cm

1



Forceps Bipolar Bayonet

1



Medicine/Specimen cup 120ml Polyethylene Transp.
Lid

1



Scalpel blade CS N°22

1



Plastic bag 20x30cm Ziplock Transp.

1



Instrument cleaner 5x5cm

1



Scalpel blade CS N°15

1



Syringe 20ml L/L 3parts Conc.

3



Light handle cover large Light green/white collar

2



Suction Catheter 12Ch 50cm 2eyes Funnel str PVC
Shore 70

1



Mayo stand cover 79x145cm Abs.65x85cm

2



Reinforced Table Cover 190 x280 wrapping

1




SA Label 1,6 x 5cm Red TEHNOMEDICA SRL

1



Tray ID	97094658
Tray name	SET CHIRURGICAL ABORD SPINAL
Colour code	
Speciality filter	NEUROSURGERY
Intervention filter	Laminectomy/Spinal
Hospital filter	Tehnomedica SRL

Description		Qty
	Adh. op sheet 240x150cm	1
	Op tape 9x49cm	1
	Adh. op sheet 175x175cm	1
	Incise drape 30x26cm Inc.30x20cm	1
	Suction And Diathermy Bag, 40x35 cm	2
	Crepe paper 60x60cm 60g/m2 White	1
	Adhesive op towel 90x75cm Reinf.	2



Hand towel 18x25cm

4



Neuro swab gauze 12,7x38,1mm 6P XRD
dot String White

10



Medicine/Specimen cup 120ml
Polyethylene Transp. Lid

1



Banded bag 140cm Circ. Elast. Transp.

1



Microscope Drape Zeiss 117 x 267cm

1



Scalpel blade CS N°11

1



Forceps Bipolar Straight 200mm

1



Scalpel blade CS N°22

1



Plastic bag 20x30cm Ziplock Transp.

1



Instrument cleaner 5x5cm

1



Syringe 20ml L/L 3parts Conc.

3



Light handle cover large Light green/white collar

2



Suction Catheter 12Ch 50cm 2eyes
Funnel str PVC Shore 70

1



Mayo stand cover 79x145cm
Abs.65x85cm

3



Reinforced Table Cover 190 x280
wrapping

1



SA Label 1,6 x 5cm Red TEHNOMEDICA
SRL

1

706100

Adhesive Absorbent Pad

Product details

Size: 55cm x 70cm

Descriptive feature: Adhesive, full length

Sterility: Sterile

Images



Delivered items

706100-09

Sales released in: Australia, Austria, Azerbaijan, Bahrain, Belgium, Czechia, Denmark, Finland, Germany, Hong Kong, Hungary, Iceland, Ireland, Italy, Luxembourg, Macedonia (the former Yugoslav Republic of), Moldova (the Republic of), Netherlands, New Zealand, Norway, Portugal, Russian Federation, Slovenia, Spain, Sweden, Switzerland, United Arab Emirates, United Kingdom of Great Britain and Northern Ireland

Country of origin: Belgium

Shelf life: 5 years

Sterilization method: Beta

Packing information: First packaging layer is a peel-open sterile barrier, plastic/plastic. Once opened the sterile barrier cannot be closed again. Second layer is a corrugated board dispenser box. Third layer is a corrugated board transport box.

Packaging level	Pack count	GS1 code
Consumer pack	1	7323190020065
Dispenser box	40	7323190019977

Packaging level	Pack count	GS1 code
Transport box	80	7323190019960
Pallet	1920	7323190019953

Material

Animal tissues:

No

Human blood derivatives:

No

Medicinal substances:

No

Phthalates:

No

Polyvinyl chloride:

No

Product Composition Drapes

Product Component	Critical Area	Less Critical Area
Drape material	Viscose/Polyester nonwoven 80 g/m ²	Viscose/Polyester nonwoven 80 g/m ²
Drape material	Polyethylene film 40 µm	Polyethylene film 40 µm
Adhesive material	Synthetic rubber based	N/A

Product Performance Drapes, Additional Tests

Characteristics	Test Method	Internal Test Method	Unit	Product Performance Critical Product Area	Product Performance Less Critical Product Area
Absorption	ISO 9073-12	T-1158	g/dm ²	3.66	3.66
Flammability	16 CFR 1610.4	N/A	s	Class 1, >3.5s	Class 1, >3.5s

Product Performance Sterile Drapes, EN 13795 High Performance

Characteristic s	Test Metho d	Interna l Test Metho d	Unit	Requireme nt Critical Product Area	Requireme nt Less Critical Product Area	Product Performanc e Critical Product Area	Product Performanc e Less Critical Product Area
Resistance to microbial penetration - Dry	ISO 22612	T-1004	CFU	Not required	≤300	-	0
Resistance to microbial penetration - Wet	ISO 22610	T-1005	BI	6.0	Not required	6.0	-
Cleanliness - Particulate Matter	EN ISO 9073-10	T-1006	IPM	≤3.5	≤3.5	1.5	1.5
Linting	EN ISO 9073-10	T-1006	Log ₁₀ (lint count)	≤4.0	≤4.0	1.8	1.8
Resistance to liquid penetration	EN 20811	T-280	cm H ₂ O	≥100	≥10	>100	>100
Bursting strength - Dry	EN ISO 13938-1	T-233 or T-1179	kPa	≥40	≥40	168	168
Bursting strength - Wet	EN ISO 13938-1	T-233 or T-1179	kPa	≥40	Not required	144	-
Tensile strength - Dry	EN 29073-3	T-229	N	≥20	≥20	36	36
Tensile strength - Wet	EN 29073-3	T-229	N	≥20	Not required	28	-

Technical

Dimension

Dimension text	Dimension value
Outer dimension	55 cm x 70 cm

Disposal instructions

Non-hazardous waste used BARRIER products and sterility barriers should, in the majority of cases, be classified as non-hazardous waste. They contain high amounts of energy and are well suited for incineration. BARRIER products do not contain any hazardous substances that can leach out if the products are land filled. Transport boxes are designed to fit existing recovery systems. The new BARRIER packaging system complies with the Packaging Waste Directive of the European Union.

Storage instructions

Mölnlycke Health Care recommends that BARRIER products are stored under normal storage conditions. All layers of packaging should be kept intact until access to the underlying layers is needed. Storage facilities for products only protected by the sterility barrier should be kept under conditions where low level of particulate air contamination prevail, so that it would not constitute a risk to the patient when the package is opened and the product is used.

Classification

Regulation type	MDD Class IS	Locally Regulated	Unregulated
MDD Classification Rule:	1		
CE Certificate Number:	CE 01966		
Notified body medical devices/PPE:	BSI (0086)		
Intended use MDD:	Surgical drapes, when sterilised, are intended to minimize the spread of micro-organisms, in order to reduce the risk for post operative wound infection.		
Sales released in:	Austria, Belgium, Czechia, Denmark, Finland, Germany, Hungary, Iceland, Ireland, Italy, Luxembourg, Netherlands, Norway, Portugal, Slovenia, Spain, Sweden, Switzerland, United Kingdom of Great Britain and Northern Ireland	Australia, Azerbaijan, Bahrain, Moldova (the Republic of), New Zealand, Russian Federation, United Arab Emirates	Hong Kong, Macedonia (the former Yugoslav Republic of)

Applied standards: The standards presented below is a selection of the most essential standards that are adhered to.

EN 1041, EN 556-1, EN 13795, EN 62366, EN ISO 9001, EN ISO 13485, EN ISO 10993-1, EN ISO 11607-1, EN ISO 11607-2, EN ISO 15223-1, ISO 15223-2, ISO 14001

Removable Label

Yes

GMDN Code (Global Medical Device Nomenclature)

47783