

Order Information

| Disposable Blood Line | | | |
|-----------------------|-------------------------------|---------------|--|
| Product Code | Pump Segment (IDxODxLength)mm | Chamber OD mm | Descriptions |
| FB-A101V11 | 6.4x9.75x350 | 22 | B.Braun(HD-SECURA) |
| FB-A001V05 | 7.9x12.1x350 | 22 | Gambro (except AK95), Fresenius (Venous Chamber 22mm) Nikkiso, Baxter, Nipro Toray, JMS. B. Braun. Bellco. |
| FB-A001V02 | 7.9x12.1x350 | 30 | Fresenius (Venous Chamber 30mm) |
| FB-A001V05 | 7.9x12.1x350 | 20 | Gambro AK95, Baxter, Nikkiso, Nipro Toray, JMS, B. Braun, Bellco. |

| Disposable A.V.Fistula Needle | | | | | |
|-------------------------------|-----------------|-------|-----------|-------------------|-----------------|
| Code(Single Pack) | Code(Twin Pack) | Gauge | Wing | Needle Length(mm) | Tube Length(mm) |
| FN-1511S | FN-1511D | 15G | Fixed | 25±2.0mm | 150±2.0mm |
| FN-1611S | FN-1611D | 16G | Fixed | 25±2.0mm | 150±2.0mm |
| FN-1711S | FN-1711D | 17G | Fixed | 25±2.0mm | 150±2.0mm |
| FN-1512S | FN-1512D | 15G | Fixed | 25±2.0mm | 300±2.0mm |
| FN-1612S | FN-1612D | 16G | Fixed | 25±2.0mm | 300±2.0mm |
| FN-1712S | FN-1712D | 17G | Fixed | 25±2.0mm | 300±2.0mm |
| FN-1621S | FN-1621D | 16G | Fixed | 32±2.0mm | 150±2.0mm |
| FN-1721S | FN-1721D | 17G | Fixed | 32±2.0mm | 150±2.0mm |
| FN-1622S | FN-1622D | 16G | Fixed | 32±2.0mm | 300±2.0mm |
| FN-1722S | FN-1722D | 17G | Rotatable | 32±2.0mm | 300±2.0mm |
| FN-1511ZS | FN-1511ZD | 15G | Rotatable | 25±2.0mm | 150±2.0mm |
| FN-1611ZS | FN-1611ZD | 16G | Rotatable | 25±2.0mm | 150±2.0mm |
| FN-1711ZS | FN-1711ZD | 17G | Rotatable | 25±2.0mm | 150±2.0mm |
| FN-1512ZS | FN-1512ZD | 15G | Rotatable | 25±2.0mm | 300±2.0mm |
| FN-1612ZS | FN-1612ZD | 16G | Rotatable | 25±2.0mm | 300±2.0mm |
| FN-1712ZS | FN-1712ZD | 17G | Rotatable | 25±2.0mm | 300±2.0mm |
| FN-1621ZS | FN-1621ZD | 16G | Rotatable | 32±2.0mm | 150±2.0mm |
| FN-1721ZS | FN-1721ZD | 17G | Rotatable | 32±2.0mm | 150±2.0mm |
| FN-1622ZS | FN-1622ZD | 16G | Rotatable | 32±2.0mm | 300±2.0mm |
| FN-1722ZS | FN-1722ZD | 17G | Rotatable | 32±2.0mm | 300±2.0mm |

ABLE[®] HAEMODIALYSIS SERIES PRODUCTS



GUANGDONG BAIHE MEDICAL TECHNOLOGY CO., LTD.
 No. 89, Taoyuan East Road, Nanhai, Foshan, Guangdong Province, P.R.C.
 Tel: +86-757-81207300 +86-757-81207302
 Fax: +86-757-81207311
 Email: able@baihemedical.com sales@baihemedical.com
 Website: www.baihemedical.com




Version: 202103

GUANGDONG BAIHE MEDICAL TECHNOLOGY CO., LTD.

BLOOD LINE

Application:

- This product is used to establish an extracorporeal circulation channel for blood dialysis therapy.




Features:


- 01 Pump Tube**
 - With high elasticity and medical grade PVC, the shape of tube remain the same after continuous pressing of 10 hours.
- 02 Drip Chamber**
 - Available various diameter sizes.
- 03 Dialysis Connector**
 - Extra large designed dialyzer connector is easy to operate.
- 04 Clamp**
 - Clamp is made of hard plastic and designed larger and thicker to guarantee sufficient stop.
- 05 Infusion Set**
 - It is convenient to install and uninstall, which ensure precision infusion and safe priming.
- 06 Drainage Bag**
 - Closed priming to meet the requirements of quality control, singleway drainage bag and double way drainage bag available.
- 07 Customized Designed**
 - Can be customized with all orders.

SUBSTITUTION LINE

L02



Y02




A.V. FISTULA NEEDLE



Features:

- Fine polishing process on blade to puncture easily and smoothly.
- Extra thin wall with back eye to get larger blood flow.
- Rotatable wing and fixed wing available.



Dull Needle

Sharp Needle