

ACCESSORIES AND SPARE PARTS CATALOGUE

CAT-ACS-2020-03-24/EX












TABLE OF CONTENTS:






Gas cylinder connections for CRYO-S/CRYO-T devices	p. 1
Accessories for CRYO-S/CRYO-T devices	pp. 2-9
Accessories for cryoprobes/cryoneedles	pp. 10-11
Accessories for ligator LC	p. 12
Accessories for CRYO-T Elephant devices	pp. 13-15
Accessories for SALSA and SWING lasers	pp. 16-17
Accessories for ATO-3 device	pp. 18-20
Accessories for BOA devices	pp. 21-23
Accessories for CRYO-T Cooler device	p. 24
Electrical power cords and fuses	p. 25

EXPLANATION OF ABBREVIATIONS:



CSM	CRYO-S Mini
CSC	CRYO-S Classic
CSEA	CRYO-S Electric II
CSEP	CRYO-S Painless
AK	previous cryosurgery device
CT2	CRYO-T 2
CTD	CRYO-T Duo
Cryo-T/NR-2	previous cryotherapy devices
CTC	CRYO-T Cooler
AOM	ATO-3
CTEH	CRYO-T Elephant
CTEHG	CRYO-T Elephant G
CTEHm	CRYO-T Elephant Mini M
CTEHgm	CRYO-T Elephant Mini GM
CTEHgs	CRYO-T Elephant Mini GS
PTC	BOA Mini+ (5 and 10 outlets version)
PTX2	BOA Max2
SALSA	all SALSA series lasers
SWING	all SWING series lasers
TWIST	TWIST Diode laser
CO2	all CO2 lasers offered by METRUM CRYOFLEX

REF PN01	N2O connection Poland		
PA01N300006	application: CSEP, CSEA, CSC, CSM	G3/8" - 19 TPI	
REF WC01	C02 connection Kazakhstan		
003835	application: CSEP, CSEA, CSC, CSM	G3/4"	
REF WC02	Medical C02 connection Italy		
PA01N300007	application: CSEP, CSEA, CSC, CSM, CTD, CT2	W 27x2	
REF WN01	N2O connection UK, ZEA		
PA01N300008	application: CSEP, CSEA, CSC, CSM	W 11/16" - 20TPI	
REF PC01	Connection: - CO2 Poland - N2O Lithuania, Belarus		
PA02N000001	application: CSEP, CSEA, CSC, CSM	W 21.8 - 14TPI	
REF WN02	N2O connection Mexico		
004186	application: CSEP, CSEA, CSC, CSM	G1/2" - 14TPI	
REF WC03	C02 connection with filter Korea		
002195	application: CT2	W 21.8 - 14TPI	
REF PC02	C02 connection with filter Poland, Mexico		
006003	application: CT2	W 21.8 - 14TPI	
REF WC04	C02 pin index connection		WARNING! Engraving "C02" on the connection. Pin spacing ≈ 14mm
004224	application: CSEP, CSEA, CSC, CSM		
REF WN03	N2O pin index connection		WARNING! Engraving "N2O" on the connection. Pin spacing smaller than in WC04
004225	application: CSEP, CSEA, CSC, CSM		

REF PZP05	Pneumatic hose gas cylinder-device	
PA02N200006	length: 0.5m application: CSEP, CSEA, CSC, CSM, CTD, CT2	
REF PZP10	Pneumatic hose gas cylinder-device	
PA02N200007	length: 1m application: CSEP, CSEA, CSC, CSM, CTD, CT2	
REF PZP20	Pneumatic hose gas cylinder-device	
PA02N200008	length: 2m application: CSEP, CSEA, CSC, CSM, CTD, CT2	
REF PZP30	Pneumatic hose gas cylinder-device	
PA02N200009	length: 3m application: CSEP, CSEA, CSC, CSM, CTD, CT2	
REF ZCZ.8	Pneumatic hose plug	
000774	size: 8mm application: CSEP, CSEA, CSC, CSM, CTD, CT2	
REF UOPZP.A	O-ring for pneumatic hose	
000599	size: 8.1mm application: CSEP, CSEA, CSC, CSM, CTD, CT2	
REF PZS.B	Hose for cryoprobes	
PA01N300003	length: 2m application: CSEA	
REF PZS.N	Hose for cryoprobes with neurostimulation	
PA01N300004	length: 2m application: CSEP	
REF 郑PS.A	Cryoprobe hose plug (device side)	
002088	application: CSEP, CSEA	
REF 郑PS.B	Cryoprobe hose plug (cryoprobe side)	
003284	application: CSEP, CSEA	

REF UOPZS.A	O-ring A for cryoprobe hose and SJA	
001821	size: Ø13.1mm application: CSEP, CSEA	
REF UOPZS.B	O-ring B for cryoprobe hose	
000566	size: Ø10.1mm application: CSEP, CSEA, CSC, CSM	
REF UOPZS.C	O-ring C for cryoprobe hose	
000558	size: Ø6mm application: CSEA, CSEP, CSC, CSM	
REF PZS.A	Hose for cryoprobes	
PA01N300002	length: 2m application: CSC, CSM, AK	
REF ZCZ.10	Cryoprobe hose plug	
001392	size: 10mm application: CSC, CSM	
REF PZD.A	Cable for nozzle	
PA02N200001	length: 2m application: CT2, CTD	
REF PZD.B	Cable for nozzle	
SERWIS	length: 2m application: CRYO-T, NR-2	
REF APZ.A	Tee for pneumatic hose gas cylinder-device	
- - - - -	application: CSEP, CSEA, CSC, CSM, CTD, CT2	

REF ACSS.A	Cone for o-rings size: LARGE ø14mm application: CSEP, CSEA	~ Ø14mm 
003840		
REF ACSS.B	Cone for o-rings size: MEDIUM ø8mm application: CSEP, CSEA, CSC, CSM, CTD, CT2	~ Ø8mm 
000480		
REF ACSS.C	Cone for o-rings size: SMALL ø6mm application: CSEP, CSEA, CSC, CSM	~ Ø6mm 
000479		
REF ACSS.D	Tool for o-rings application: CSEP, CSEA, CSC, CSM, CTD, CT2	
000538		
REF ACSK.A	Flat wrench for gas cylinder connection size: 22mm application: CSC, CSM	
000629		
REF ACSK.B	Flat wrench for gas cylinder connection size: 24 / 27mm application: CSC, CSM	
000630		
REF ACSK.C	Adjustable wrench for gas cylinder connection size: max. 35mm application: CSEA, CSEP, CSC, CSM, CTD, CT2	
004606		
REF ZCZ.12	Plug for device front/back sockets size: 12mm application: CSEP, CSEA, CSC, CSM, CTD, CT2	~ Ø12mm 
000775		
REF ZCZ.24	Plug for device front socket size: 24mm application: CSEA, CSEP	~ Ø24mm 
003279		

REF PWN.E	Two-pedal footswitch	 connection XLR 3 PIN
SERWIS	cable length: 2.4m application: CSEP, CSEA	
REF ACSE.E	Pantograph for cryoprobes	
SERWIS	size: L=50cm application: CSEP, CSEA	
REF ACSE.F	Passive electrode cable	
004813	application: CSEP	
REF ACSS.G	Passive electrode	
003774	application: CSEP	
REF SJA	Single use probe adapter (set)	 -1x adapter -1x flat wrench #10mm (ASJA.A) -1x Allen wrench #0.9mm (ASJA.B) -3x o-ring 1.42x1.52 -3x o-ring 1.78x1.78 -9x o-ring 8.3x2.4 (UOPZS.A) -6x o-ring 6.0x1.5 -1x Cone for o-rings large (ACSS.A) -2x setscrew M2x5 (ASJA.E) -1x instruction of o-ring replacement
PA01N270004	application: CSEP	
REF ASJA.A	Flat wrench for SJA	
003307	size: 10mm application: SJA	
REF ASJA.B	Allen wrench for SJA	
003080	size: 0.9mm application: SJA	
REF ASJA.C	Set of o-rings red	 Set includes: - 2x 1.42x1.52mm - 2x 6.0x1.5mm
002327 003300	size: 1.42x1.52mm and 6.0x1.5mm application: SJA  (red)	
REF ASJA.D	Set of o-rings blue	 Set includes: - 2x 1.78x1.78mm - 2x 6.0x1.5mm
002970 003300	size: 1.78x1.78mm and 6.0x1.5mm application: SJA  (blue)	
REF ASJA.E	Setscrew	
004295	size: M2x5mm application: SJA	

REF ACSM.A	Device fixing screw	M6x20	
001684	size: M6x20 application: CSC, CTD		
REF ACSM.B	Device fixing screw	4mm	
000433	size: M4x25 application: CSM		
REF ACSM.C	Allen wrench for fixing screw	M4x25	
000627	size: 4mm application: CSC, CTD		
REF AZN.A	Set of accessories for Cryo-S Classic	- 1x N2O connection Poland (PN01) - 1x CO2 connection Poland (PC01) - 1x flat gasket for N2O connection (UPZB.A) - 1x flat gasket for CO2 connection (UPZB.B) - 1x Allen wrench for fixing screw (ACSM.C) - 5x o-ring 3.3x2.4mm (UOPZP.A) - 3x o-ring 2.57x1.78mm (UOPZS.C) - 2x o-ring 5.3x2.4mm (UOPZS.B) - 1x cone for o-rings medium (ACSS.B) - 1x cone for o-rings small (ACSS.C) - 1x tool for o-rings (ACSS.D) - 1x instruction of o-ring replacement	
004511	application: CSC		
REF AZN.B	Set of accessories for Cryo-S Mini	- 1x N2O connection Poland (PN01) - 1x CO2 connection Poland (PC01) - 1x flat gasket for N2O connection (UPZB.A) - 1x flat gasket for CO2 connection (UPZB.B) - 5x o-ring 3.3x2.4mm (UOPZP.A) - 3x o-ring 2.57x1.78mm (UOPZS.C) - 2x o-ring 5.3x2.4mm (UOPZS.B) - 1x cone for o-rings medium (ACSS.B) - 1x cone for o-rings small (ACSS.C) - 1x tool for o-rings (ACSS.D) - 1x instruction of o-ring replacement	
004513	application: CSM		
REF AZN.C	Set of accessories for Cryo-S Electric II/ Painless	- 1x N2O connection Poland (PN01) - 1x CO2 connection Poland (PC01) - 1x flat gasket for N2O connection (UPZB.A) - 1x flat gasket for CO2 connection (UPZB.B) - 5x o-ring 8.3x2.4mm (UOPZS.A) - 2x o-ring 5.3x2.4mm (UOPZS.B) - 3x o-ring 3.3x2.4mm (UOPZP.A) - 3x o-ring 2.57x1.78mm (UOPZS.C) - 1x cone for o-rings medium (ACSS.B) - 1x cone for o-rings small (ACSS.C) - 1x tool for o-rings (ACSS.D) - 1x instruction of o-ring replacement - 2x fuse 1.25A (B520.C) - 1x plug for device socket (ZCZ.24) - 1x plug for device socket (ZCZ.12)	
PA01N300005	application: CSEA, CSEP		



REF UPZB.A	Flat gasket for N2O connection	Ø13.7 x Ø7mm	
003679	size: 13.7x7 #3,5mm application: PN01, WN01		
REF UPZB.B	Flat gasket for CO2 connection	Ø19.5 x Ø11mm	
003570	size: 19.5x11 #4mm application: PC01, WC01		
REF UPZB.C	Gasket for CO2 connection with filter	Ø17 x Ø10mm	
003720	size: 17x10 #7mm application: PC02		
REF UOFCO.A	Black o-ring for CO2 connection with filter		
000557	size: 18x3 application: PC02, WC03		
REF UOFCO.B	Red o-ring for CO2 connection with filter		
000555	size: 18x3 application: PC02, WC03		
REF FCO.A	Paper filter for CO2 connection with filter		
RA-201	application: PC02, WC03		
REF FCO.B	Sintered filter for CO2 connection with filter		
000020	application: PC02, WC03		
REF ACT.Z	Power supply 24V DC		
005500	application: CT2		
REF	Plexi shelf for trolley		
	application: trolley for CSEP and CSEA		


REF AZN.D	Set of accessories for Cryo-T 2 / Cryo-T Duo application: CT2, CTD	- 1x CO2 connection with filter Poland (PC02) - 1x gasket 17x10mm, #7mm (UPZB.C) - 1x paper filter for CO2 connection (FCO.A) - 1x sintered filter for CO2 connection (FCO.B) - 1x black o-ring 18x3mm (UOFCO.A) - 1x red o-ring 18x3mm (UOFCO.B) - 5x o-ring 3.3x2.4 (UOPZP.A) - 1x cone for o-rings medium (ACSS.B) - 1x tool for o-rings (ACSS.D) - 1x instruction of o-ring replacement - 1x instruction of CO2 connection assembling - 2x fuse 600mA (B520.A) - 1x Allen wrench for fixing screw (ACSM.C) - CTD only
PA02N000000		
REF AZN.E	Set of accessories for Cryo-T 2 Korea application: CT2 (Korean version)	- 1x CO2 connection with filter Korea (WC03) - 1x paper filter for CO2 connection (FCO.A) - 1x sintered filter for CO2 connection (FCO.B) - 1x black o-ring 18x3mm (UOFCO.A) - 1x red o-ring 18x3mm (UOFCO.B) - 5x o-ring 5.3x2.4mm (UOPZS.B) - 5x o-ring 3.3x2.4mm (UOPZP.A) - 1x cone for o-rings medium (ACSS.B) - 1x tool for o-rings (ACSS.D) - 1x instruction of o-ring replacement - 1x instruction of korean CO2 connection assembling - 2x fuse 600mA (B520.A) - 1x wrench for connection with filter (FCO.C) - 1x angular connector (ACT.B)
004510		
REF FCO.C	Wrench for connection with filter application: WC03	
003227		
REF ACT.A	Flat wrench for gas cylinder connection size: 32mm application: CT2, CTD	
000632		
REF ACT.B	Angular connector Korea size: 20mm application: CT2 KOREA	
SERWIS		
REF PWN.M	Footswitch with pneumatic valve cable length: 1.8m application: CSM	single stage connection Ø6mm 
005087		
REF PWN.C	Footswitch with pneumatic valve cable length: 1.8m application: CSC	double stage connection Ø10/Ø6mm 
005088		
REF ZCZ.6	Footswitch hose plug size: 6mm application: CSM	~Ø6mm 
001791		

REF DT1S.A	Cryotherapy nozzle D1 application: CTD, CT2	
PA02N100001		
REF DT2S.A	Cryotherapy nozzle D2 application: CTD, CT2	
PA02N100002		
REF DT3S.A	Cryotherapy nozzle D3 application: CTD, CT2	
PA02N100003		
REF DT1K.A	Cryotherapy nozzle D1 with switch application: CT2	

REF DT1K.H	Cryotherapy hyperbaric nozzle D4 with switch application: CT2	
PA02N100007		
REF DT2K.A	Cryotherapy nozzle D2 with switch application: CT2	
PA02N100005		
REF DT3K.A	Cryotherapy nozzle D3 with switch application: CT2	
PA02N100006		

REF ASSD.A	Set of accessories for cryoprobe SD-10 application: CSC, CSM	 - 1x tool for spray type probe cleaning - 6x freezing field limiter - 10x swab for large lesions freezing - 1x tool for swab preparation
PA01N210001		
REF ASSG.A	Set of accessories for cryoprobe SG-10 application: CSC, CSM	 - 1x tool for spray type probe cleaning - 6x freezing field limiter - 10x swab for large lesions freezing - 1x tool for swab preparation
PA01N230001		
REF ASSD.B	Plug for SD-10 cryoprobe size: Ø1.6mm, L=19mm application: SD-10	~ Ø1.6x19mm 
002053		
REF ASSG.B	Plug for SG-10 cryoprobe size: Ø1.8mm, L=25.5mm application: SG-10	~ Ø1.8x25.5mm 
001607		
REF ASSD.C	Tip cover for SD-10 cryoprobe size: Ø8mm/10mm, L=12mm application: SD-10	Ø8/ Ø10mm L=12mm 
000824		
REF ASK.A	Plug for cryoprobe application: cryoprobes for CSC and CSM	M16x1 
000429		
REF ASK.B	Plug for cryoprobe with RFID application: cryoprobes for CSEA and CSEP	M16x1 
PA01N270003		
REF ASK.C	Cannula for cryoprobe 2mm, long application: cryoprobes for CSEP	 
TA01S100006		




REF ASK.D TA01N100003	Sterilization box large size: 290x80x55mm application: cryoprobes for CSEP, CSEA, CSC, CSM	L=290mm D=80mm H=55mm	
REF ASK.E TA01N100004	Sterilization box small size: 180x80x51mm application: cryoprobes for CSEP, CSEA, CSC, CSM	L=180mm D=80mm H=51mm	

REF ALP.Z	Proctological ligator LC application: CSC	
PA01N100001		
REF ALP.A	Cone for ligator size: ø14mm, L=44mm application: LC	~ Ø 14 L=44mm 
004075		
REF ALP.B	Allen key for ligator size: 2mm application: LC	#2mm 

REF ALP.C	Tool for ligation rings application: LC	
001089		
REF ALP.D	O-ring for tool for ligation rings size: 8x2mm application: LC	8.0x2.0 NBR 70Sh 
000573		
REF ALP.E	Silicone hose size: 3.5x6mm application: LC	
000206		
REF ALP.F	HEPA filter application: LC	

REF ALP.G	Luer connection CT 61.1 application: LC	
000982		
REF ALP.H	Luer connection CT 65.1 application: LC	
000983		
REF ALP.I	Container for ligator application: LC	
001213		


REF ACTE.A	Dewar clamp		
000093	application: CTEH, CTEHG, CTEHm, CTEHgm, CTEHgs		
REF ACTE.B	Head / handpiece holder		
000392 005238	application: CTEH, CTEHG, CTEHm, CTEHgm, CTEHgs		
REF ACTE.C	Protective gloves		
001560	application: CTEH, CTEHG, CTEHm, CTEHgm, CTEHgs		
REF ACTE.D	Protective goggles		
001389	application: CTEH, CTEHG, CTEHm, CTEHgm, CTEHgs		
REF ACTE.E	Handpiece holder		
004208	application: CTEH, CTEHG		
REF ACTE.F	Distancer for handpiece holder		
005414	application: CTEH, CTEHG		
REF ACTE.G	Silicone tank cover		
004455	application: CTEH, CTEHG, CTEHm, CTEHgm, CTEHgs		
REF ACTE.H	Dewar cap		
002121	application: CTEH, CTEHG, CTEHm, CTEHgm, CTEHgs		
REF ACTE.I	Cryogenic hose	L=2.4m	
PA02N200002	with application handpiece mini – handpiece without temperature measurement application: CTEHm, CTEHgm, CTEHgs		
REF ACTE.J	Cryogenic hose	L=2.4m	
PA02N200003	with application handpiece laser – handpiece with temperature measurement application: CTEH, CTEHG		

REF SEE TABLE	Head with heater Mini G application: CTEHgm, CTEHgs	L (mm) 	<table><tr><td></td><td></td><td>230V</td><td>110V</td></tr><tr><td>TR11</td><td>L=526</td><td>AGE.A 004561</td><td>AGE.A2 004562</td></tr><tr><td>TR26</td><td>L=606</td><td>AGE.B 004560</td><td>AGE.B2 -----</td></tr><tr><td>TR35</td><td>L=576</td><td>AGE.C 004691</td><td>AGE.C2 004692</td></tr></table>			230V	110V	TR11	L=526	AGE.A 004561	AGE.A2 004562	TR26	L=606	AGE.B 004560	AGE.B2 -----	TR35	L=576	AGE.C 004691	AGE.C2 004692
			230V	110V															
TR11	L=526	AGE.A 004561	AGE.A2 004562																
TR26	L=606	AGE.B 004560	AGE.B2 -----																
TR35	L=576	AGE.C 004691	AGE.C2 004692																
SEE TABLE																			
REF SEE TABLE	Head with heater G application: CTEHG	L (mm) 	<table><tr><td></td><td></td><td>230V</td><td>110V</td></tr><tr><td>TR35</td><td>L=576</td><td>AGE.D 004558</td><td>AGE.D2 004559</td></tr><tr><td>TR60</td><td>L=781</td><td>AGE.E 004556</td><td>AGE.E2 004557</td></tr></table>			230V	110V	TR35	L=576	AGE.D 004558	AGE.D2 004559	TR60	L=781	AGE.E 004556	AGE.E2 004557				
			230V	110V															
TR35	L=576	AGE.D 004558	AGE.D2 004559																
TR60	L=781	AGE.E 004556	AGE.E2 004557																
SEE TABLE																			
REF SEE TABLE	Head P / P Mini application: CTEH / CTEHm	L (mm)  ELEPHANT P  ELEPHANT P Mini	<table><tr><td></td><td></td><td>P MINI</td><td>P</td></tr><tr><td>TR26</td><td>L=600</td><td>APE.A 004697</td><td></td></tr><tr><td>TR35</td><td>L=570</td><td>APE.B 004697</td><td>APE.C 004566</td></tr><tr><td>TR60</td><td>L=765</td><td></td><td>APE.D 004566</td></tr></table>			P MINI	P	TR26	L=600	APE.A 004697		TR35	L=570	APE.B 004697	APE.C 004566	TR60	L=765		APE.D 004566
			P MINI	P															
TR26	L=600	APE.A 004697																	
TR35	L=570	APE.B 004697	APE.C 004566																
TR60	L=765		APE.D 004566																
SEE TABLE																			

REF ACTE.K	O-ring for cryogenic hose size: 21x3mm application: ACTE.I, ACTE.J	21x3mm NBR 70	
000583			
REF ACTE.L	O-ring for head P air connection size: 13.95x2.62mm application: Head P, P MINI	13.95x2.62mm SILICON	
001286			
REF ACTE.M	O-ring for head G and P tank size: 59.99x2.62mm application: Head P, P MINI, G, G MINI	59.99x2.62mm SILICON	
000585			
REF ACTE.N	Ring for head / handpiece holder application: ACTE.B	Ø 65mm	
005238			

REF ALD.A	Protective goggles FG1 1470nm application: SWING lasers	
TA05N500004		
REF ALD.B	Protective goggles YG3 940-980nm application: SALSA lasers	
TA05N500005		
REF ALD.C	Protective goggles EC2 10600nm application: CO2 lasers	
TA05N500003		
REF ALD.D	Footswitch with protective shield application: SALSA and SWING lasers	
SERWIS		
REF ALD.E	Interlock connection cover application: SALSA and SWING lasers	

REF ALD.F	Diode SMA connection cover application: SALSA and SWING lasers	UNC 1/4" - 36 Ø 10, L=25mm 

REF ALL.A	Set for laser lipolysis application: SWING lasers	- ULL handpiece - spare grip for optical fiber - cannula 20cm - cannula 15cm - cannula 10cm - case for the set	
005229 TA05N900005			
REF ALL.B	ULL handpiece for laser lipolysis application: ALL.A		

REF ALL.C	Grip for optical fiber application: ALL.A		

REF ALL.D	Cannula for optical fiber size: 0.8-1mm, L=20cm application: ALL.A		

REF ALL.E	Cannula for optical fiber size: 0.8-1mm, L=15cm application: ALL.A		

REF ALL.F	Cannula for optical fiber size: 0.8-1mm, L=10cm application: ALL.A		

REF AATO.A	Oxygen reducer Farum or Korgiel	
000973	application: ATO-3	
REF AATO.B	Connector device inlet-outlet	
-----	application: ATO-3	
REF UOATO.A	O-ring device inlet-outlet	
000569	size: ~Ø11mm application: ATO-3	
REF AATO.C	Hose reducer-device	
002261	length: ~4m application: ATO-3	
REF AATO.D	Hose oxygen concentrator-device	
-----	length: ~1.5m application: ATO-3	
REF AATO.M	Safety container green	
SERVICE	application: ATO-3	
REF AATO.N	Safety container red	
SERVICE	application: ATO-3	
REF UOATO.B	O-ring for safety container	
SERVICE	size: 34x2mm application: AATO.E, AATO.F	

REF AATO.G	Luer connection CT67	
000984	application: ATO-3	
REF AATO.H	Antibacterial filter 0.45µm	
610-06308	application: ATO-3	
REF AATO.I	3D stopcock	
002545	Luer lock extender for infusion pump application: ATO-3	
REF AATO.J	LC-1 line	
005055	application: ATO-3	
REF AATO.K	OP-3 line	
005054	application: ATO-3	
REF AATO.L	LR-1 line	
005056	application: ATO-3	
REF AATO.E	SO-3 line	
005488	application: ATO-3	
REF AATO.F	LR-1 sleeve connector	
001304	application: ATO-3	



REF AAOZ.A	Container for fluid ozonation	
-----	application: ATO-3	
REF AAOZ.B	S0 ozone pistol	
-----	application: ATO-3	
REF AAOZ.C	Ozone container cap	
-----	application: ATO-3	
REF UOATO.C	O-ring for container cap	24.99x3.53mm SILICON
-----	application: ATO-3	
REF AAOZ.E	S0 container	
-----	application: ATO-3	
REF UOATO.D	Set of o-rings for S0 container	85x3mm SILICON 94x3mm SILICON
-----	size: 85x3mm and 94x3mm application: ATO-3	

REF ABOA.A	Plug for BOA device	
001253	application: BOA Max 2	
REF ABOA.B	12-pin connector	
002335 002341 002340 002388 002340 002872	application: MBX.A, MBX.B, MBX.C, MBX.D, MBX.E, MBX.F	
REF ABOA.C	Plug for device outlet	Ø 5x15mm 
SERWIS	application: BOA Mini+	
REF MBC.A	5-chamber arm sleeve left arm	5 chambers 
PA04N120004	application: BOA Mini+	
REF MBC.B	5-chamber arm sleeve right arm	
PA04N120005	application: BOA Mini+	
REF MBC.C	5-chamber leg sleeve left leg	5 chambers 
PA04N120002	application: BOA Mini+	
REF MBC.D	5-chamber leg sleeve right leg	
PA04N120003	application: BOA Mini+	
REF MBC.E	10-chamber pants type sleeve	10 chambers 
PA04N120001	application: BOA mini+ (10-outlets version)	
REF MBC.F	Leg sleeve expander	
TA04N100010 (left) TA04N100016 (right)	application: MBC.C, MBC.D	
REF MBC.G	Arm sleeve expander	
TA04N100011	application: BOA Mini+	

REF MBC.H	Set of expanders for 10-chamber pants type sleeve	
TA04N100012	2 pieces in the set application: MBC.E	
REF MBX.A	12-chamber arm sleeve left arm	12 chambers 
PA04N110005	application: BOA Max 2	
REF MBX.B	12-chamber arm sleeve right arm	12 chambers 
PA04N110006	application: BOA Max 2	
REF MBX.C	12-chamber leg sleeve left leg	12 chambers 
PA04N110003	application: BOA Max 2	
REF MBX.D	12-chamber leg sleeve right leg	12 chambers 
PA04N110004	application: BOA Max 2	
REF MBX.I	Leg sleeve expander	
TA04N100013 (left) TA04N100016 (right)	application: MBX.C, MBX.D	
REF MBX.J	Arm sleeve expander	
TA04N100019 (left) TA04N100018 (right)	application: MBX.A, MBX.B	
REF MBX.E	24-chamber jacket type sleeve	24 chambers 
PA04N110001	application: BOA Max 2	
REF MBX.F	24-chamber pants type sleeve	24 chambers 
PA04N110002	application: BOA Max 2	
REF MBX.G	Set of expanders for jacket type sleeve	
TA04N100009	3 pieces in the set application: MBX.E	
REF MBX.H	Set of expanders for pants type sleeve	
PA04N110008	2 pieces in the set application: MBX.F	

REF AMB.A	Disposable insert for jacket type sleeve	illustrative photo 
TA04N100001	application: MBX.E	
REF AMB.B	Disposable insert covered with foil for jacket type sleeve	
TA04N100002	application: MBX.E	
REF AMB.C	Disposable insert for pants type sleeve	
TA04N100007	application: MBX.F, MBC.E	
REF AMB.D	Disposable insert covered with foil for pants type sleeve	
TA04N100008	application: MBX.F, MBC.E	
REF AMB.E	Disposable insert for leg sleeve	
TA04N100003	application: MBX.C, MBX.D, MBC.C, MBC.D	
REF AMB.F	Disposable insert covered with foil for leg sleeve	
TA04N100004	application: MBX.C, MBX.D, MBC.C, MBC.D	
REF AMB.G	Disposable insert for arm sleeve	
TA04N100005	application: MBX.A, MBX.B, MBC.A, MBC.B	
REF AMB.H	Disposable insert covered with foil for arm sleeve	
TA04N100006	application: MBX.A, MBX.B, MBC.A, MBC.B	

REF ACTC.A	Condensate bottle	
000004	capacity: 500ml application: CRYO-T Cooler	
REF ACTC.B	Cryogenic hose	L=1.8m 
PA02N200004	with application handpiece handpiece without temperature measurement application: CRYO-T Cooler	
REF ACTC.C	Cryogenic hose	L=1.8m 
PA02N200005	with application handpiece handpiece with temperature measurement application: CRYO-T Cooler	
REF ACTC.D	Hose holder for pantograph	
000625	application: CRYO-T Cooler	
	Pantograph	
	application: CRYO-T Cooler	

REF B520.A	Fuse 160mA T	5x20
000362	size: 5x20 application: CT2, CTD	
REF B520.B	Fuse 1A T	
004170	size: 5x20 application: PTC, PTX2	
REF B520.C	Fuse 1.25A T	
001767	size: 5x20 application: ATO-3, CSEA, CSEP, CTEH (110V)	
REF B520.D	Fuse 800mA T	
003078	size: 5x20 application: CTEHm	
REF B520.E	Fuse 630mA T	
001567	size: 5x20 application: CTEH, CTC	
REF B520.F	Fuse 2.5A T	
002215	size: 5x20 application: CTEHG, CTEHgm, CTEHgs	
REF B520.G	Fuse 5A T	
002785	size: 5x20 application: CTEHG (110V), CTEHgm/S (110V), SWING, SALSA (except SALSA 250)	
REF B520.H	Fuse 4A T	
006074	size: 5x20 application: SWING, TWIST Diode	
REF B520.I	Fuse 6.3A T	
	size: 5x20 application: Salsa 250	
REF PZE.E	Electrical cord type E/F	
002661	L=3m, standard: PL, EU application: all devices except CSM and CSC	
REF PZE.B	Electrical cord type B	
002746	L=5m, standard: US, South America application: all devices except CSM and CSC	
REF PZE.G	Electrical cord type G	
003539	L=5m, standard: UK, Asia application: all devices except CSM and CSC	



METRUM CRYOFLEX

Manufacture

Zielna 29
05-082 Blizne Łaszczyńskiego
POLAND, EU

Headoffice

Kolejowa 16A
05-092 Łomianki
POLAND, EU

Tel: +48 22 33 13 750
+48 22 33 13 830
Fax: +48 22 33 13 766

export@metrum.com.pl
www.metrum.com.pl