


Disposable AV Fistula Needle Sets



VITAL HEALTHCARE SDN.BHD.

Add: Lot 3, Jalan Sultan Mohamed 3, Bandar Sultan Sulaiman,
42000 Pelabuhan Klang, Selangor Darul Ehsan, Malaysia
Tel: +603-3179 0722 **Fax:** +603-3179-0822

| Gauge | Needle Length | Wing | Tube Spec | Sterilization |
|-------|---------------|-----------------|-------------|---------------|
| 14G | 20±2mm | fixed/rotatable | ID3.2-OD5.0 | EO/Gamma |
| 14G | 25±2mm | fixed/rotatable | ID3.2-OD5.0 | EO/Gamma |
| 14G | 32±2mm | fixed/rotatable | ID3.2-OD5.0 | EO/Gamma |
| 15G | 20±2mm | fixed/rotatable | ID3.2-OD5.0 | EO/Gamma |
| 15G | 25±2mm | fixed/rotatable | ID3.2-OD5.0 | EO/Gamma |
| 15G | 32±2mm | fixed/rotatable | ID3.2-OD5.0 | EO/Gamma |
| 16G | 20±2mm | fixed/rotatable | ID3.2-OD5.0 | EO/Gamma |
| 16G | 25±2mm | fixed/rotatable | ID3.2-OD5.0 | EO/Gamma |
| 16G | 32±2mm | fixed/rotatable | ID3.2-OD5.0 | EO/Gamma |
| 17G | 20±2mm | fixed/rotatable | ID3.2-OD5.0 | EO/Gamma |
| 17G | 25±2mm | fixed/rotatable | ID3.2-OD5.0 | EO/Gamma |
| 17G | 32±2mm | fixed/rotatable | ID3.2-OD5.0 | EO/Gamma |

- Ultra thin and siliconized wall for smooth penetration
- Special process on blade reduces pain and skin damage
- Medical grade PVC tubing minimizes side effect to health