

RADIAL ARTERY INTRODUCER

Set with hemostasis valve to introduce and exchange catheters and endocavitary electrodes



Set components:

- sheath with a valve
- angiographic needle
- guidewire
- dilator

Lot 165

The set contains a guidewire of special construction. The soft distal part of the guidewire enables its introduction to the vessel, and the stiff proximal part provides an excellent scaffold for the dilator, along with preventing vessel perforation. The distal part of the sheath and dilator is finished in a way that enables an easy and atraumatic introduction of the catheter to the vessel.

Reference:	Sheath size	Sheath length [mm]	Colour identification	Dilator size	Angiographic needle size	Guidewire size [inch]	Guidewire length [cm]	Guidewire type
<u>INT4FKR</u>	4F	70	red	4F	21G	<u>.018</u>	45	"J"
INT4FKRP						.018		STRAIGHT
<u>INT5FKR</u>	5F	110	grey	5F		<u>.018</u>	45	"J"
INT5FKRP						.018		STRAIGHT
INT5F21KRP						.021		STRAIGHT
<u>INT5FKR23</u>		230				<u>.018</u>	60	"J"
INT5FKR23P						.018		STRAIGHT
INT5F21KR23P						.021		STRAIGHT
<u>INT6FKR</u>	6F	110	green	6F		<u>.018</u>	45	"J"
INT6FKRP						.018		STRAIGHT
INT6F21KRP						.021		STRAIGHT
<u>INT6FKR23</u>		230				<u>.018</u>	60	"J"
INT6FKR23P						.018		STRAIGHT
INT6F21KR23P						.021		STRAIGHT
INT7FKR	7F	110	orange	7F		.018	45	"J"
INT7FKRP						.018		STRAIGHT
INT7F21KRP						.021		STRAIGHT
INT7FKR23		230				.018	60	"J"
INT7FKR23P						.018		STRAIGHT
INT7F21KR23P						.021		STRAIGHT

When ordering use Ref. Number, e.g.: INT5FKR.

Special order can contain introducer with radiopaque marker

Model reference number: INT7FKR**M** – introducer with a marker.

Hydrophilic coating is available on request.

Model reference number: INT7FKR**H** – introducer with a hydrophilic surface protection.

