

Clot Activator

Clot Activator Tubes are covered with clot activator powder to accelerate blood coagulation.

Silica (clot activator) is used as the activator.

Serum is a substance freed from the fibrinogen but containing fragments of thrombocyte and product of the metabolism.

Usage: to research serum in the clinical chemistry, immunology, serology, protein electrophoresis.

Specimen: serum.

Clotting time: from 10 to 30 minutes.

Conditions of centrifugation: 1500-1800 g for 10 minutes.

REF	Tube Size (mm)	Draw Volume (mL)	Tube Materia	Pack/Carton
C104035			PET	100PCS/1200PCS
C304035	13x75	3.5mL	GLASS	100PCS/1200PCS
C104040			PET	100PCS/1200PCS
C304040	13x75	4mL	GLASS	100PCS/1200PCS
C104050			PET	100PCS/1200PCS
C304050	13x100	5mL	GLASS	100PCS/1200PCS
C104060			PET	100PCS/1200PCS
C304060	13x100	6mL	GLASS	100PCS/1200PCS
C104070			PET	100PCS/1200PCS
C304070	13x100	7mL	GLASS	100PCS/1200PCS
C104080			PET	100PCS/1200PCS
C304080	16x100	8mL	GLASS	100PCS/1200PCS
C104090			PET	100PCS/1200PCS
C304090	16x100	9mL	GLASS	100PCS/1200PCS
C104100			PET	100PCS/1200PCS
C304100	16x100	10mL	GLASS	100PCS/1200PCS



Gel & Clot Activator

As the activator it is used silica dioxide and biologically inert olefin gel.

Gel provides the separation of serum and clot up to 48 hours without re- centrifugation.

Gel is the special material intended for formation of barrier between cellular components of blood and serum during centrifugation.

Usage: to research serum in the clinical chemistry, immunology, serology, protein electrophoresis.

Specimen: serum.

Clotting time: from 5 to 30 minutes.

Conditions of centrifugation: 1800-2000 g for 10 minutes.

REF	Tube Size (mm)	Draw Volume (mL)	Tube Materia	Pack/Carton
C105035			PET	100PCS/1200PCS
C305035	13x75	3.5mL	GLASS	100PCS/1200PCS
C305040	13x75	4mL	GLASS	100PCS/1200PCS
C105050			PET	100PCS/1200PCS
C305050	13x100	5mL	GLASS	100PCS/1200PCS
C105060			PET	100PCS/1200PCS
C305060	13x100	6mL	GLASS	100PCS/1200PCS
C105070			PET	100PCS/1200PCS
C305070	13x100	7mL	GLASS	100PCS/1200PCS
C105080			PET	100PCS/1200PCS
C305080	16x100	8mL	GLASS	100PCS/1200PCS
C105090			PET	100PCS/1200PCS
C305090	16x100	9mL	GLASS	100PCS/1200PCS
C105100			PET	100PCS/1200PCS
C305100	16x100	10mL	GLASS	100PCS/1200PCS

