

No Additive (Tube with Red cap)



- This tube is a plain containing no anticoagulant.
- It mainly used for collection of Serum for selected chemistry tests

Ex. Hepatitis C Virus (HCV Ab), Hormones and blood chemistry Tests.

2 ml

3 ml

4 ml

5 ml

6 ml

Clot Activator (Tube with Red cap)

- This tube contains BCA as coagulation Catalyst.
- It mainly used for collection of serum for various laboratory test.

Ex. Serology tests and blood chemistry

2 ml

3 ml

4 ml

5 ml

6 ml



EDTA K2 (Tube with lavender cap)



- This tube contains Dipotassium Ethyleneimine tetra acetic acid as anticoagulant.
- It mainly used for hematological tests.

Ex. Complete Blood Picture (CBC)

2 ml

3 ml

4 ml

5 ml

6 ml

