



SW-500

Rebound Tonometer



New Technologies And Powerful Functions
Simple operation interface provides
with accurate operation experience.

SHANGHAI LINK INSTRUMENTS

Features

- With vertical and horizontal two working modes, wireless output print data. The equipment is used to measure intraocular pressure, using the principle of: the probe hits the surfaces of different hardness at a certain speed, has different reaction when the probe rebounds.
- Be of advantages of high accuracy, portable, without anesthesia, without the cross-infection etc.



Specification

Measuring Range: 3mmHg~70mmHg

Measuring Error: $\pm 1.5\text{mmHg}$ ($3\text{mmHg} \leq \text{intraocular pressure} \leq 25\text{mmHg}$)

$\pm 2.5\text{mmHg}$ ($25\text{mmHg} \leq \text{intraocular pressure} \leq 70\text{mmHg}$)

Measurement can go along both in standing and lying posture

Wireless data printing.

Easy to learn, easy to use

Lightweight and handy

No anesthesia, avoiding discomfort reaction



Add.No.1101 Huyi Highway, Shanghai, 201802, China
Tel. +86-21-69512871
Fax. +86-21-69512872
E-mail. sales@linkcsh.com
Http://www.cnlink-group.com

SHANGHAI LINK INSTRUMENTS CO., LTD.