


ADJUSTMENT OF DRAWING

# BLOOD TUBING SET.



Form No.BL-PI-020-02

MODEL:NO.	A290R-V801R	TYPE	(A)
-----------	-------------	------	-----



**REMARK**  
 GAMMA STERILE (CAUTION)  
 4.3x6.8 (65)  
 3.5x5.5 (65)  
 1.0x2.1 (57)  
 8.0x12.0 (87)  
 2.0x3.65 (57)

\* Pump tube cut length 352 mm.

Packing Method	
Item	GAMMA Sterilization
Individual Package	CE GAMMA Package
Outer Box	CE GAMMA
Attach sticker	Attach sticker barcode on outer box
IFU	CE GAMMA
Packing Q'ty	24 Sets.
Special Labelling M/C Name	FOR B. BRAUN DIALOG DN
PRIMING VOLUME A290R-V801R	151 ML

Change symbol "/" to "-"  
 A290R/V801R → A290R-V801R

ESTABLISHED DATE	DRAWING	CHECK	APPROVE
10 Jun.'11	Ms. Yanisa	<i>Yanisa</i> 4 Apr.'17	<i>[Signature]</i> 24 Apr.'17 SEE BACK SIDE FOR APPROVED BY NIPRO

RECORD OF CHANGE			
REVISED DATE	REVISION	APPR.NO.	DATE OF RECEIPT
05 Apr.'14	8	DC-BL-14D-008	24 Apr. 2017
09 Jun.'16	9	DC-BL-16F-009	
31 Mar.'17	10	DC-BL-17D-011	FILE NO. 252-833
			REGISTERED <i>[Signature]</i> 24 Apr. 17

BASE FROM MODEL X2186R

Form executed Date on 01 Apr.2014

CONFIDENTIAL

NIPRO