




Material. Design. Optics.



IOL Portfolio



www.medicontur.com

About us

Our company is constantly striving for innovation in order to help patients with various vision related problems.

Since its founding in 1989 Medicontur has developed into one of the major players in the field of cataract and refractive surgery.

Our lens portfolio provides a solution for a range of eye conditions including cataract, presbyopia, astigmatism, dry AMD, and colour vision deficiency.

30+

33 years
of experience

7+m

intraocular lenses
implanted

over 60

international
markets

Active presence in




Medicontur presence worldwide

TRIFOCAL IOLS

PRELOADED

NON-PRELOADED



Trifocal
LibERT7
Accuject PRO

677CMY


Liberty Preloaded



Trifocal
Toric
LibERT7
Accuject PRO

677CMTY

Liberty Toric Preloaded



Q | Trifocal
FLEX
Accuject PRO

640CMY

Q-Flex Trifocal Preloaded



Toric
LibERT7

677MY

Liberty

Standard Powers	+8.0 D → +30.0 D (0.5 D steps)		+8.0 D → +30.0 D (0.5 D steps)		+8.0 D → +30.0 D (0.5 D steps)		+8.0 D → +30.0 D (0.5 D steps)	
Extreme Powers	+30.5 D → +35.0 D (0.5 D steps)		+30.5 D → +35.0 D (0.5 D steps)		+30.5 D → +35.0 D (0.5 D steps)		+30.5 D → +35.0 D (0.5 D steps)	
Cylinders	-		1.0 D 1.5 D 2.0 D 2.5 D 3.0 D 3.5 D 4.0 D 4.5 D 525D* 6.0 D*		-		-	
Addition	+3.5 D near +1.75 D intermediate		+3.5 D near +1.75 D intermediate		+3.5 D near +1.75 D intermediate		+3.5 D near +1.75 D intermediate	
Material	Copolymer of hydrophobic and hydrophilic monomers with 25% water content, UV blocker and blue light filter		Copolymer of hydrophobic and hydrophilic monomers with 25% water content, UV blocker and blue light filter		Copolymer of hydrophobic and hydrophilic monomers with 25% water content, UV blocker and blue light filter		Copolymer of hydrophobic and hydrophilic monomers with 25% water content, UV blocker and blue light filter	
RI Abbe Number	1.46	58	1.46	58	1.46	58	1.46	58
Optic Design	Aspheric - Aberration Neutral Biconvex (+8.0 D → +35.0 D)		Aspheric - Aberration Neutral Biconvex (+8.0 D → +35.0 D)		Aspheric - Aberration Neutral Biconvex (+8.0 D → +35.0 D)		Aspheric - Aberration Neutral Biconvex (+8.0 D → +35.0 D)	
Diffraction or EDOF Zone	3 mm diameter diffractive array on the anterior surface utilizing EPS technology		3 mm diameter diffractive array on the anterior surface utilizing EPS technology		3 mm diameter diffractive array on the anterior surface utilizing EPS technology		3 mm diameter diffractive array on the anterior surface utilizing EPS technology	
Optic Diameter	6 mm		6 mm		6 mm		6 mm	
Overall Diameter	13 mm		13 mm		11.0 mm (8.0 D → +15.0 D) 10.7 mm (+15.5 D → +35.0 D)		13 mm	
Haptic Design	Posterior vaulting fenestrated C-loops with 0° angulation		Posterior vaulting fenestrated C-loops with 0° angulation		Posterior vaulting 4 closed loops with 0° angulation		Posterior vaulting fenestrated C-loops with 0° angulation	
PCO Prevention	360° sharp edge		360° sharp edge		360° sharp edge		360° sharp edge	
Shelf Life	3 years after sterilization		3 years after sterilization		3 years after sterilization		5 years after sterilization	

*available above +10.0 D SEQ

EDOF IOL

PRELOADED



Trifocal
Toric
LIBERTY7

677MTY

Liberty Toric

+8.0 D → +30.0 D (0.5 D steps)

+30.5 D → +35.0 D (0.5 D steps)

1.0 D 1.5 D 2.0 D 2.5 D 3.0 D
3.5 D 4.0 D 4.5 D 5.25 D* 6.0 D*

+3.5 D near
+1.75 D intermediate

Copolymer of hydrophobic and
hydrophilic monomers with 25%
water content, UV blocker and blue
light filter

1.46

58

Aspheric - Aberration Neutral
Biconvex (+8.0 D → +35.0 D)

3 mm diameter diffractive array on
the anterior surface utilizing
EPS technology

6 mm
13 mm

Posterior vaulting fenestrated
C-loops with 0° angulation

360° sharp edge

5 years after sterilization



Trifocal
Q FLEX

640MY

Q-Flex Trifocal

+8.0 D → +30.0 D (0.5 D steps)

+30.5 D → +35.0 D (0.5 D steps)

-

+3.5 D near
+1.75 D intermediate

Copolymer of hydrophobic and
hydrophilic monomers with 25%
water content, UV blocker and blue
light filter

1.46

58

Aspheric - Aberration Neutral
Biconvex (+8.0 D → +35.0 D)

3 mm diameter diffractive array on
the anterior surface utilizing
EPS technology

6 mm
11.0 mm (8.0 D → +15.0 D)
10.7 mm (+15.5 D → +35.0 D)

Posterior vaulting 4 closed loops
with 0° angulation

360° sharp edge

5 years after sterilization

EDOF by ELON

Bi-Flex platform



3! **elON**
FLEX | POB-MA HYDROPHOBIC PRELOADED

877PEY

Bi-Flex ELON POB-MA

+8.0 D → +30.0 D (0.5 D steps)

+31.0 D → +35.0 D (1.0 D steps)

-

-

Hydrophobic acrylic with UV
blocker and blue light filter

1.47

58

Aspheric - Aberration Neutral
Biconvex (+8.0 D → +35.0 D)

Refraction-based central
Wavefront Linking zones

6 mm
13 mm

Posterior vaulting fenestrated
C-loops with 0° angulation


360° sharp edge

2.5 years after sterilization

MONOFOCAL IOLS

PRELOADED


NON-PRELOADED



3! POB-MA
FLEX | HYDROPHOBIC PRELOADED

877PA | 877PAY


Bi-Flex POB-MA



3! PIL-MA
FLEX | HYDROPHILIC PRELOADED

677P | 677PY


Bi-Flex Preloaded



Q | PIL-MA
FLEX | HYDROPHILIC PRELOADED

640P | 640PY

Q-Flex Preloaded




3! FLEX HB

877FAB | 877FABY

Bi-Flex HB

Standard Powers	0.0 D → +9.0 D (1.0 D steps) +10.0 D → +30.0 D (0.5 D steps)		0.0 D → +30.0 D (0.5 D steps)		0.0 D → +9.0 D (1.0 D steps) +10.0 D → +30.0 D (0.5 D steps)		0.0 D → +9.0 D (1.0 D steps) +10.0 D → +30.0 D (0.5 D steps)	
Extreme Powers	-10.0 D → -1.0 D (1.0 D steps) +31.0 D → +35.0 D (1.0 D steps)		-10.0 D → -1.0 D (1.0 D steps) +31.0 D → +35.0 D (1.0 D steps)		+31.0 D → +35.0 D (1.0 D steps)		+31.0 D → +35.0 D (1.0 D steps)	
Cylinders	-		-		-		-	
Addition	-		-		-		-	
Material	Hydrophobic acrylic with UV blocker and optional blue light filter		Copolymer of hydrophobic and hydrophilic monomers with 25% water content, UV blocker and optional blue light filter		Copolymer of hydrophobic and hydrophilic monomers with 25% water content, UV blocker and optional blue light filter		Hydrophobic acrylic with UV blocker and optional blue light filter	
RI Abbe Number	1.47	58	1.46	58	1.46	58	1.47	58
Optic Design	Aspheric - Aberration Neutral Biconcave (-10.0 D → -1.0 D) Biconvex (0.0 D → +35.0 D)		Aspheric - Aberration Neutral Convex-concave (-10.0 D → -1.0 D) Biconvex (0.0 D → +35.0 D)		Aspheric - Aberration Neutral Biconvex (0.0 D → +35.0 D)		Aspheric - Aberration Neutral Biconvex (0.0 D → +35.0 D)	
Diffraction or EDOF Zone	-		-		-		-	
Optic Diameter	6 mm		6 mm		6 mm		6 mm	
Overall Diameter	13 mm		13 mm		11.0 mm (0.0 D → +15.0 D) 10.7 mm (+15.5 D → +35.0 D)		13 mm	
Haptic Design	Posterior vaulting fenestrated C-loops with 0° angulation		Posterior vaulting fenestrated C-loops with 0° angulation		Posterior vaulting 4 closed loops with 0° angulation		Posterior vaulting fenestrated C-loops with 0° angulation	
PCO Prevention	360° sharp edge		360° sharp edge		360° sharp edge		360° sharp edge	
Shelf Life	2.5 years after sterilization		3 years after sterilization		3 years after sterilization		5 years after sterilization	



3! FLEX HL

677AD | 677ADY

Bi-Flex

0.0 D → +9.0 D (1.0 D steps)
+10.0 D → +30.0 D (0.5 D steps)

-10.0 D → -1.0 D (1.0 D steps)
+31.0 D → +45.0 D (1.0 D steps)

-

-

Copolymer of hydrophobic and hydrophilic monomers with 25% water content, UV blocker and optional blue light filter

1.46

58

Aspheric - Aberration Neutral
Convex-concave (-10.0 D → -1.0 D)
Biconvex (0.0 D → +45.0 D)


-

6 mm
13 mm

Posterior vaulting fenestrated
C-loops with 0° angulation

360° sharp edge

5 years after sterilization



Q FLEX

640AD | 640ADY

Q-Flex

0.0 D → +9.0 D (1.0 D steps)
+10.0 D → +30.0 D (0.5 D steps)

+31.0 D → +35.0 D (1.0 D steps)

-

-

Copolymer of hydrophobic and hydrophilic monomers with 25% water content, UV blocker and optional blue light filter

1.46

58

Aspheric - Aberration Neutral
Biconvex (0.0 D → +35.0 D)


-

6 mm
11.0 mm (0.0 D → +15.0 D)
10.7 mm (+15.5 D → +35.0 D)

Posterior vaulting 4 closed loops
with 0° angulation

360° sharp edge

5 years after sterilization



3! FLEX T
Accuject PRO

677CTA | 677CTAY

Bi-Flex Toric Preloaded

8.0 D → +30.0 D (0.5 D steps)

+31.0 D → +35.0 D (1.0 D steps)

1.0 D 1.5 D 2.25 D 3.0 D
3.75 D 4.5 D 5.25 D 6.0 D

-

Copolymer of hydrophobic and hydrophilic monomers with 25% water content, UV blocker and optional blue light filter

1.46

58

Aspheric - Aberration Neutral
Biconvex (8.0 D → +35.0 D)

-

6 mm
13 mm

Posterior vaulting fenestrated
C-loops with 0° angulation

360° sharp edge

3 years after sterilization



3! FLEX T

677TA | 677TAY

Bi-Flex Toric*

0.0 D → +30.0 D (0.5 D steps)

-10.0 D → -1.0 D (1.0 D steps)
+31.0 D → +35.0 D (1.0 D steps)

1.0 D 1.5 D 2.25 D 3.0 D
3.75 D 4.5 D 5.25 D 6.0 D

-

Copolymer of hydrophobic and hydrophilic monomers with 25% water content, UV blocker and optional blue light filter

1.46

58

Aspheric - Aberration Neutral
Convex-concave (-10.0 D → +5.0 D)
Biconvex (+5.5 D → +35.0 D)

-

6 mm
13 mm

Posterior vaulting fenestrated
C-loops with 0° angulation

360° sharp edge

5 years after sterilization

*Extreme cylinders available: +6.75 D, +7.5 D, +8.25 D, +9.0 D, +10.0 D

SECONDARY SULCUS IOLS

NON-PRELOADED



 1stQ AddOn®

A46R

1stQ AddOn® Refractive

 1stQ AddOn®

A45RT

1stQ AddOn® Refractive Toric

 1stQ AddOn®

A45RD2

1stQ AddOn® Trifocal

 1stQ AddOn®

A45DT

1stQ AddOn® Trifocal Toric

Standard Powers	-10.0 D → +10.0 D (0.25 D steps)		-10.0 D → +10.0 D (0.25 D steps)		-5.0 D → +5.0 D (0.25 D steps)		-3.0 D → +3.0 D (0.5 D steps)	
Extreme Powers	-		-		-		-	
Cylinders	-		+1.0 D +1.5 D → +9.0 D (0.75 D steps) +10.0 D, +11.0 D*		-		+1.0 D → +4.5 D (0.5 D steps)	
Addition	-		-		+3.0 D near +1.5 D intermediate		+3.0 D near +1.5 D intermediate	
Material	Copolymer of hydrophobic and hydrophilic monomers with 25% water content and UV blocker		Copolymer of hydrophobic and hydrophilic monomers with 25% water content and UV blocker		Copolymer of hydrophobic and hydrophilic monomers with 25% water content and UV blocker		Copolymer of hydrophobic and hydrophilic monomers with 25% water content and UV blocker	
RI Abbe Number	1.46	58	1.46	58	1.46	58	1.46	58
Optic Design	Aspheric - Aberration Neutral Convex-concave (-10.0 D → +10.0 D)		Aspheric - Aberration Neutral Convex-concave (-10.0 D → +10.0 D)		Aspheric - Aberration Neutral Convex-concave (-5.0 D → +5.0 D)		Aspheric - Aberration Neutral Convex-concave (-3.0 D → +3.0 D)	
Diffraction or EDOF Zone	-		-		3 mm diameter diffractive array on the anterior surface utilizing EPS technology		3 mm diameter diffractive array on the anterior surface utilizing EPS technology	
Optic Diameter	6 mm		6 mm		6 mm		6 mm	
Overall Diameter	13 mm		13 mm		13 mm		13 mm	
Haptic Design	Special 4 closed loops with 0° angulation		Special 4 closed loops with 0° angulation		Special 4 closed loops with 0° angulation		Special 4 closed loops with 0° angulation	
PCO Prevention	-		-		-		-	
Shelf Life	5 years after sterilization		5 years after sterilization		5 years after sterilization		5 years after sterilization	

*available for SEQ range: -3.0 D → +8.0 D



1st Q AddOn®

SML | SCHARIOTH
MACULA LENS

A45SML

SML

-4.0 D → +4.0 D (0.5 D steps)

-

-

+10 D

Copolymer of hydrophobic and hydrophilic monomers with 25% water content and UV blocker

1.46

58

Aspheric - Aberration Neutral
Convex-concave (-4.0 D → +4.0 D)

-


6 mm

13 mm

Special 4 closed loops with
0° angulation

-

5 years after sterilization




Constants table

		Nominal A-constant	SRK/T A-constant	Haigis (a ₀)	Haigis (a ₁)	Haigis (a ₂)	Hoffer Q (pACaD)	Holladay I (SF)	Holladay II (ACD)**	Holladay II (SF)**	Barrett Universal II (Lens Factor)**
Liberty Preloaded	677CMY	118,9	118.828*	0.190*	0.192*	0.173*	5.431*	1.682*	5,45	1,69	1,79
Liberty Toric Preloaded	677CMTY	118,9	118.828*	0.190*	0.192*	0.173*	5.431*	1.682*	5,45	1,69	1,79
Q-Flex Trifocal Preloaded	640CMY	118,9	118,9	1,243	0.400	0.100	5,46	1,67	5,49	1,73	1,83
Liberty	677MY	118,9	118.828*	0.190*	0.192*	0.173*	5.431*	1.682*	5,45	1,69	1,79
Liberty Toric	677MTY	118,9	118.828*	0.190*	0.192*	0.173*	5.431*	1.682*	5,45	1,69	1,79
Q-Flex Trifocal	640MY	118,9	118,9	1,243	0.400	0.100	5,46	1,67	5,49	1,73	1,83
Bi-Flex ELON POB-MA	877PEY	118,9	118,9	1.320	0.400	0.100	5,46	1,7	5,49	1,73	1,83
Bi-Flex POB-MA	877PA(Y)	118,9	118,9	1.320	0.400	0.100	5,46	1,7	5,49	1,73	1,83
Bi-Flex Preloaded	677P(Y)	118,9	118.828*	0.190*	0.192*	0.173*	5.431*	1.682*	5,45	1,69	1,79
Q-Flex Preloaded	640P(Y)	118,9	118,9	1,243	0.400	0.100	5,46	1,67	5,49	1,73	1,83
Bi-Flex HB	877FAB(Y)	118,9	118,9	1.320	0.400	0.100	5,46	1,7	5,49	1,73	1,83
Bi-Flex	677AD(Y)	118,9	118.828	0.190	0.192	0.173	5.431	1.682	5,45	1,69	1,79
Q-Flex	640AD(Y)	118,9	118,9	1,243	0.400	0.100	5,46	1,67	5,49	1,73	1,83
Bi-Flex Toric Preloaded	677CTA(Y)	118,9	118.828*	0.190*	0.192*	0.173*	5.431*	1.682*	5,45	1,69	1,79
Bi-Flex Toric	677TA(Y)	118,9	118.828*	0.190*	0.192*	0.173*	5.431*	1.682*	5,45	1,69	1,79

IOL-Injector-OVD compatibility chart available: <https://medicontur.com/professionals/compatibility>

* Optimized IOL constants: n=400, date: 2018.
** Barrett Universal II and Holladay II constants were calculated with https://calc.apacs.org/barrett_universal2105/ and <http://www.hicsoap.com> online calculators.



Focusing on your
vision since 1989.

Some of the products shown in this IOL chart may not be registered in your country and may not be available there.
Please contact our local distributor for availability of our products in your market.

The trademarks and logos used in this document are the property of MEDICONTUR Medical Engineering Ltd.
©2020 MEDICONTUR Medical Engineering Ltd. All rights reserved.



Material. Design. Optics.

Head Office - Zsámbék

HU-2072 Herceghalmi Road 1
Tel: +36 23 56 55 55
mc@medicontur.hu

Sales Office - Budapest

Invoice: HU-1113 Karolina Street 65
Mailing: HU-1113 Daróczi Street 80
Tel: +36 1 214 20 33
iol@medicontur.com

International Office - Plan-les-Ouates

CH-1228 Chemin du Champ-des-Filles 36
export@medicontur.com