


HEALTHCARE PROFESSIONALS

Sentrant Introducer Sheath

EVAR and TEVAR Procedures



REQUEST MORE INFORMATION

CONTACT US

Icons: Print, Email, Link

AORTIC CATALOGUE

Choose from a complete portfolio of market-leading aortic vascular therapies to treat aneurysm disease.

DOWNLOAD CATALOGUE

RELATED PRODUCTS

Aortic Stent Graft Products

Cardiovascular Products

OVERVIEW

Achieve improved seal in both straightforward and complex anatomies with Sentrant™ introducer sheath.\*

INDICATIONS

See the device manual for detailed information regarding the instructions for use, indications, contraindications, warnings, precautions, and potential adverse events. For further information, contact your local Medtronic representative and/or consult the [Medtronic Europe website](#).

The Medtronic Sentrant™ Introducer Sheaths with Hydrophilic Coating are intended to provide a conduit for the insertion of diagnostic or endovascular devices into the vasculature and to minimise blood loss associated with such insertions.



SENTRANT INTRODUCER SHEATH FEATURES

The Sentrant introducer sheath is the choice for superior hemostasis.\*

1. Smooth dilator to sheath transition
2. Radiopaque marker band
3. Reinforced coiled tubing
4. Two working lengths: 28 cm and 64 cm
5. Hydrophilic coating
6. Locking mechanism on dilator handle
7. Guidewire diameter: 0.035" / 0.89 mm

Compatible with Endurant™ and Valiant Navion™ stent graft delivery systems

ENGINEERED TO DELIVER PROCEDURAL CONFIDENCE

- EnsureSeal Technology delivers superior leak resistance versus competitors\*
- Coil-reinforced tubing for added stability and kink resistance
- Maintains lubricity after multiple insertions
- Radiopaque dilator shaft and sheath tip for accurate visualization and guidance
- 64 cm configuration available for broader anatomical applicability
- Compatible with a wide range of endovascular portfolios

SUPERIOR LEAK RESISTANCE VERSUS COOK CHECK-FLO PERFORMER™ INTRODUCER SHEATH AND GORE DRYSEAL™ FLEX INTRODUCER SHEATH\*

Comparison of leak resistance for the Medtronic Sentrant introducer sheath, Cook Check Flo Performer introducer, and Gore DrySeal Flex sheath.\* The Sentrant introducer sheath had zero leakage on 0.035" / 0.89 mm guidewire.



MANUALS AND TECHNICAL GUIDES

Instructions for Use

Find this technical manual in the [Medtronic Manual Library](#) or in the product labeling supplied with each device.

MODEL SPECIFICATIONS

ORDERING INFORMATION

Table with 3 columns: Model Number, Sheath Size (F), Working Length (cm)

\* Leak Resistance Bench Test Data on file at Medtronic. Test data not indicative of clinical performance. Bench Test compared Cook Check-Flo Performer™™ Introducer Sheath 16 Fr and Gore DrySeal™™ Flex 16 Fr to Sentrant™ 16 Fr. Graph shows average leakage over an extrapolated 2.5 hour procedure.

™™ Third party brands are trademarks of their respective owners.

Disclaimer: 64 cm working length sheath may not be available in your region or country. Please consult the approved indications for use. Content on specific Medtronic products is not intended for users in markets that do not have authorisation for use.

Medtronic

Engineering the extraordinary

Privacy Statement | Terms of Use | Contact

UC201801391b EE Last Updated May 2023 © 2023 Medtronic