



Angiographic Syringes

- The angiographic syringe of SCW Medicath is used to inject contrast material into the patient's vascular system to diagnostic and therapeutic the vessels by x-ray photography.
- Angiographic syringe is compatible with high-pressure injection equipment during CT, DSA and MRI, and also compatible with a broad range of different brand machines.
- Widely used in cardiology, radiology, and interventional procedures.

[INQUIRY](#)

FEATURES

- **Competitive Pricing:** As a professional [medical supplies manufacturer](#), we leverage supply chain advantages to offer products of the same high quality at more competitive prices compared to peers. [Quote now.](#)

resonance, and imaging. The contrast media used is the gold standard for detecting cancerous tissues, delivering high-definition lesion images.

- **Wide Compatibility:** SCW Medicath's syringes come in a wide variety and are compatible with injectors from nearly all major manufacturers on the market.
 - **High Pressure Resistance:** Designed to handle maximum pressures of up to 350 Pa or 1200 Pa, ensuring reliable performance during procedures.
 - **Customizable Connection Tubes:** Exported products are equipped with a single check valve by default, with two check valves available for special configurations.
 - **Durable Materials:** Made from high-quality medical grade materials resistant to high temperatures and pressures. This angiographic injector and syringe has good permeability.
 - **Enhanced Visibility:** The clear syringe body provides excellent clarity and glass-like finish, allowing easy detection of air bubbles and contrast media.
 - **Export Business:** We provide sample delivery, wholesale, distribution and OEM/ODM services. For business cooperation, our minimum order amount is 1 carton.
- [Quote](#) for excellent products and prices now.

USE METHOD

This video tutorial shows you how to use our angiographic syringe step by step -

Angiographic Syringe Tutorial | Step-by-Step Guide | SCW Medicath



PARAMETER

Product No.		Volume(ml)	Description	
Single	Dual		Used for Medrad	Sin
42.16.10000	42.16.10001	200	Used for CT: MCT & plus; OP-100; VCT610; ECT710	CT-200-
42.16.10002	42.16.10003	200	Used for CT: STELANT	CT-200-
42.16.10004		150	Used for DSA: Mark V&Mark V Provis	DSA-15
42.16.10005		200	Used for DSA: Mark V& Mark V Provis	DSA-20
42.16.10007		130	Used for DSA: Mark III & Mark IV	DSA-13
	42.16.10008	65	Used for MRI: Spectris	
	42.16.10009	115	Used for MRI: Spectris	
	42.16.10010	65/115	Used for MRI: Spectris	
42.16.10011	42.16.10012	190	Used for CT: SALENT	CT-190-
42.16.10016		150	Used for DSA: Mark 7	DSA-15

Product No.		Volume(ml)	Description	
Single	Dual		Used for Libel-Flarsheim	Single
42.16.20000	42.16.20001	200	Used for CT: CT 9000& CT 9000 ADV	CT-200-
42.16.20002		150	Used for DSA: Angiomat 6000	DSA-150
42.16.20003		150	Used for CT & DSA: Angiomat Illumena	CT/DSA-150
42.16.20005	42.16.20006	60	Used for MRI: Optistar	MRI-60-

Product No.		Volume(ml)	Description	
Single	Dual		Used for EZEM	Single
42.16.50000	42.16.50001	200	Used for CT: Empower	CT-200-

Product No.		Volume(ml)	Description	
-------------	--	------------	-------------	--

42.16.30000	42.16.30001	100	Used for CT: A-25	CT-100-
42.16.30002	42.16.30003	200	Used for CT: A-25 A-60	CT-200-
42.16.30005		120	Used for DSA: Nemoto 120S	DSA-120
42.16.30006	42.16.30007	60	Nemoto Dual Shot. Sonic Shot 50	MRI-60-
	42.16.30008	100/200	Used for CT.A-25	
	42.16.30012	60/200	Used for CT: A-25.A-60	

Product No.		Volume(ml)	Description	
Single	Dual		Used for Metron	Single
42.16.40000	42.16.40001	200	Used for CT: ACCUTRON	CT-200-M
42.16.40002	42.16.40003	65	Used for MRI: ACCUTRON	MRI-65-M
	42.16.40004	65/200	Used for MRI: ACCUTRON	

CE CERTIFICATION

This is a CE certified product. Here is the picture for your reference:

Manufacturer:

SCW MEDICATH LTD.

NO.4, Baolong 6th Road, Baolong Industrial
Town, Longgang District, Shenzhen,
518116, Guangdong, P.R. China

We, the manufacturer, herewith declare that the products
Angiographic Syringes

Models: DSA-60-A1, DSA-150-A1, CT-100-X1, CT-200-A1, MRI-115-A1,

DSA-100-B1, DSA-150-B1, CT-200-B1, MRI-65-A1, CT-200-A2,
DSA-200-A1, DSA-130-A1, CT/DSA-150-B1, MRI-60-B1, CT-200-NE,
DSA-90-NE, DSA-120-NE, CT-100-NE, CT-190-A1, CT-190/190-A1,
CT-200/200-A2, CT-200/200-B1, CT-100/100-NE, CT-100/200-NE,
CT-200/200-NE, CT-200-EM, CT-200/200-EM, CT-200-MED,
CT-200/200-MED, MRI-65/115-A1, MRI-65/65-A1, MRI-115/115-A1,
MRI-60/60-B1, MRI-60/60-NE, MRI-65-MED, MRI-65/200-MED

meet the provisions of Directive 93/42/EEC which apply to them.

The medical device has been assigned to class IIa according to Annex IX Rule 2 of the
Directive 93/42/EEC. It bears the mark

CE 0197

The product concerned has been manufactured under a quality management system
according to Annex II of Directive 93/42/EEC.

Compliance of the designated product with the Directive 93/42/EEC has been assured via
assessment of the quality management system by the Notified Body

TÜV Rheinland LGA Products GmbH
Tillystraße 2
90431 Nürnberg
Deutschland

Certificate No.: HD 60144232 0001

Issue date: 26.05.2020

Expiry date: 26.05.2024

following the procedure relating to the EC Declaration of Conformity set out in Annex II of
Directive 93/42/EEC.

This Declaration of conformity is valid in connection with the release document for the
respective batch of produced devices.

The above mentioned declaration of conformity is exclusively under the responsibility of

SCW MEDICATH LTD.
NO.4, Baolong 6th Road, Baolong Industrial Town,
Longgang District, Shenzhen, 518116, Guangdong, P.R. China

Shenzhen, 2020/07/16

Place, date

Miriam Xia, RA

Legally binding signature, Function

EC Declaration of Conformity
SCW-MDTE-AS-DOC A.7

PRODUCT RECOMMENDATION



Angiographic Syringes
Accessories

What else do you want to know? Contact us right now!

*Your Name

Country

*Your Email

Your Phone

*Message

*Verify Code

20K7

Submit

We'd Love to Hear From You

Curious about product features, price quote, registration, delivery time, etc? We're ready to answer any and all questions within one business day.

Contact Us

Contact

Email: info@scwmed.com

Mobile Phone: +8613691980679

Telephone: 86-755-89926266 ext 4

Fax: 86-755-89668199 ext 827

NO.4 Baolong 6th Road, Baolong industrial Town,
Longgang District, Shenzhen, P.R. China

