

MEDIVATORS™ ISA™

Endoscope Reprocessor | **Reprocess**





THE COMPLETE CIRCLE OF PROTECTION

As the global vanguard in infection prevention, Cantel Medical delivers the Complete Circle of Protection, a full-value, proactive partnership dedicated to helping you remove risk, streamline operational efficiencies and optimize your success.

Evolutionary process, Revolutionary solution

From more than 20 years, Medivators is dedicated to developing and providing innovative solutions in infection prevention and control worldwide.

As a result of this experience, the MEDIVATORS™ ISA™ Endoscope Reprocessor represents the unique new generation response. It is a rapid, efficient and flexible reprocessing solution designed for endoscopy units demanding the most advanced technology and the highest quality standards.



REPROCESS

High-level disinfection is the cornerstone of infection prevention. Reprocessing technologies from Cantel are designed to ensure patient safety by adhering to the strictest of standards for high level disinfection, optimize workflow efficiency, and deliver versatility by supporting a wide range of scopes.

MEDIVATORS™ ISA™ Endoscope Reprocessor – The Fast And Flexible Reprocessing Solution


Fast & Simple Operation for Improved Efficiency

- Full cycle in 20 minutes with channel monitoring
- Simple, user-friendly handling and operation
- Insertion of endoscopes and operation information through RFID system – helps to avoid errors and saves time




Compact & Ergonomic Design for maximum user comfort

- Large and easy to load basin accommodates all types of flexible endoscopes
- Intuitive interface assures easy and fast access to all the cycle information
- Compact size contributes to easy installation, operation and maintenance



Ensuring safety for operators, patients & endoscopes

- Hands free operation and RFID – based technology for complete infection control
- Full traceability through printed and electronic cycle documentation
- Validated process fully compliant to the standards UNI EN ISO 15883-1/4 and UNI CEN ISO/TS 15883-5 (including Austrian, German and French test soil)



The right formula for the best performance & result

- Dedicated room temperature chemicals compatible with all types of flexible endoscopes
- Improved water and chemical consumption
- Designed specifically for use with MEDIVATORS ISA Endoscope Reprocessor



Product Specifications MEDIVATORS™ ISA™ Endoscope Reprocessor



Dimensions

700 mm (l) x 1400 mm (h) x 650 mm (d) - monitor included

Weight

75 kg

Electrical Requirements

100 V-240 V, 50-60Hz

Operating Room Temperature

5° - 40° C

Humidity

20 - 80%

Reprocessing Time

20 minutes (full cycle), 12 minutes only (HLD cycle)

Process Capacity

1 endoscope per cycle

High Level Disinfectant/Sterilant

ISASPOR™ Single Shot Disinfectant/Sterilant, two parts, Peracetic acid and ISAZONE™ Activator.

RAPICIDE™ PA Single Shot Disinfectant/Sterilant, two part, Peracetic acid.

Detergent

ISACLEAN™ Detergent, multienzymatic detergent

INTERCEPT™ PLUS Detergent, Non-Enzymatic detergent

Ordering Information

| Part Number | Description | Quantity |
|-------------|---|------------------|
| 500187 | ISASPOR Single Shot High-Level Disinfectant/Sterilant | 1 x 10L Sol. A+B |
| 500188 | ISASPOR Single Shot High-Level Disinfectant/Sterilant | 1 x 5L Sol. A+B |
| ML02-0158 | RAPICIDE PA Single Shot High-Level Disinfectant/Sterilant | 1 x 5L A+B |
| ML02-0115 | RAPICIDE PA Single Shot High-Level Disinfectant/Sterilant | 2 x 5L AA+BB |
| 500189 | ISACLEAN Detergent | 2 x 5L |
| ML02-0145 | INTERCEPT PLUS Detergent | 2 x 5L |
| 93250-222 | ISA CPC Filters | 0,45µm + 0,1µm |
| ISA4400 | ISA Air Filter | 1 |
| ISAROTS2 | ISA Printer Roll | 1 |

ISAT™, ISASPOR™, ISAZONE™ and ISACLEAN™ are trademarks of Cantel Medical (Italy) S.r.l.

MEDIVATORS™, RAPICIDE™ and INTERCEPT™ PLUS are trademarks of Medivators Inc.



Cantel Medical (Italy) S.r.l.
Via Laurentina, 169
00071 Pomezia (RM) Italia
Tel.: +39 06 9145399
Fax: +39 06 9146099
www.cantelmedical.it
www.cantelmedical.eu

Medivators Inc.
14605 28th Avenue North
Minneapolis, MN
55447-4822 USA
Tel: +1.763.553.3300
Fax: +1.763.553.3387

Medivators BV
Sourethweg 11
6422 PC Heerlen
The Netherlands
Tel: +31.45.5.471.471
Fax: +31.45.5.429.595

Cantel (UK) Limited
Campfield Road
Shoeburyness, Essex
SS3 9BX,
United Kingdom
Tel: +44 (0) 1702.291878
Fax: +44 (0) 1702.290013

Cantel Medical Asia/Pacific Pte. Ltd.
1A International Business Park
#05-01 Singapore 609933
Tel: +65.6227.9698
Fax: +65.6225.6848

Cantel Medical Devices (China) Co. Ltd.
Unit 804-805, Innov Tower Block A,
Hongmei Road, Xuhui 200233 Shanghai
Tel: +86 21 60161380
Fax: +86 21 61210913