
## **Clear-Retriever**



#### (1) Soft Mesh Net

• Less possibility of damaging tissues.

#### (2) Wide Net Loop

• Enables to hold piecemeal polyps.

# Retrieval of polyp and multiple polyps







### SPECIFICATION

Model	Loop Size	Type of Retation	Loop Shape	Sheath 0.D	Working Length	Number of Use	Validity of Sterilzation
FM-ER0008	30mm	Rotational	Hexagonal	¢2.5mm	2400mm	Single Use Only	3 years

# **ClearCap Distal Attachment**



#### FEATURE & BENEFITS

#### (1) Maintaining the Clear View

• This ClearCap Distal Attachment helps maintain the clear view by keeping a moderate distance from the tissues.

### (2) Protecting the tissues

• By adopting the soft raw material, it helps users protect the

#### (3) Two draining holes

• These holes make the discharge of foreign materials and saline irrigation more easily.





#### SPECIFICATION

Model	Appilcable Lesion	0.D	I.D	Quantity of Hole for Draining	Number of Use	Validity of Sterilzation
FM-EC0001	Lower GI	¢15mm	¢12.6mm			
FM-EC0002	Upper GI	¢11.8mm	¢9.4mm			
FM-EC0003	Upper GI	¢11.35mm	¢8.9mm			
FM-EC0004	Upper GI	¢12.4mm	¢10mm	2	Single Use Only	3years
FM-EC0005	Lower GI	¢13.4mm	¢11mm		,	
FM-EC0006	Lower GI	¢14mm	¢11.6mm			
FM-EC0007	Lower GI	¢15.7mm	¢13.3mm			

Imagine What a Woderful difference You have 9