

WEGO Blood Bags



OUR COMMITMENT


The Weigao Group was established in 1988, and has since become a top medical supplies manufacturer in China with a workforce of over 10000; principally involved in the research and development ("R&D"), production and sales of branded disposable medical devices.

The Blood bag is one of Weigao Group's most experienced production lines, which has been manufactured since its foundation in 1992. Weigao Group has since introduced certified blood bags, production and quality control to market complete with CE marks, ISO9001, ISO13485 and GMP certificates.

All manufacturing processes are quality controlled and maintained to the highest standards; Weigao's production process includes the use of imported high quality welding equipment from the Netherlands and extrusion equipment from Italy. Blood management systems include all solutions such as CPDA-1, CPD+SAGM allowing platelet storage times of 3-5 days.

Wego is committed to delivering a quality service and supply for blood transfusion and blood storage products. The responsible management of blood and its components are important issues in blood transfusion and storage. Wego works with the co-operation of blood banks and donation centres in order to ensure the continued delivery of new and fresh blood to where ever it is needed; with Wego's advanced safe blood management system, more lives are saved every day.

SINGLE BLOOD BAG



Single bag is used for collection of whole blood with anticoagulant CPDA-1 solutions USP.

This single system with a primary bag containing anti-coagulant CPDA-1 (citrate-phosphate-Dextrose-Adenine) which preserves red blood cells up to 35 days at 2°C-6°C.

Sizes and Specification:




Code Number	Bag Capacity	Satellite	Needle	Anticoagulant		Al foil pack /carton
SC-100W	100ml		16G	CPDA-1	14ml	5/100
SC-250W	250ml		16G	CPDA-1	35ml	5/100
SC-350W	350ml		16G	CPDA-1	49ml	5/100
SC-450W	450ml		16G	CPDA-1	63ml	5/100
SC-500W	500ml		16G	CPDA-1	70ml	5/100

Code Number	Bag Capacity	Satellite	Needle	Anticoagulant		Al foil pack /carton
SC-100E	100ml		16G	CPDA-1	14ml	5/100
SC-250E	250ml		16G	CPDA-1	35ml	5/100
SC-350E	350ml		16G	CPDA-1	49ml	5/100
SC-450E	450ml		16G	CPDA-1	63ml	5/100
SC-500E	500ml		16G	CPDA-1	70ml	5/100



DOUBLE BLOOD BAG



This system is used for separation of two components from whole blood. This double system includes one primary bag with anticoagulant CPDA-1 Solutions USP and one empty satellite bag.


Sizes and Specification:



Code Number	Bag Capacity	Satellite	Needle	Anticoagulant		Al foil pack /carton
DC-250W	250ml	1x200ml	16G	CPDA-1	35ml	5/80
DC-350W	350ml	1x300ml	16G	CPDA-1	49ml	5/80
DC-450W	450ml	1x300ml	16G	CPDA-1	63ml	5/80
DC-500W	500ml	1x300ml	16G	CPDA-1	70ml	5/80

Code Number	Bag Capacity	Satellite	Needle	Anticoagulant		Al foil pack /carton
DC-250E	250ml	1x200ml	16G	CPDA-1	35ml	5/80
DC-350E	350ml	1x300ml	16G	CPDA-1	49ml	5/80
DC-450E	450ml	1x300ml	16G	CPDA-1	63ml	5/80
DC-500E	500ml	1x300ml	16G	CPDA-1	70ml	5/80

TRIPLE BLOOD BAG



This system is used for separation of three components from whole blood. The triple system includes one primary bag with anticoagulant CPDA-1 solution BP and two empty satellite bags.

The triple SAGM system includes one primary bag with anticoagulant CPD Solution USP, SAGM solution in another bag and one empty satellite bag.


Sizes and Specification:




Code Number	Bag Capacity	Satellite	Needle	Anticoagulant		Al foil pack /carton
TC-250W	250ml	2x200ml	16G	CPDA-1	35ml	2/50
TC-350W	350ml	2x300ml	16G	CPDA-1	49ml	2/50
TC-450W	450ml	2x300ml	16G	CPDA-1	63ml	2/50
TC-500W	500ml	2x300ml	16G	CPDA-1	70ml	2/50
TS-350W	350ml	2x300ml	16G	CPD+SAGM	49ml+80ml	2/50
TS-450W	450ml	2x300ml	16G	CPD+SAGM	63ml+ 102ml	2/50
TS-500W	500ml	2x300ml	16G	CPD+SAGM	70ml+ 114ml	2/50

Code Number	Bag Capacity	Satellite	Needle	Anticoagulant		Al foil pack /carton
TC-250E	250ml	2x200ml	16G	CPDA-1	35ml	2/50
TC-350E	350ml	2x300ml	16G	CPDA-1	49ml	2/50
TC-450E	450ml	2x300ml	16G	CPDA-1	63ml	2/50
TC-500E	500ml	2x300ml	16G	CPDA-1	70ml	2/50
TS-350E	350ml	2x300ml	16G	CPD+SAGM	49ml+80ml	2/50
TS-450E	450ml	2x300ml	16G	CPD+SAGM	63ml+ 102ml	2/50
TS-500E	500ml	2x300ml	16G	CPD+SAGM	70ml+ 114ml	2/50





QUADRUPLE BLOOD BAG



This system is used for separation of four components from whole human blood. The quadruple system includes one primary bag with anticoagulant CPDA-1 solution USP and three empty satellite bags.

The Quadruple SAGM Blood Bag system includes one primary bag with anticoagulant CPD solution USP, SAGM solution in another bag and two empty satellite bags.


Sizes and Specification:



Code Number	Bag Capacity	Satellite	Needle	Anticoagulant		Al foil pack /carton
QC-350W	350ml	3x300ml	16G	CPDA-1	49ml	2/40
QC-450W	450ml	3x300ml	16G	CPDA-1	63ml	2/40
QC-500W	500ml	3x300ml	16G	CPDA-1	70ml	2/40
QS-350W	350ml	3x300ml	16G	CPD+SAGM	49ml+80ml	2/40
QS-450W	450ml	3x300ml	16G	CPD+SAGM	63ml+102ml	2/40
QS-500W	500ml	3x300ml	16G	CPD+SAGM	70ml+114ml	2/40

Code Number	Bag Capacity	Satellite	Needle	Anticoagulant		Al foil pack /carton
QC-350E	350ml	3x300ml	16G	CPDA-1	49ml	2/40
QC-450E	450ml	3x300ml	16G	CPDA-1	63ml	2/40
QC-500E	500ml	3x300ml	16G	CPDA-1	70ml	2/40
QS-350E	350ml	3x300ml	16G	CPD+SAGM	49ml+80ml	2/40
QS-450E	450ml	3x300ml	16G	CPD+SAGM	63ml+102ml	2/40
QS-500E	500ml	3x300ml	16G	CPD+SAGM	70ml+114ml	2/40

QUALITY FEATURES



1. 16G sharp pointed Japanese needle that is siliconised with ultra-thin wall. 17G needle is also available.
2. Excellent break-off needle cover makes the needle non-reusable.
3. Provided with standard donor tubing and code number on the tube surface.
4. Tamper proof, safe and easy-to-open Port covers are provided to avoid contamination.
5. Rounded shape of the bag minimizes loss of blood components during Transfer and Transfusion.
6. Convenient hanger slits and holes are provided for use during Blood Collection and Blood Transfusion. This also allows easy suspension of bag in a vertical position.
7. Particle-free and high quality transparent Medical grade PVC sheet is provided for accurate and easy blood monitoring during Collection, Transfer and Transfusion.

AVAILABLE OPTIONS

1. Blood bag types available : CPDA -1 / CPD / SAGM.
2. With Safety Needle Shield.
3. With Sampling bag and Vacuum Blood Collection Tube Holder.
4. High quality film suitable for extended storage of viable platelets for approximately 5 days.
5. Blood bag with leukoreduction filter.
6. Transfer empty bag is also available from 150ml to 2000ml for separating blood components from whole blood.





WEGO 威高

YOUR HEALTH OUR CARE

Weigao Group Corporation

Head Office:

312 Shichang Road
Weihai, Shandong
264209

P.R. of China

Email: info@weigagroup.com
exp@weigagroup.com

Tel: +86 631 5660699

Fax: +86 631 5622563

Wego Medical Europe Limited

Email: info@weigao.co.uk

Tel: 01702 328398

Fax: 01702 328387

Registered Office

Colechurch House, 1 London Bridge Walk
London SE1 2SX

Registered in England Number: 5878396