






AMPLITUDE[®]

COMETE[®]
ACL femoral cortical fixation

COMETE® : Device features



- Titanium plate TA6VELI
- Polyethylene terephthalate loop



- Traction wire (polyethylene terephthalate) Ø 0.85 mm, rigid and fixed tip facilitating introduction into the eye



- Cannulated reamer Ø 5 mm



- Graduated drill pin with eye Ø 2.4 mm, 300 mm length

COMETE® : implants instrumentation

CE 0120
CE 0124



Making by Cousin Biotech

Making by Amplitude

COMETE® ACL femoral cortical fixation 15 mm length
COMETE® ACL femoral cortical fixation 20 mm length
COMETE® ACL femoral cortical fixation 25 mm length
COMETE® ACL femoral cortical fixation 30 mm length
COMETE® ACL femoral cortical fixation 35 mm length
COMETE® ACL femoral cortical fixation 40 mm length

OAMGEFX15U
OAMGEFX20U
OAMGEFX25U
OAMGEFX30U
OAMGEFX35U
OAMGEFX40U

Cannulated reamer Ø 5 mm
Graduated drill pin with eye Ø 2.4 mm, 300 mm length

2-0405201
2-0404700

AMPLITUDE

Customer Service France : Porte du Grand Lyon. 01700 Neyron. France. Tél. : +33 (0)4 37 85 19 19. Fax : +33 (0)4 37 85 19 18
Customer Service Export : 11, cours Jacques Offenbach. Zone Mozart 2. 26000 Valence. France. Tél. : +33 (0)4 75 41 87 41. Fax : +33 (0)4 75 41 87 42
E-mail : amplitude@amplitude-ortho.com. Internet : www.amplitude-ortho.com