

OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION
SALMONELLA SELECTIVE SUPPLEMENT SR0194E
SALMONELLA SELECTIVE SUPPLEMENT
SR0194E
Formula

Vial contents (each vial is sufficient to supplement 500ml of medium)

Novobiocin	2.5 mg
Cefsulodin	6.0 mg

Description

 A selective supplement for the isolation of *Salmonella* spp.

Directions

Aseptically add 2ml of distilled water to one vial to dissolve. Add the vial contents to 500ml of *Brilliance*[™] Salmonella Agar Base (CM1092). Mix well and bring to the boil with frequent agitation. DO NOT AUTOCLAVE. Cool to 50°C, mix well and pour into sterile Petri dishes.

Physical Characteristics

White or pale yellow pellet

Microbiological Tests Using Optimum Inoculum Dilution

Control Medium: Tryptone Soya Agar

Reactions after incubation at 37°C for 24 ± 2 hours

 Tested in *Brilliance*[™] Salmonella Agar Base CM1092

Medium is challenged with 10-100 colony-forming units

<i>Salmonella typhimurium</i>	ATCC®14028	1-2mm purple/pink colonies
<i>Salmonella enteritidis</i>	ATCC®13076	1-2mm purple/pink colonies
<i>Salmonella arizonae</i>	ATCC®13314	1-4mm purple/pink colonies

A satisfactory result is represented by recovery of positive strains equal to or greater than 70% of the control medium.

	Document Owner Department: QC	BT-SPEC-0531
		Page 2 of 3
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
SALMONELLA SELECTIVE SUPPLEMENT SR0194E		

Medium is challenged with 1E+02 to 1E+04 colony-forming units

<i>Escherichia coli</i>	ATCC®25922	No growth or 1-2mm cream colonies
-------------------------	------------	-----------------------------------

Medium is challenged with greater than 1E+04 colony-forming units

<i>Pseudomonas aeruginosa</i>	ATCC®27853	No growth or 1-2mm pink colonies
-------------------------------	------------	----------------------------------

Negative strains are inhibited or shall produce at least a 2 log (10) reduction when compared to the control medium.

Medium is challenged with 1E+04 to 1E+06 colony-forming units

<i>Klebsiella pneumoniae</i>	ATCC®29665	1-3mm mucoid, blue colonies
<i>Enterococcus faecalis</i>	ATCC®29212	No growth
<i>Proteus mirabilis</i>	NCTC10975	No growth or ppt-0.5mm straw colonies

Negative strains are inhibited or shall produce a negative diagnostic reaction.

	Document Owner Department: QC	BT-SPEC-0531
		Page 3 of 3
OXOID QUALITY ASSURANCE PRODUCT SPECIFICATION		
SALMONELLA SELECTIVE SUPPLEMENT SR0194E		

Revision History

Section / Step	Description of Change	Reason for Change	Reference
<i>Proteus mirabilis</i>	Update to testing specification	Change control	BT-CC-1859
Entire document	Correction of typographical/minor errors. Addition of control media + result criteria.	Change control	BT-CC-1918
Directions	Remove 'completely'	Change control	BT-CC-2725