



Cardiac Markers



ELISA Kits



Infectious Diseases



Cancer Markers



Drugs of Abuse

LumiQuick

DIAGNOSTICS, INC.

In Vitro Diagnostic Devices

CATALOG

www.LumiQuick.com



www.youtube.com/lumiquick





ISO 13485

cGMP



Introduction

LumiQuick Diagnostics Inc., located in the heart of Silicon Valley, California, USA, develops, manufactures and markets high quality point of care tests and other immunoassay kits for the world-wide in-vitro diagnostic market. Through years of continuous and diligent efforts, we are proud to offer extensive product lines of rapid tests, ELISA, CLIA and raw materials. These product lines provide the immunoassays in various formats to detect infectious diseases, cardiac markers, cancer markers, drugs of abuse, fertility hormones and other disease markers. Many of the products have been CE marked and registered in different countries.

Our mission is to provide the high quality products with affordable cost to help people around the world to fight human diseases and illicit substance abuse. We are committed to providing exceptional quality, customer service and flexibility to meet our customer's need.

We appreciate very much your interests in our products and look forward to be your partner in IVD business!




Table of Contents

Infectious Diseases	3 - 5
Drugs of Abuse	6 - 7
Cardiac Markers	8
Fertility Hormones	9
Cancer Markers	10
Enzyme Immunoassay Tests	11 - 12
Urine Test Strip	13
Raw Materials	13
Contract Services	13
Product Highlights	14 - 25
Product Specifications and Order Information	26


Infectious Diseases

Human Immunodeficiency Virus (HIV)



Cat. No.	Product Name	No. per Kit
71001	HIV I&II Test Strip	50
71002	HIV I&II Test Card	25
71068	HIV I & II Tri-line Test Card	25

Hepatitis B Virus (HBV)



Cat. No.	Product Name	No. per Kit
71003	HBsAg Test Strip	50
71004	HBsAg Test Card	25
71005	HBsAb Test Strip	50
71006	HBsAb Test Card	25
71007	HBcAb Test Strip	50
71008	HBcAb Test Card	25
71009	HBeAg Test Strip	50
71010	HBeAg Test Card	25
71011	HBeAb Test Strip	50
71012	HBeAb Test Card	25
71013-1	HBV-3 Panel Test Card (HBsAg/HBeAg/HBcAb)	25
71013-2	HBV-3 Panel Test Card (HBsAg/HBsAb/HBcAb)	25
71014	HBV-5 Panel Test Card (HBsAg/HBsAb/HBeAg/HBeAb/HBcAb)	25
71069-1	HBV-2 Panel Test Card (HBsAg/HBsAb)	25

Tropical Diseases



Cat. No.	Product Name	No. per Kit
71019	Dengue IgG/IgM Combo Test Card	25
71028	Salmonella Typhi Antigen Test Card	25
71047	Dengue IgG/IgM Combo Test Strip	50
71049	Malaria pf Antigen Test Card	25
71050	Malaria pf/pv Antigen Test Card	25
71053	Malaria pan Antigen Test Card	25
71054	Salmonella Typhi/Paratyphi Antigen Test Card	25
71055	Salmonella Typhi IgG/IgM Duo Test Card	25
71066	Dengue NS1 Antigen Test Card*	25

* Products that are under development. Please contact us for the availability.

Infectious Diseases

Hepatitis C Virus (HCV)




Cat. No.	Product Name	No. per Kit
71027	HCV Antibody Test Strip	50
71030	HCV Antibody Test Card	25

Gastro-Intestinal Diseases



Cat. No.	Product Name	No. per Kit
71020	H. Pylori Stool Antigen Test Card	25
71024	H. Pylori Antibody Test, Serum/Plasma/Whole Blood	25
71029	Rotavirus Stool Antigen Test Card	25
71032	Adenovirus Stool Antigen Test Card	25
71033	Adeno/Rota Virus Stool antigen Combo Test Card	25
71046	H. Pylori Antibody Test Card, Serum/Plasma	25
71051	Adeno/Rota Virus Stool antigen Combo Test Strip	50
71064	H. Pylori Antibody Test Strip, Serum/Plasma	50

Respiratory Diseases



Cat. No.	Product Name	No. per Kit
71034	Legionella Pneumonia Urinary Antigen Test Card	25
71035	Legionella Pneumonia Urinary Antigen Test Strip	50
71048	Strep A Test	25
71062	Tuberculosis Antibody Test Card	25
71065	Influenza A Test Card (Flu A)	25
71067	Influenza A+ Influenza B Test Card (Flu A+B)	25

Infectious Diseases

ToRCH



Cat. No.	Product Name	No. per Kit
71036	Toxo IgM Test*	25
71037	Rubella IgG Test*	25
71038	Rubella IgM Test*	25
71039	CMV IgG Test*	25
71040	CMV IgM Test*	25
71041	HSV I IgG Test*	25
71042	HSV I IgM Test*	25
71043	HSV II IgG Test*	25
71044	HSV II IgM Test*	25
71075	Toxo IgG Test*	25
71076	Toxo IgG/IgM Test	25

Sexually Transmitted Diseases (STD)



Cat. No.	Product Name	No. per Kit
71015	Syphilis Antibody Test Strip, Serum/Plasma	50
71016	Syphilis Antibody Test Card, Serum/Plasma/Whole Blood	25
71045	Chlamydia Test Card*	25
71057	Syphilis Antibody Test Card, Serum/Plasma	25

Other Infectious Diseases





Cat. No.	Product Name	No. per Kit
71017	Rickettsia IgG/IgM Test Card	25
71018	Hantavirus IgG/IgM (HFRS) Test Card	25

* Products that are under development. Please contact us for the availability.

Drugs of Abuse

Single



Cat. No.	Product Name	No. per Kit
74001	Saliva Alcohol Test	50
74002	Tramadol Test Card	25
74003	Tramadol Test Strip	50
74013	Amphetamine Strip	50
74014	Amphetamine Test Card	25
74015	Barbiturate Test Strip	50
74016	Barbiturate Test Card	25
74017	Benzodiazepine Test Strip	50
74018	Benzodiazepine Test Card	25
74019	Cocaine Test Strip	50
74020	Cocaine Test Card	25
74021	EDDP Test Strip	50
74022	EDDP Test Card	25
74023	MDMA Test Strip	50
74024	MDMA Test Card	25
74025	Methadone Test Strip	50
74026	Methadone Test Card	25
74027	Methamphetamine Test Strip	50
74028	Methamphetamine Test Card	25
74029	Opiate Test Strip, 300 ng/ml	50
74030	Opiate Test Card, 300 ng/ml	25
74031	Opiate Test Strip, 2000 ng/ml	50
74032	Opiate Test Card, 2000 ng/ml	25
74033	Phencyclidine Test Strip	50
74034	Phencyclidine Test Card	25
74035	Marijuana (THC) Test Strip	50
74036	Marijuana (THC) Test Card	25
74037	TCA Test Strip	50
74038	TCA Test Card	25
74039	Ketamine Test Strip	50
74040	Ketamine Test Card	25
74041	Buprenorphine Test Strip	50
74042	Buprenorphine Test Card	25
74043	Oxycodone Test Strip	50
74044	Oxycodone Test Card	25
74049	Urine Alcohol Test	50
74051	Propoxyphene Test Strip	50
74052	Propoxyphene Test Card	25
74055	Methylphenidate (Ritalin) Test Strip	50
74056	Methylphenidate (Ritalin) Test Card	25
74057	Fentanyl Test Strip	50
74058	Fentanyl Test Card	25
74059	Clonazepam Test Strip	50
74060	Clonazepam Test Card	25
74061	Cotinine Test Strip	50
74062	Cotinine Test Card	25
74063	Clenbuterol Test Strip	50
74064	Clenbuterol Test Card	25

Drugs of Abuse

Multiple Drug Panel Test

Cat. No.	Product Name	No. per Kit
74004	DOA-2 Panel Test, Dip Device	25
74004-TC	DOA-2 Panel Test Card, Cassette Device	25
74005	DOA-3 Panel Test, Dip Device	25
74005-TC	DOA-3 Panel Test Card, Cassette Device	25
74006	DOA-4 Panel Test, Dip Device	25
74006-TC	DOA-4 Panel Test Card, Cassette Device	25
74007	DOA-5 Panel Test, Dip Device	25
74007-TC	DOA-5 Panel Test Card, Cassette Device	25
74008	DOA-6 Panel Test, Dip Device	25
74008-TC	DOA-6 Panel Test Card, Cassette Device	25
74009	DOA-7 Panel Test, Dip Device	25
74009-TC	DOA-7 Panel Test Card, Cassette Device	25
74010	DOA-8 Panel Test, Dip Device	25
74010-TC	DOA-8 Panel Test Card, Cassette Device	25
74011	DOA-9 Panel Test, Dip Device	25
74011-TC	DOA-9 Panel Test Card, Cassette Device	25
74012	DOA-10 Panel Test, Dip Device	25
74012-TC	DOA-10 Panel Test Card, Cassette Device	25
74045-A	Methamphetamine/Marijuana (THC) Test	25
74045-B	Methamphetamine/Opiates Test	25
74046	DOA Panel Test Cup, up to 12 drugs	25
74053	DOA-11 Panel Test, Dip Device	25
74053-TC	DOA-11 Panel Test Card, Cassette Device	25
74054	DOA-12 Panel Test, Dip Device	25
74054-TC	DOA-12 Panel Test Card, Cassette Device	25




Cardiac Markers

Panel



Cat. No.	Product Name	No. per Kit
75003	Cardiac Panel Test Card, Serum/Plasma	25
75004	Cardiac Panel Test Card, Serum/Plasma/Whole Blood	25

Troponin I



Cat. No.	Product Name	No. per Kit
75001	Troponin I Test Card, Serum/Plasma	25
75002	Troponin I Test Card, Serum/Plasma/Whole Blood	25

Myoglobin



Cat. No.	Product Name	No. per Kit
75005	Myoglobin Test Card, Serum/Plasma	25
75006	Myoglobin Test Card, Serum/Plasma/Whole Blood	25

CK-MB



Cat. No.	Product Name	No. per Kit
75007	CK-MB Test Card, Serum/Plasma	25
75008	CK-MB Test Card, Serum/Plasma/Whole Blood	25

Fertility Hormones

HCG



Cat. No.	Product Name	No. per Kit
73001	HCG Midstream Test	25
73002	HCG Pregnancy Urine Test Card	25
73005	HCG Pregnancy Serum/Urine Test Card	25
73006	HCG Pregnancy Urine Test Strip	50
73011	HCG Serum/Urine Test Card	50

LH



Cat. No.	Product Name	No. per Kit
73003	LH Ovulation Midstream Test	25
73004	LH Test Card	25
73007	LH Test Strip	50

TSH



Cat. No.	Product Name	No. per Kit
73008	TSH Test Strip*	50
73009	TSH Test Card*	25

* Products that are under development. Please contact us for the availability.

Cancer Markers

FOB



Cat. No.	Product Name	No. per Kit
72001	Fecal Occult Blood Test Card	25
72006	Fecal Occult Blood Test Strip	50

PSA



Cat. No.	Product Name	No. per Kit
72003	PSA Test Card	25
72004	PSA Whole Blood Test Card	25
72005	PSA Test Strip	50

AFP



Cat. No.	Product Name	No. per Kit
72002	AFP Test Card	25
72007	AFP Test Strip	50

CEA



Cat. No.	Product Name	No. per Kit
72009	CEA Test Card	25
72010	CEA Test Strip	50

Enzyme Immunoassay Tests

ELISA Kits



Cat. No.	Product Name	No. per Kit
78001	HAV-IgM	96
78002	HAV-IgG	96
78003	HBsAg	96
78004	Anti-HBs	96
78005	HBeAg	96
78006	Anti-HBe	96
78007	Anti-HBc	96
78008	Anti-HBc-IgM	96
78009	HCV	96
78010	HIV	96
78011	HEV-IgM	96
78012	HEV-IgG	96
78014	HTLV	96
78015	Treponema Pallidum	96
78016	AFP	96
78017	CEA	96
78018	CA 125	96
78019	CA 15-3	96
78020	CA 19-9	96
78021	PSA	96
78023	Free Beta-HCG	96
78025	TSH	96
78026	Ultra-sensitive TSH	96
78027	T3	96
78028	Free T3	96
78029	T4	96
78030	Free T4	96
78031	Rubella IgG	96
78032	Rubella IgM	96
78033	Toxoplasma IgG	96
78034	Toxoplasma IgM	96
78035	CMV IgG	96
78036	CMV IgM	96
78037	HSV I IgG	96
78038	HSV I IgM	96
78039	HSV II IgG	96
78040	HSV II IgM	96
78041	H. Pylori IgG Qualitative	96
78042	H. Pylori IgG Quantitative	96
78043	H. Pylori IgM	96

Enzyme Immunoassay Tests

ELISA Kits

Cat. No.	Product Name	No. per Kit
78044	HCG	96
78045	FSH	96
78046	LH	96
78047	Prolactin	96
78048	Visual hCG Pregnancy	96
78049	Ferritin	96
78050	HGH	96
78051	Ig E	96
78052	Estradiol	96
78053	Progesterone	96
78054	Testosterone	96
78055	Troponin I	96
78056	Myoglobin	96
78057	High sensitivity-CRP	96
78058	Free-PSA	96
78059	HDV-IgG	96
78060	HDV-IgM	96
78061	HDV Ag	96



Urine Test Strips

Urine Test Strip



Cat. No.	Product Name	No. per Kit
76001	10 Parameter Test Strip	100
76002	7 Parameter Adulteration Test Strip	25

Raw Materials

Raw Materials

Well paired materials and reagents for infectious disease, drugs of abuse, cancer markers and cardiac marker tests

Contract Services

Contract Services

We provide contract services for feasibility study, product development and manufacturing.

Product Highlights

Product Highlights

This section contains only selected products' flyers of LumiQuick Diagnostics, Inc. Please contact us or our local distributor if you are interested in any of the flyers for your marketing program.

COLORECTAL CANCER

Fecal Occult Blood Rapid Test

QuickView™ Fecal Occult Blood Test is a qualitative test that detects human hemoglobin in human fecal specimens. The test is a visual one step, in vitro assay. It is intended for professional use to help diagnose gastrointestinal bleeding.



IVD

LumiQuick Diagnostics, Inc. Tel: 1-408-855-0061 Fax: 1-408-855-0062
Email: info@lumiquick.com www.lumiquick.com

ACQUIRED IMMUNODEFICIENCY SYNDROME (AIDS)

HIV Type 1/2 Rapid Test

QuickView™ HIV 1/2 Test Card is a single use, rapid test device, double antigen sandwich immunoassay for qualitative detection of antibodies to HIV-1 and HIV-2 in serum, plasma or whole blood. It is intended for medical institutions as the aid for the diagnosis of possible HIV infection. This product may be used for the quick screening of blood donors and blood products.



IVD

LumiQuick Diagnostics, Inc. Tel: 1-408-855-0061 Fax: 1-408-855-0062
Email: info@lumiquick.com www.lumiquick.com

Hepatitis B virus Infection

Hepatitis B virus Antigen Rapid Test

QuickView™ HBsAg Test Card is a double antibody sandwich immunoassay. Colloidal gold conjugated anti-HBsAg antibody complexes are dry-immobilized in the test device. When the sample is added, it migrates by capillary diffusion through the strip rehydrating the gold conjugate complexes. If present, HBsAg will react with the gold conjugate complexes forming particles. These particles will continue to migrate along the strip until the Test Zone (T) where they are captured by anti-HBsAg antibodies immobilized there and a visible red line appears. If there is no HBsAg in sample, no red line will appear in the test zone (T). The gold conjugate complexes will continue to migrate along until they are captured in the Control Zone (C) by immobilized goat anti-mouse IgG antibody aggregating a red line, which indicates the validity of the test.



IVD

LumiQuick Diagnostics, Inc. Tel: 1-408-855-0061 Fax: 1-408-855-0062
Email: info@lumiquick.com www.lumiquick.com

H. PYLORI INFECTION

H. Pylori Antibody Rapid Test

QuickView™ H. Pylori Antibody Test is a rapid immunochromatographic assay for the detection of antibodies to H. Pylori in human whole blood, serum or plasma. The assay is to provide an aid in diagnosis of H. pylori infection.



IVD

LumiQuick Diagnostics, Inc. Tel: 1-408-855-0061 Fax: 1-408-855-0062
Email: info@lumiquick.com www.lumiquick.com

H. PYLORI INFECTION

H. Pylori Antigen Rapid Test

QuickView™ H. pylori Antigen Test is an in vitro qualitative immunochromatographic assay for the rapid detection of Helicobacter pylori antigens in human stool specimen. The test results are intended to aid in the diagnosis of H. pylori infection, to monitor the effectiveness of therapeutic treatment and to confirm the eradication of H. pylori in peptic ulcer patients.



IVD

LumiQuick Diagnostics, Inc. Tel: 1-408-855-0061 Fax: 1-408-855-0062
Email: info@lumiquick.com www.lumiquick.com

TYPHOID FEVER

Salmonella typhi IgG/IgM Antibody Duo Rapid Test

QuickView™ Salmonella IgG/IgM Antibody Test is an in vitro qualitative immunochromatographic assay for the rapid detection of anti-S. Typhi antibodies in human whole blood/plasma/serum specimens. The test results are intended to help in the diagnosis of S. Typhi infection.



IVD

LumiQuick Diagnostics, Inc. Tel: 1-408-855-0061 Fax: 1-408-855-0062
Email: info@lumiquick.com www.lumiquick.com

GROUP A STREPTOCOCCAL (GAS) DISEASE

Strep A Rapid Test

Group A Streptococcus is a bacterium often found in the throat and on the skin. People may carry group A streptococci in the throat or on the skin and have no symptoms of disease. Most GAS infections are relatively mild illnesses such as "strep throat" or impetigo. Occasionally these bacteria can cause severe and even life-threatening diseases.

QuickView™ Strep A Test is a qualitative test that detects group A streptococcal antigen. The test is a rapid, in vitro assay. It is intended for professional use to help diagnose Strep A infection in hospital and clinical laboratories.



IVD

LumiQuick Diagnostics, Inc. Tel: 1-408-855-0061 Fax: 1-408-855-0062
Email: info@lumiquick.com www.lumiquick.com

TYPHOID FEVER

Salmonella Typhi Antigen Rapid Test

QuickView™ S. Typhi Antigen Test is an in vitro qualitative immunochromatographic assay for the rapid detection of S. Typhi antigens in human stool or blood specimen. The test results are intended to help the diagnosis of S. Typhi infection and to monitor the effectiveness of therapeutic treatment.



IVD

LumiQuick Diagnostics, Inc. Tel: 1-408-855-0061 Fax: 1-408-855-0062
Email: info@lumiquick.com www.lumiquick.com

DENGUE FEVER

Dengue IgG/IgM Rapid Test

QuickView™ Dengue Test is a rapid immunochromatographic assay for the simultaneous detection of IgG and IgM antibodies to Dengue virus in human whole blood, serum or plasma. The assay is used as a screening test for Dengue viral infection and as an aid for differential diagnosis of primary and secondary infections in conjunction with other criteria.



IVD

LumiQuick Diagnostics, Inc. Tel: 1-408-855-0061 Fax: 1-408-855-0062
Email: info@lumiquick.com www.lumiquick.com

Product Highlights

Product Highlights

This section contains only selected products' flyers of LumiQuick Diagnostics, Inc. Please contact us or our local distributor if you are interested in any of the flyers for your marketing program.

SYPHILIS

Syphilis (T. Pallidum) Antibody Test

QuickView™ Syphilis Antibody Test is a rapid immunochromatographic assay for the rapid qualitative determination of antibodies to *Trichomonas pallidum* in human whole blood, serum or plasma. The assay is used as a screening test for *T. pallidum* infection (also known as Syphilis).



IVD

LumiQuick Diagnostics, Inc. Tel: 1-408-855-0061 Fax: 1-408-855-0062
Email: info@lumiquick.com www.lumiquick.com

MALARIA

Malaria pan-LDH Antigen Test

QuickView™ Malaria pan-LDH Antigen Test is a rapid immunochromatographic assay for the rapid qualitative determination of one or more of the known Malaria species: *P. falciparum*, *P. vivax*, *P. ovale*, and/or *P. malariae* by detecting lactate dehydrogenase (LDH) in human blood. This test is an aid in the diagnosis of Malaria infection.



IVD

LumiQuick Diagnostics, Inc. Tel: 1-408-855-0061 Fax: 1-408-855-0062
Email: info@lumiquick.com www.lumiquick.com

INFLUENZA

Influenza A / B Rapid Test

QuickView™ Flu Rapid Test Devices are in vitro diagnostic tests for the qualitative detection of influenza type A and B nucleoprotein antigens in nasopharyngeal swabs and nasal aspirate specimens, using immunochromatographic method. The identification is based on the monoclonal antibodies specific for the nucleoprotein of either influenza virus A or B. They provide information to help clinical professionals to diagnose and treat influenza infection.



IVD

LumiQuick Diagnostics, Inc. Tel: 1-408-855-0061 Fax: 1-408-855-0062
Email: info@lumiquick.com www.lumiquick.com

DRUGS OF ABUSE

Drugs of Abuse Tests

Available 18 drugs in a variety of formats:

- Alcohol
- Benzodiazepines
- Buprenorphine
- COCAINE
- Ecstasy
- Marijuana
- MDMA
- Opioids
- Phencyclidine
- Propoxyphene
- THC
- Tramadol
- Amphetamines
- Marijuana
- Ecstasy
- Marijuana
- Ecstasy
- Marijuana
- Ecstasy



IVD

LumiQuick Diagnostics, Inc. Tel: 1-408-855-0061 Fax: 1-408-855-0062
Email: info@lumiquick.com www.lumiquick.com

DRUG OF ABUSE

Methamphetamine / Marijuana (THC) Test

QuickView™ Methamphetamine / Marijuana (THC) Test Card is an immunochromatographic based one step in vitro test. It is designed for qualitative determination of Methamphetamine, THC, and their metabolites in human urine specimens. This assay may be used in the point of care setting.



IVD

LumiQuick Diagnostics, Inc. Tel: 1-408-855-0061 Fax: 1-408-855-0062
Email: info@lumiquick.com www.lumiquick.com

CARCINOEMBRYONIC ANTIGEN

CEA Rapid Test

QuickView™ CEA Test is a rapid immunochromatographic assay for the qualitative detection of carcinoembryonic antigen (CEA) in human whole blood, serum or plasma. This device is intended as an aid in monitoring patients for disease progress or response to therapy or for the detection of recurrent or residual disease.



IVD

LumiQuick Diagnostics, Inc. Tel: 1-408-855-0061 Fax: 1-408-855-0062
Email: info@lumiquick.com www.lumiquick.com

ACUTE MYOCARDIAL INFARCTION

Cardiac Panel Rapid Test

QuickView™ Cardiac Panel Test is an immunochromatographic based one step in vitro test. It is designed for qualitative determination of cardiac troponin (cTnI), CK-MB and Myoglobin in human serum, plasma or whole blood specimens as an aid in the diagnosis of myocardial infarction.



IVD

LumiQuick Diagnostics, Inc. Tel: 1-408-855-0061 Fax: 1-408-855-0062
Email: info@lumiquick.com www.lumiquick.com

ACUTE MYOCARDIAL INFARCTION

Troponin I Rapid Test

QuickView™ Troponin I Test is an immunochromatographic based one step in vitro test. It is designed for qualitative determination of cardiac troponin (cTnI) in human serum, whole blood or plasma specimens as an aid in the diagnosis of myocardial infarction.



IVD

LumiQuick Diagnostics, Inc. Tel: 1-408-855-0061 Fax: 1-408-855-0062
Email: info@lumiquick.com www.lumiquick.com

PREGNANCY

HCG Pregnancy Rapid Test

QuickView™ Pregnancy HCG Test Card is a rapid test to detect the presence of hCG in urine specimens in a qualitative format sensitive to 25 mIU/mL. The test utilizes monoclonal antibodies to selectively detect elevated levels of hCG in urine. The immunological specificity of the test kit virtually eliminates cross reactivity interference from the structurally related glycoprotein hormones hLH, hH and hTSH in physiological levels.



IVD


LumiQuick Diagnostics, Inc. Tel: 1-408-855-0061 Fax: 1-408-855-0062
Email: info@lumiquick.com www.lumiquick.com

ACQUIRED IMMUNODEFICIENCY SYNDROME (AIDS)



HIV Type 1/2 Rapid Test

QuickView™ HIV I&II Test Card is a single use, rapid test device, double antigen sandwich immunoassay for qualitative detection of antibodies Immunodeficiency Viruses or HIV I&II in serum, plasma or whole blood. It is intended for medical institution as the aid for the diagnosis of possible HIV infection. The product may be used for the quick screening of blood donors and blood products.



IVD



LumiQuick Diagnostics, Inc.

2946 Scott Blvd. Santa Clara, CA 95054 USA

Tel: 1-408-855-0061 Fax: 1-408-855-0063


Email: info@LumiQuick.com www.LumiQuick.com



Preparation



Bring the kit components to room temperature before testing.

Remove the test card from the sealed foil pouch.



Test Operation


Use micropipetter to transfer 150 µL of sample, or place the transfer pipette supplied with the device in the specimen and depress the bulb to withdraw a sample.




Hold the pipette in a vertical position over the sample well of the test card and deliver 3 drops (120-150 µL) of sample into the sample well.



Read the result at 20 minutes.



Test Result Appearance and Interpretation



Ordering Information

Catalog No. **Description**

71002 QuickView™ HIV Type I & II Test Card 25 tests

More infectious disease tests are available. Please contact us or your agent for details.



LumiQuick Diagnostics, Inc.

2946 Scott Blvd. Santa Clara, CA 95054 USA

Tel: 1-408-855-0061 Fax: 1-408-855-0063

Email: info@LumiQuick.com

www.LumiQuick.com

DENGUE FEVER



Dengue IgG/IgM Rapid Test

QuickView™ Dengue Test is a rapid immunochromatographic assay for the simultaneous detection of IgG and IgM antibodies to Dengue virus in human whole blood, serum or plasma. The assay is used as a screening test for Dengue viral infection and as an aid for differential diagnosis of primary and secondary infections in conjunction with other criteria.



IVD



LumiQuick Diagnostics, Inc.

2946 Scott Blvd. Santa Clara, CA 95054 USA

Tel: 1-408-855-0061 Fax: 1-408-855-0063

Email: Info@LumiQuick.com www.LumiQuick.com




Preparation

Bring the kit components to room temperature before testing.



Test Operation


Apply 5µL of serum, plasma, or whole blood to the "S1" area indicated by an arrow mark





Add 2 drops of sample buffer to well marked as "S"



At the end of 20 minutes read the results.



Test Result Appearance and Interpretation



Ordering Information

Catalog No.	Description	
71019	QuickView™ Dengue IgG/IgM Combo Test Card	25 tests
71047	QuickView™ Dengue IgG/IgM Combo Test Strip	50 tests



LumiQuick Diagnostics, Inc.

2946 Scott Blvd. Santa Clara, CA 95054 USA

Tel: 1-408-855-0061 Fax: 1-408-855-0063

Email: info@LumiQuick.com


www.LumiQuick.com

H. PYLORI INFECTION



H. Pylori Antigen Rapid Test

QuickView™ H. pylori Antigen Test is an in vitro qualitative immuno-chromatographic assay for the rapid detection of Helicobacter pylori antigens in human stool specimen. The test results are intended to aid in the diagnosis of H. pylori infection, to monitor the effectiveness of therapeutic treatment and to confirm the eradication of H. pylori in peptic ulcer patients.



IVD



LumiQuick Diagnostics, Inc.

2946 Scott Blvd. Santa Clara, CA 95054 USA

Tel: 1-408-855-0061 Fax: 1-408-855-0063

Email: info@LumiQuick.com www.LumiQuick.com



Preparation

Bring the kit components to room temperature before testing.

Stool specimens can be collected at any time of the day.

Collect a random sample of feces in a clean, dry receptacle.



Test Operation

Unscrew the cap of the collection tube and remove the applicator stick.

Insert the stick into the fecal specimen at several different sites.



Insert the sampled applicator back to the tube and tighten the cap securely.


Shake the tubes vigorously for about 5 seconds to release and disperse the stool sample into the collection buffer.



Hold the tube upright with tip pointed toward the direction away from the test performer,
Snap off the tip.




Test Result Appearance and Interpretation



Hold the tube in a vertical position over the sample well of the test card and deliver 3 drops (120-150 μ L) of sample into the sample well marked as "S" on the cassette.



Read the result within 10 to 15 minutes



Ordering Information

Catalog No.	Description
71020	QuickView™ H. Pylori Antigen Test 25 tests

More infectious disease tests are available. Please contact us or your agent for details.



LumiQuick Diagnostics, Inc.

2946 Scott Blvd. Santa Clara, CA 95054 USA

Tel: 1-408-855-0061 Fax: 1-408-855-0063

Email: info@LumiQuick.com

www.LumiQuick.com

TYPHOID FEVER



Salmonella typhi IgG/IgM Antibody Duo Rapid Test

QuickView™ Salmonella IgG/IgM Antibody Test is an in vitro qualitative immunochromatographic assay for the rapid detection of anti-S.Typhi antibodies in human whole blood/plasma/serum specimen.

The test results are intended to help in the diagnosis of S. Typhi infection.



IVD



LumiQuick Diagnostics, Inc.

2946 Scott Blvd. Santa Clara, CA 95054 USA

Tel: 1-408-855-0061 Fax: 1-408-855-0063

Email: Info@LumiQuick.com www.LumiQuick.com



Movies: www.youtube.com/LumiQuick





Preparation

Bring the kit components to room temperature before testing.
Remove the test card from the sealed foil pouch.



Test Operation

Apply 10 µL of serum, plasma, or whole blood to each sample well.





Add 2 drops of sample buffer to each sample well, marked as "S".



Read the result at 20 minutes.



Test Result Appearance and Interpretation



Ordering Information

Catalog No.	Description	
71055	QuickView™ Salmonella typhi IgG/IgM Antibody Duo Test Card	25 tests
71028	QuickView™ Salmonella typhi Antigen Test Card	25 tests



LumiQuick Diagnostics, Inc.

2946 Scott Blvd. Santa Clara, CA 95054 USA

Tel: 1-408-855-0061 Fax: 1-408-855-0063

Email: info@LumiQuick.com www.LumiQuick.com

NEW YouTube Movies: www.youtube.com/LumiQuick


DRUGS OF ABUSE



Drugs of Abuse Tests

Available 18 drugs in a variety of formats:

Alcohol	Amphetamine
Barbiturate	Benzodiazepine
Buprenorphine	Cocaine
EDDP	Ketamine
Methadone	Methamphetamine (Ecstasy)
MDMA	Opiate
Oxycodone	Phencyclidine
Cannabinoid (THC)	Propoxyphene
Tramadol	Tricyclic antidepressant (TCA)



IVD



LumiQuick Diagnostics, Inc.

2946 Scott Blvd. Santa Clara, CA 95054 USA

Tel: 1-408-855-0061 Fax: 1-408-855-0063

Email: Info@LumiQuick.com www.LumiQuick.com

NEW

You Tube

Movies: www.youtube.com/LumiQuick



Preparation

CUP





DIP DEVICES



CARD DEVICES



Test Operation



Test Result Appearance and Interpretation



Ordering Information

Description

Quick View™ Single Drug Strip or Card Test
Quick Profile™ Multiple Drug Test Dip Device or Card Device
LumiCup™ Multiple Drug Test Cup
Quick Check™ Adulteration Strip Test

Please contact us or your agent for details.

 **LumiQuick Diagnostics, Inc.**
2946 Scott Blvd. Santa Clara, CA 95054 USA
Tel: 1-408-855-0061 Fax: 1-408-855-0063
Email: info@LumiQuick.com www.LumiQuick.com

NEW You Tube Movies: www.youtube.com/LumiQuick

Product Specification

LumiQuick presents its products with well designed QuickView™ and QuickProfile™ package. The standard rapid test package contains 25 test cards or 50 test strips. The standard ELISA kit contains 96 tests.

Besides LumiQuick's own package, we also provide OEM and private label service for customers. The products can be in finished kit or bulk package format.

LumiQuick's products are for in-vitro diagnostic use only. All the products and reagents are warranted until the date of expiration.

Order Information

When you place the order, please include at least the following information.

1. Complete shipping address
2. Complete billing address
3. Purchase order number
4. Correct catalog number, product name and quantity
5. Any special labeling and packaging instructions
6. Shipping instructions

For prompt service, please send your order by fax to 1-408-855-0063, email to purchasing@lumiquick.com or your sales representative.

All orders are F.O.B. from our factories. Freight charges and insurance charges are the responsibility of the customer.

LumiQuick



DIAGNOSTICS, INC.


**2946 Scott Blvd.
Santa Clara, California
95054 USA**

www.LumiQuick.com



www.youtube.com/lumiquick

You Tube



**Questions? Please contact us at
info@lumiquick.com**

Or, your local agent

For Authorized Distributors Use