

Product catalog

This document reflects your selection of products from the Hamilton Medical online product catalog. The content is based on the following parameters:

Creation Date: **20 January 2022**
Country settings: **International**

Your selection:

- Accessories & consumables**
 - Flow sensors
- Ventilators**
 - HAMILTON-C1
- Patient groups**
 - Adult

Search keywords

The online product catalog is subject to change without notice. The latest and most comprehensive product information is available online.

Legend

 Single use

 Adult



The packaging unit of each single product is marked within this box icon.

 Reusable

 Pediatric

 Autoclavable

 Neonatal

HAMILTON
MEDICAL

Flow sensors

Product

Description



281637

Flow sensor



OD15 x ID15 / OD22 connectors



For ventilator: HAMILTON-C6/S1/G5/C3/C2/C1/T1/MR1
GALILEO, RAPHAEL



GTIN: 07630002800860

Dimensions: L: 188 cm / 74 in

Includes: [PN 279937](#) Calibration adapter for flow sensor



282049

Flow sensor



OD15 x ID15 / OD22 connectors



For ventilator: HAMILTON-C6/S1/G5/C3/C2/C1/T1/MR1
GALILEO, RAPHAEL



GTIN: 07630002800877

Dimensions: L: 260 cm / 102 in

Includes: [PN 279937](#) Calibration adapter for flow sensor



282051

Flow sensor



OD15 x ID15 / OD22 connectors



For ventilator: HAMILTON-C6/S1/G5/C3/C2/C1/T1/MR1
GALILEO, RAPHAEL



GTIN: 07630002800891

Dimensions: L: 330 cm / 130 in

Includes: [PN 279937](#) Calibration adapter for flow sensor



279937

Calibration adapter for flow sensor



For use with flow sensor PN281637, PN282049, PN282051



For ventilator: HAMILTON-C6/S1/G5/C3/C2/C1/T1/MR1
GALILEO



GTIN: 07630002800914



950185

Flow sensor



OD15 x ID15 x OD22 connectors

Includes 3 tubing attachment straps



For ventilator: HAMILTON-C6/S1/G5/C3/C2/C1/T1/MR1
GALILEO, RAPHAEL



GTIN: 07630002802529

Dimensions: L: 188 cm / 74 in

Notes: The accuracy of the sensor is reduced to +/- 15% or 15 ml whichever is greater, even after calibration when used with the RAPHAEL ventilator.

Includes: [PN 282323](#) Calibration adapter for flow sensor



Hamilton Medical AG

Via Crusch 8, 7402 Bonaduz, Switzerland

☎ +41 (0)58 610 10 20

info@hamilton-medical.com

www.hamilton-medical.com

Specifications are subject to change without notice. Some features are options. Not all features are available in all markets. INTELLiVENT-ASV is not available in the US. For all proprietary trademarks and third-party trademarks used by Hamilton Medical AG, see www.hamilton-medical.com/trademarks. © 2021 Hamilton Medical AG. All rights reserved.