



MAIN APPLICATIONS

Emergency Dept.

General Surgery

Gynecology

Liposuction

Neurosurgery

Obstetrics

Operating Theatre

Dental Practice

NEW HOSPIVAC 350 has been designed for professional aspiration of bodily fluids, tissues or bones of the patient during or after surgery. The state-of-the-art oilless pump guarantees high performances with excellent suction capacities and max vacuum reachable within a few seconds, with no need for maintenance. A clear dashboard along with a full range of accessories makes it the ideal device for surgical suction and liposuction.

| AVAILABLE MODELS AND ACCESSORIES INCLUDED | BASIC | FS | FULL |
|---|-------|----|------|
| Liquid Collection Jar with overflow valve system (different options, see below) | o | o | o |
| 8x14mm Silicone Tubes, length 150cm (autoclavable) + Conical Connector ø 10-11-12mm | o | o | o |
| Antibacterial & Hydrophobic Filter (single-patient) | o | o | o |
| Power Cord with Schuko plug | o | o | o |
| Footswitch (may be ordered at a later time) | | o | o |
| Electronic Change-Over System from jar to jar | | | o |

TECHNICAL FEATURES

| | |
|----------------------------|---|
| Motor | Oilless and maintenance-free piston pump |
| Power Feeding | 230V-50Hz - other configurations on request |
| ISO 10079-1 Classification | HIGH VACUUM / HIGH FLOW |
| Max Vacuum (adjustable) | -0.90 bar -90 kPa -675 mmHg |
| Max free air flow rate | 60 l/min |
| Noise Level | 51,7 dBA |
| Duty cycle | Non-stop operation |
| Weight | kg 13 |
| Size | cm 46 x 85 x 42 |
| Years of Warranty | 2 |
| Shipping carton | 1 |
| Place of Manufacturing | Italy |

Italian
Medical
Touch
MADE IN ITALY

NEW HOSPIVAC CONFIGURATIONS AND JARS INTERCHANGE

HOSPIVAC series offer a great versatility to users when it comes to change type of jars for liquid collection, regardless in which configuration the suction unit has been supplied in origin. A full line of accessories complete the range. The tables in the next pages provide all the information for ordering.

CE 0123

NEW HOSPIVAC 350 - AVAILABLE MODELS and ACCESSORIES INCLUDED

| | | AUTOCLAVABLE JARS | | | | DISPOSABLE LINERS | | | |
|--------------|----------------|----------------------------|-----------|-------------------|-----------|-------------------|-----------|-------------|-------------|
| | | Makrolon® 2858 (max 121°C) | | Apec® (max 143°C) | | Flovac® systems | | | |
| REF. | MODEL | 2x2000 ml | 2x5000 ml | 2x2000 ml | 2x5000 ml | 2x2000 ml | 2x3000 ml | FOOT SWITCH | CHANGE OVER |
| RE 410356 | BASIC 2 | | | | | | | | |
| RE 410356/06 | FS 2 | | | | | | | | |
| RE 410356/01 | FULL 2 | | | | | | | | |
| RE 410356/39 | BASIC 5 | | | | | | | | |
| RE 410356/40 | FS 5 | | | | | | | | |
| RE 410356/41 | FULL 5 | | | | | | | | |
| upon request | AP BASIC 2 | | | | | | | | |
| upon request | AP FS 2 | | | | | | | | |
| upon request | AP FULL 2 | | | | | | | | |
| upon request | AP BASIC 5 | | | | | | | | |
| upon request | AP FS 5 | | | | | | | | |
| upon request | AP FULL 5 | | | | | | | | |
| RE 410356/27 | FLOVAC BASIC 2 | | | | | | | | |
| RE 410356/29 | FLOVAC FS 2 | | | | | | | | |
| RE 410356/28 | FLOVAC FULL 2 | | | | | | | | |
| RE 410356/02 | FLOVAC BASIC 3 | | | | | | | | |
| RE 410356/09 | FLOVAC FS 3 | | | | | | | | |
| RE 410356/30 | FLOVAC FULL 3 | | | | | | | | |