

DATA & ESSENTIALS SHEET

SST VACUUM BLOOD COLLECTION TUBES WITH SERUM SEPERATION GEL & CLOT ACTIVATOR



1.	General Information	
1.1.	Description:	<p>TTT WORLD® SST Vacuum blood collection tubes. Manufactured by TTT SAĞLIK HİZMETLERİ EĞİT. İNŞ. SAN. VE TİC. A.Ş.</p> <p>Frequently used for clinical biochemistry and immunolgy serogy testing.</p> <p>The tubes are produced in various volumes.</p> <p>With the gel, a barrier is created between the serum formed upon centrifugation. The gel forms a barrier in between.</p> <p>The inner surface of the tube is coated with the sprayed and dried clot activator to accelerate coagulation.</p>
1.2.	Advantages for user:	<p>TTT WORLD® Blood collection tubes are produced in accordance to the international quality standards as with appropriate additive, accurate vacuum capacity and proper no-leakage caps.</p> <p>TTT WORLD® Blood collection Tubes are comfortable, ergonomic and useful for both patient and practitioner.</p>
1.3.	Safety:	<p>Shall be dispose according to the medical waste rules. For Single Use Only.</p> <p>Please check the usage information on the carton and IFU provided.</p>



1.4.	How to Use:	TTT WORLD® Blood collection tubes should be used in accordance with national and international phlebotomy blood collection guidelines.
2	Design and Components	
2.1.	Composed of:	PET Tube, Rubber Gasket, Cap.
2.2.	Additive:	Serum Separation Gel Clot Activator
2.3.	Volume:	From 0.5 ml pediatric tubes to 9 ml tubes, all volumes are available.
2.4.	Tube Sizes:	0,5 ml-1 ml: Micro Tubes 1 ml- 4 ml: 13x75 mm Tube 4 ml- 5.5 ml: 13x100 mm Tube 6-9 ml: 16x100 mm Tube The volumes can be adjusted.
2.6.	Shelf Life:	Indicated on the individual label.
3	Complied Standards and Certificates	
3.1.	CE Marking	Yes
3.2.	Directive	98/79 / EC IVD Medical Devices Directive
3.3.	Class	IVD
3.4.	Other Certificates	ISO 13485
3.5.	GMDN	41128

