

Biliary Stents – Teflon

Used for drainage of obstructed biliary ducts
Material TEFLON (PTFE)



Art. No.	Size (Fr/CH)	Length (cm)
BSST 7/5	7	5
BSST 7/7	7	7
BSST 7/10	7	10
BSST 7/12	7	12
BSST 7/15	7	15
BSST 8.5/5	8.5	5
BSST 8.5/7	8.5	7
BSST 8.5/10	8.5	10
BSST 8.5/12	8.5	12
BSST 8.5/15	8.5	15

Art. No.	Size (Fr/CH)	Length (cm)
BSST 10/5	10	5
BSST 10/7	10	7
BSST 10/10	10	10
BSST 10/12	10	12
BSST 10/15	10	15
BSST 11.5/5	11.5	5
BSST 11.5/7	11.5	7
BSST 11.5/10	11.5	10
BSST 11.5/12	11.5	12
BSST 11.5/15	11.5	15

Biliary Stents – Double Pigtail

Used for drainage of obstructed biliary ducts
Material PEBAX



Art. No.	Size (Fr/CH)	Length (cm)
BS-DP 5/4	5	4
BS-DP 5/5	5	5
BS-DP 5/7	5	7
BS-DP 5/10	5	10
BS-DP 5/12	5	12
BS-DP 5/15	5	15
BS-DP 7/4	7	4
BS-DP 7/5	7	5
BS-DP 7/7	7	7
BS-DP 7/9	7	9
BS-DP 7/10	7	10
BS-DP 7/12	7	12
BS-DP 7/15	7	15
BS-DP 8.5/4	8.5	4
BS-DP 8.5/5	8.5	5
BS-DP 8.5/7	8.5	7
BS-DP 8.5/10	8.5	10
BS-DP 8.5/12	8.5	12
BS-DP 8.5/15	8.5	15

Art. No.	Size (Fr/CH)	Length (cm)
BS-DP 10/4	10	4
BS-DP 10/5	10	5
BS-DP 10/7	10	7
BS-DP 10/10	10	10
BS-DP 10/12	10	12
BS-DP 10/15	10	15
BS-DP 11.5/4	11.5	4
BS-DP 11.5/5	11.5	5
BS-DP 11.5/7	11.5	7
BS-DP 11.5/10	11.5	10
BS-DP 11.5/12	11.5	12
BS-DP 11.5/15	11.5	15

Biliary Stents – Standard

Used for drainage of obstructed biliary ducts
Material PEBAX



Art. No.	Size (Fr/CH)	Length (cm)
BSS 7/5	7	5
BSS 7/7	7	7
BSS 7/9	7	9
BSS 7/10	7	10
BSS 7/12	7	12
BSS 7/15	7	15
BSS 8.5/5	8.5	5
BSS 8.5/7	8.5	7
BSS 8.5/9	8.5	9
BSS 8.5/10	8.5	10
BSS 8.5/12	8.5	12
BSS 8.5/15	8.5	15

Art. No.	Size (Fr/CH)	Length (cm)
BSS 10/5	10	5
BSS 10/7	10	7
BSS 10/8	10	8
BSS 10/9	10	9
BSS 10/10	10	10
BSS 10/12	10	12
BSS 10/15	10	15
BSS 11.5/5	11.5	5
BSS 11.5/7	11.5	7
BSS 11.5/9	11.5	9
BSS 11.5/10	11.5	10
BSS 11.5/12	11.5	12
BSS 11.5/15	11.5	15

Biliary Stents – Amsterdam

Used for drainage of obstructed biliary ducts
Material PEBAX



Art. No.	Size (Fr/CH)	Length (cm)
BSA 7/5	7	5
BSA 7/7	7	7
BSA 7/10	7	10
BSA 7/12	7	12
BSA 7/15	7	15
BSA 8.5/5	8.5	5
BSA 8.5/7	8.5	7
BSA 8.5/10	8.5	10
BSA 8.5/12	8.5	12
BSA 8.5/15	8.5	15

Art. No.	Size (Fr/CH)	Length (cm)
BSA 10/5	10	5
BSA 10/7	10	7
BSA 10/9	10	9
BSA 10/10	10	10
BSA 10/12	10	12
BSA 10/15	10	15
BSA 11.5/5	11.5	5
BSA 11.5/7	11.5	7
BSA 11.5/10	11.5	10
BSA 11.5/12	11.5	12
BSA 11.5/15	11.5	15

Pancreatic Stents – Curved

Used for drainage of obstructed biliary ducts
Material PEBAX + PTFE



Art. No.	Size (Fr/CH)	Length (cm)
PS-C 5/3	5	3
PS-C 5/5	5	5
PS-C 5/7	5	7
PS-C 5/9	5	9
PS-C 5/10	5	10
PS-C 5/12	5	12
PS-C 5/15	5	15

Art. No.	Size (Fr/CH)	Length (cm)
PS-C 8.5/3	8.5	3
PS-C 8.5/5	8.5	5
PS-C 8.5/7	8.5	7
PS-C 8.5/9	8.5	9
PS-C 8.5/10	8.5	10
PS-C 8.5/12	8.5	12
PS-C 8.5/15	8.5	15

PS-C 7/3	7	3
PS-C 7/5	7	5
PS-C 7/7	7	7
PS-C 7/9	7	9
PS-C 7/10	7	10
PS-C 7/12	7	12
PS-C 7/15	7	15

PS-C 11.5/3	11.5	3
PS-C 11.5/5	11.5	5
PS-C 11.5/7	11.5	7
PS-C 11.5/9	11.5	9
PS-C 11.5/10	11.5	10
PS-C 11.5/12	11.5	12
PS-C 11.5/15	11.5	15

Pancreatic Stents – Pigtail

Used for drainage of obstructed biliary ducts
Material PEBAX



Art. No.	Size (Fr/CH)	Length (cm)
PS-P 5/3	5	3
PS-P 5/5	5	5
PS-P 5/7	5	7
PS-P 5/9	5	9
PS-P 5/10	5	10
PS-P 5/12	5	12
PS-P 5/15	5	15

Art. No.	Size (Fr/CH)	Length (cm)
PS-P 8.5/3	8.5	3
PS-P 8.5/5	8.5	5
PS-P 8.5/7	8.5	7
PS-P 8.5/9	8.5	9
PS-P 8.5/10	8.5	10
PS-P 8.5/12	8.5	12
PS-P 8.5/15	8.5	15


PS-P 7/3	7	3
PS-P 7/5	7	5
PS-P 7/7	7	7
PS-P 7/9	7	9
PS-P 7/10	7	10
PS-P 7/12	7	12
PS-P 7/15	7	15

PS-P 11.5/3	11.5	3
PS-P 11.5/5	11.5	5
PS-P 11.5/7	11.5	7
PS-P 11.5/9	11.5	9
PS-P 11.5/10	11.5	10
PS-P 11.5/12	11.5	12
PS-P 11.5/15	11.5	15

Biliary Stent Introducer System

Used for endoscopically biliary stent placement

→ Radiopaque marker over guiding catheter



Art. No.	Productname	Size (Fr/CH)	Required Guidewire	Length (cm)
SSIS-8.5	Biliary Stent Introducer System	8.5	0.038"	220
SSIS-10	Biliary Stent Introducer System	10	0.038"	220
SSIS-11.5	Biliary Stent Introducer System	11.5	0.038"	220


Supplied without Guidewire and Stent

Pusher only

BSI 5	Pusher for 5Fr Biliary Stents	170
BSI 7	Pusher for 7Fr Biliary Stents	170

NEG Nitinol Guidewire

Used to gain access through difficult biliary and pancreatic strictures and to facilitate exchange of compatible wire guided accessories to the desired biliary or pancreatic ducts.




Flexible tip

Art. No.	Guidewire Ø	Tip	Length (cm)
NEG-025-S	0.025"	Straight	450
NEG-025-A	0.025"	Angled	450
NEG-032-S	0.032"	Straight	450
NEG-032-A	0.032"	Angled	450
NEG-035-S	0.035"	Straight	450
NEG-035-A	0.035"	Angled	450
NEG-038-S	0.038"	Straight	450
NEG-038-A	0.038"	Angled	450

Hydra NEG Nitinol Guidewire

Used to gain access through difficult biliary and pancreatic strictures and to facilitate exchange of compatible wire guided accessories to the desired biliary or pancreatic ducts



With Extra Filiform Wire

Super flexible tip

Art. No.	Guidewire Ø	Tip	Length (cm)
HYNEG-032-S	0.032"	Straight	450
HYNEG-032-A	0.032"	Angled	450
HYNEG-035-S	0.035"	Straight	450
HYNEG-035-A	0.035"	Angled	450
HYNEG-038-S	0.038"	Straight	450
HYNEG-038-A	0.038"	Angled	450