


3.1 Composition of Study List

3.1.1 Database



Item	Description
1	Database Ribbon Toolbar Displays a toolbar available in the Database tab of the Study List window.
2	Study List Displays a study list searched from the database tab.
3	Thumbnail List Displays thumbnails of the selected study.
4	User-defined screen layout Sets the view screen configuration for the selected study.
5	Report Displays report information of the selected study.
6	Past Study List Displays past study list for selected patients.
7	Status Bar Displays functions and their status used from the Database tab. (Refer to <3.1.3 Status Bar>)

3.1.2 User-Defined Screen Layout

Image	Description
	<ul style="list-style-type: none"> Configures the view screen with the specified layout. Resets the layout.