

Patient Cables

Connection cables for ablation catheters



Product Highlights

- **Flexibility**
All BIOTRONIK ablation catheters connect to all major RF ablation generator

Patient cables can be resterilized.
Amount of resterilization cycles and procedure according to technical manual.

Ordering Information

Model	Length	Order number
PK-109	2.5 m	330 970
PK-111	2.0 m	330 080
PK-112	2.0 m	330 081
PK-142	2.5 m	362 442
PK-143	2.5 m	371 597
PK-185	2.5 m	345 945

Patient Cables

Technical Data

Patient cable	Cable length	Suitable for connection to	Order number
PK-109	2.5 m	Stockert EP-Shuttle	330 079
PK-111	2.0 m	Osyпка HAT 300 Smart	330 080
PK-112	2.0 m	Medtronic Atakr, Radionics RFG 3-E	330 081
PK-142	2.5 m	IBI 1500	362 442
PK-143	2.5 m	Boston Scientific Maestro 3000	371 597
PK-185	2.5 m	BIOTRONIK AbControl MDS	345 954

Generator Connection Side

Ablation Catheter Connection Side



Patient Cables

Connection cables for diagnostic catheters



Product Highlights

- **Flexibility**

All BIOTRONIK diagnostic catheters connect to Stimulators and EP recording systems

MPK-R cables can be resterilized.
Amount of resterilization cycles and procedure according to technical manual.

Ordering Information

Model	Length	Order number
MPK-4-R	2.0 m	353 177
MPK-10-R	2.0 m	353 173

Patient Cables

Technical Data

Patient cable	Cable length	Suitable for connection to	Order number
MPK-4-R	2m	MultiCath: 2 and 4 poles Woxx: 4 poles ViaCath: 4 poles	353 177
MPK-10-R	2m	MultiCath: 5, 6 and 10 poles Woxx: 10 poles ViaCath: 8 and 10 poles Lexx: 10 poles	353 173

MultiCath: 2 and 4 poles, Woxx: 4 poles, ViaCath: 4 poles

MPK-4-R



MultiCath: 5, 6 and 10 poles, Woxx: 10 poles, ViaCath: 8 and 10 poles, Lexx: 10 poles

MPK-10-R

