

# AKASPRAY TÜCHER

Disinfection Wipes for Surfaces – Alcohol Based



## USAGE AREAS:

- Used for surface disinfection in every place where health services are provided such as operating rooms, intensive care, dental units, doctor's office, veterinary clinic, laboratory and public areas such as home, office, factory, school etc.

## PRODUCT PROPERTIES:

- Ready-to-use product.
- Does not contain phenol and aldehyde.
- Effective in 30 seconds
- Dries quickly and is pleasantly fragrant.
- High material compatibility
- Biocidal product license dated 08.06.2018 and numbered 2018/60, obtained from the Turkish Ministry of Health

## COMPOSITION:

- Ethanol, 2-propanol, didecyldimethylammonium propionate

## MICROBIOLOGICAL ACTIVITY:

- Bactericidal, fungicidal, virucidal, tuberculocidal in 30 seconds

## USAGE:

- Sufficient amount of wipes are taken according to the size of the equipment to be wiped. It is applied by wiping method as a surface disinfectant. It is waited for 30 seconds to take effect.

## PACKAGE:

| QUANTITY | PACKAGING TYPE       | BARCODE       |
|----------|----------------------|---------------|
| 100 pcs  | 12 x<br>CYLINDER BOX | 8680097103574 |



SURFACE  
DISINFECTANTS

**detrox**<sup>®</sup>  
Integrated Disinfection Solutions