


fixLOCK Anterolateral Distal Tibia Plate, 3.5 mm

Item Code: 936

Indications

- Comminuted fracture distal end tibia- periarticular & intrarticular
- Osteotomy around distal tibia
- Management of osteoporotic fractures
- Fracture non-union distal end tibia



Features

- Anatomically shaped- right or left
- Low profile construct with minimum implant prominence or soft tissue irritation
- Helical twist ensures well contoured plate for laterally placed shaft holes and anteriorly placed metaphyseal segment holes
- fixLOCK screws at the curvature to ensure strength to the construct even if the screw purchase is compromised
- Tapered tip for submuscular insertion with minimal soft tissue stripping
- Uniform stiffness of all segments, prolonging fatigue life of the implant
- Smaller area of post-fixation avascularity
- Available in both Titanium & Stainless Steel

Sizes Available

No. of Holes	Stainless Steel		Titanium	
	Code- Left	Code- Right	Code- Left	Code- Right
4	936.04L	936.04R	936.04TL	936.04TR
6	936.06L	936.06R	936.06TL	936.06TR
8	936.08L	936.08R	936.08TL	936.08TR
10	936.10L	936.10R	936.10TL	936.10TR
12	936.12L	936.12R	936.12TL	936.12TR
14	936.14L	936.14R	936.14TL	936.14TR