


# EMR ClearGrasp Snare



## FEATURE & BENEFITS

### (1) Built-in saline Irrigation port

- This product has the saline irrigation port. Using it, the users do not need to exchange accessories to flush the saline or water to detect the bleeding site.

### (2) Looping rotation with the handle at 1:1 ratio

- The loop of this product can be rotated at 1:1 ratio with the handle, and the users can perform the EMR procedures more easily.

### (3) Various shapes and sizes of loop (6mm-30mm)

- The users can select the shape and diameter of loop according to their needs.
- Cold Snare Polypectomy can be performed by using our small size snare (6-10mm)


### (4) Cold & Hot Snare




- Shield shape enables effective cold snare Polypectomy.
- Covering both Hot and Cold Snare Polypectomy.
- 6mm loop size facilitates the resection of diminutive polyps (size ≤5mm)
- Irrigation Port helps control bleeding during polypectomy.




### OVAL 1 TYPE




### CRAB TYPE



### CRESCENT TYPE




### COLD TYPE



**SPECIFICATION**

Type	Model	Type of Retation	Loop Size	Loop Stiffness	Sheath O.D	Working Length	Number of Use	Validity of Sterilization
Oval1	FM-ESN004(2400)(Soft)	Non-Rotation	6mm	Soft	ø2.4mm	2,400mm	Single Use Only	3 years
	FM-ESN006(2400)(Stiff)			Stiff				
	FM-ESN006(2400)(Moderato)			Moderato				
	FM-ES0006(2400)(Soft)	Rotation	6mm	Soft				
	FM-ES0006(2400)(Stiff)			Stiff				
	FM-ES0006(2400)(Moderato)			Moderato				
	FM-ESN004(2400)(Soft)	Non-Rotation	10mm	Soft				
	FM-ESN004(2400)(Stiff)			Stiff				
	FM-ESN004(2400)(Moderato)			Moderato				
	FM-ES0004(2400)(Soft)	Rotation	10mm	Soft				
	FM-ES0004(2400)(Stiff)			Stiff				
	FM-ES0004(2400)(Moderato)			Moderato				
	FM-ESN005(2400)(Soft)	Non-Rotation	15mm	Soft				
	FM-ESN005(2400)(Stiff)			Stiff				
	FM-ESN005(2400)(Moderato)			Moderato				
	FM-ES0005(2400)(Soft)	Rotation	15mm	Soft				
	FM-ES0005(2400)(Stiff)			Stiff				
	FM-ES0005(2400)(Moderato)			Moderato				
FM-ESN002(2400)(Soft)	Non-Rotation	24mm	Soft					
FM-ESN002(2400)(Stiff)			Stiff					
FM-ESN002(2400)(Moderato)			Moderato					
FM-ES0002(2400)(Soft)	Rotation	24mm	Soft					
FM-ES0002(2400)(Stiff)			Stiff					
FM-ES0002(2400)(Moderato)			Moderato					
FM-ESN003(2400)(Soft)	Non-Rotation	30mm	Soft					
FM-ESN003(2400)(Stiff)			Stiff					
FM-ESN003(2400)(Moderato)			Moderato					
FM-ES0003(2400)(Soft)	Rotation	30mm	Soft					
FM-ES0003(2400)(Stiff)			Stiff					
FM-ES0003(2400)(Moderato)			Moderato					
Crab	FM-ESN006(Crab)(2400)(Stiff)	Non-Rotation	6mm	Stiff				
	FM-ES0006(Crab)(2400)(Stiff)	Rotation	6mm					
	FM-ESN004(Crab)(2400)(Stiff)	Non-Rotation	10mm					
	FM-ES0004(Crab)(2400)(Stiff)	Rotation	10mm					
	FM-ESN005(Crab)(2400)(Stiff)	Non-Rotation	15mm					
	FM-ES0005(Crab)(2400)(Stiff)	Rotation	15mm					
	FM-ESN002(Crab)(2400)(Stiff)	Non-Rotation	24mm					
	FM-ES0002(Crab)(2400)(Stiff)	Rotation	24mm					
FM-ESN003(Crab)(2400)(Stiff)	Non-Rotation	30mm						
FM-ES0003(Crab)(2400)(Stiff)	Rotation	30mm						
Crescent	FM-ES0004(Crescent)(2400)	Non-Rotation	10mm	Soft				
	FM-ES0004(Crescent)(2400)		15mm	Stiff				
	FM-ES0005(Crescent)(2400)		15mm	Soft				
	FM-ES0005(Crescent)(2400)		24mm	Stiff				
	FM-ES0002(Crescent)(2400)		24mm	Soft				
Cold	FM-ESN006(Cold)(2400)	Non-Rotation	6mm	Soft (0.2mm)				
	FM-ES0006(Cold)(2400)	Rotation	6mm					
	FM-ESN004(Cold)(2400)	Non-Rotation	10mm					
	FM-ES0004(Cold)(2400)	Rotation	10mm					

# Clear-CoaBite



**FEATURE & BENEFITS**

**(1) Built-in Saline Irrigation Port**

• This product has the saline Irrigation port, This could let the users don't need to change snare to flush saline or water for detecting bleeding site.

**(2) Improving usability**

• By shortening the conductive length at the distal end of sheath from 9mm to 4mm, operators do not need to focus out in order to take the angle again.



**SPECIFICATION**

Model	Shape of Cup	Cup(Jaw)O.D	Sheath O.D	Working Length	Compatible Endoscopes		Number of Use	Validity of Sterilization
					Length	Channal Inner Diameter		
FM-EB0001	Oval	ø1.72mm	ø2.3mm	1,600mm	Max. 1,100mm	Min. ø2.8mm	Single Use Only	3 years
FM-EB0002				2,400mm	Max. 1,800mm			
FM-EB0003				1,600mm	Max. 1,100mm			
FM-EB0004				2,400mm	Max. 1,800mm			