


Flexor® Ureteral Access Sheath

Used to establish a conduit during endoscopic urological procedures, facilitating the passage of endoscopes and other instruments into the urinary tract. The working channel protects the ureter during repeated instrument exchanges. The sheath also protects delicate instruments and smaller flexible scopes from damage. The sheath has a hydrophilic coating to create a low-friction surface as the ureter is dilated.



Order Number	Reference Part Number	Sheath ID Fr	Length cm
G32853	FUS-095013*	9.5	13
G19177	FUS-095020	9.5	20
G19176	FUS-095028	9.5	28
G19174	FUS-095035	9.5	35
G46167	FUS-095045	9.5	45
G19175	FUS-095055	9.5	55
G34876	FUS-107013*	10.7	13
G34877	FUS-107020	10.7	20
G34878	FUS-107028	10.7	28
G34879	FUS-107035	10.7	35
G34880	FUS-107045	10.7	45
G34881	FUS-107055	10.7	55

Order Number	Reference Part Number	Sheath ID Fr	Length cm
G32852	FUS-120013*	12.0	13
G19173	FUS-120020	12.0	20
G19172	FUS-120028	12.0	28
G19168	FUS-120035	12.0	35
G46168	FUS-120045	12.0	45
G19169	FUS-120055	12.0	55
G30391	FUS-140020	14.0	20
G30393	FUS-140028	14.0	28
G30389	FUS-140035	14.0	35
G46169	FUS-140045	14.0	45
G30390	FUS-140055	14.0	55

* For paediatric use in children who are two years old or older.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Support & Delivery for details.

Please see product risk information in the IFU at cookmedical.eu.