



K3EDTA


Intended to be used in hematology laboratory for withdrawing Whole blood or EDTA Plasma in order to examine various kinds of blood cells.

Tubes manufactured using transparent PET or Glass. Caps manufactured using polyethylene.

Specimen: whole blood.

After taking blood into the tube it should be mixed by turning 6-8 times. Lack of mixing can also cause platelet aggregation, coagulation or formation of micro clots.

■ In addition, we can provide EDTA blood collection tubes with Na₂EDTA as additive. If you need to order, please contact our sales manager.



| REF | Additive | Cap Colour | Tube Size (mm) | Draw Volume (mL) | Tube Material | Pack/Carton | Label |
|--------------|----------|------------|----------------|------------------|---------------|----------------|-------------------|
| RJ1803101-0P | K3EDTA | Purple | 13x75 | 1mL | PET | 100PCS/1200PCS | Paper/Transparent |
| RJ1803101-0G | K3EDTA | Purple | 13x75 | 1mL | GLASS | 100PCS/1200PCS | Paper/Transparent |
| RJ1803102-0P | K3EDTA | Purple | 13x75 | 2mL | PET | 100PCS/1200PCS | Paper/Transparent |
| RJ1803102-0G | K3EDTA | Purple | 13x75 | 2mL | GLASS | 100PCS/1200PCS | Paper/Transparent |
| RJ1803103-0P | K3EDTA | Purple | 13x75 | 3mL | PET | 100PCS/1200PCS | Paper/Transparent |
| RJ1803103-0G | K3EDTA | Purple | 13x75 | 3mL | GLASS | 100PCS/1200PCS | Paper/Transparent |
| RJ1803104-0P | K3EDTA | Purple | 13x75 | 4mL | PET | 100PCS/1200PCS | Paper/Transparent |
| RJ1803104-0G | K3EDTA | Purple | 13x75 | 4mL | GLASS | 100PCS/1200PCS | Paper/Transparent |
| RJ1803205-0P | K3EDTA | Purple | 13x100 | 5mL | PET | 100PCS/1200PCS | Paper/Transparent |
| RJ1803205-0G | K3EDTA | Purple | 13x100 | 5mL | GLASS | 100PCS/1200PCS | Paper/Transparent |
| RJ1803206-0P | K3EDTA | Purple | 13x100 | 6mL | PET | 100PCS/1200PCS | Paper/Transparent |
| RJ1803206-0G | K3EDTA | Purple | 13x100 | 6mL | GLASS | 100PCS/1200PCS | Paper/Transparent |
| RJ1803207-0P | K3EDTA | Purple | 13x100 | 7mL | PET | 100PCS/1200PCS | Paper/Transparent |
| RJ1803207-0G | K3EDTA | Purple | 13x100 | 7mL | GLASS | 100PCS/1200PCS | Paper/Transparent |
| RJ1803308-0P | K3EDTA | Purple | 16x100 | 8mL | PET | 100PCS/1200PCS | Paper/Transparent |
| RJ1803308-0G | K3EDTA | Purple | 16x100 | 8mL | GLASS | 100PCS/1200PCS | Paper/Transparent |
| RJ1803309-0P | K3EDTA | Purple | 16x100 | 9mL | PET | 100PCS/1200PCS | Paper/Transparent |
| RJ1803309-0G | K3EDTA | Purple | 16x100 | 9mL | GLASS | 100PCS/1200PCS | Paper/Transparent |
| RJ1803310-0P | K3EDTA | Purple | 16x100 | 10mL | PET | 100PCS/1200PCS | Paper/Transparent |
| RJ1803310-0G | K3EDTA | Purple | 16x100 | 10mL | GLASS | 100PCS/1200PCS | Paper/Transparent |