



		Lunghezza del cavo (cm) / Cable length (cm):		F7981	F7981P
GE	RESPONDER 3000 (attacco diretto- <i>direct connection</i>)	●	●	200	200
	RESPONDER 2000 (attacco diretto- <i>direct connection</i>)	●	●		
	RESPONDER 1100 (attacco diretto- <i>direct connection</i>)	●	●		
	CARDIOSERV (attacco diretto- <i>direct connection</i>)	●	●		
	DASH RESPONDER (attacco diretto- <i>direct connection</i>)	●	●		




		Lunghezza del cavo (cm) / Cable length (cm):						
		F7986	F7986W	F7986P	F7986PW	F7986UW	F7986UW-D700	F7786UW-D700
Mediana	D500	●	●	●	●	-	-	-
	HEARTON AED A15	-	-	-	-	●	-	-
	D700	-	-	-	-	-	●	●
Tecno-Gaz	TECNOHEART PLUS	-	-	-	-	●	-	-



F7986 - F7986W - F7986P - F7986PW - F7986UW-D700 - F7786UW-D700 : Non commercializzabili in Italia - *Not for sale in Italy*

		Lunghezza del cavo (cm) / Cable length (cm):		F7987W	F7987PW
Zoll Medical Corp	AED PLUS	●	●	120	120
	AED PRO	●	●		
	ZOLL AED 3	●	●		
	R Series (#)	●	●		
	X Series (*)	●	●		



F7987W-F7987PW: senza funzione CPR (risposta alle compressioni e metronomo) - *without CPR function (feedback on compression and pace metronome)*

(#) Con cavo terapia CPR-D: 8000-000863-40 o adattatore: 8009-0020 (per cavi terapia: 8009-0749 o 8009-0750) - *With therapy cable CPR-D: 8000-000863-40 or adapter: 8009-0020 (for therapy cables: 8009-0749 or 8009-0750)*

(*) Con cavo terapia CPR-D: 8300-000831-40 o adattatore: 8009-0370 (per cavo terapia: 8300-000676) - *with therapy cable CPR-D: 8300-000831-40 or adapter: 8009-0370 (for therapy cable: 8300-000676)*