

Cleaning brush and valve cleaning brush

Reference	Designation	Working channel diameter (mm)	Catheter	Brush (mm)		
			Length and Colour (cm)	Single / Double / Triple	Length	Diameter
Endoscope small diameter : bronchial and pediatric						
CJ-EDB-100-02	Ureteroscope	1.2	100 White	●—●	6 / 6	2
CJ-EDB-120-02	Bronchial pediatric endoscope	1.5	120 White	●—●	12 / 12	2
CJ-EDB-180-02	Echo-endoscope	1.25	180 White	●—●	10 / 10	2
CJ-EDB-230-02	For Air/Water channel	1.2 - 1.5	230 White	●—●	12 / 12	2
CJ-EDB-120-03	Bronchial endoscope	≥ 2.0	120 Yellow	●—●	20 / 20	3
CJ-EDB-230-03	Pediatric endoscope	≥ 2.0	230 Yellow	●—●	20 / 20	3
CJ-EDB-300-03*	Pediatric Enteroscope	≥ 2.2	300 Yellow	●—●	20 / 20	3
Adult : Endoscope, enteroscope and cystoscope						
CJ-ESB-020-06	For beginning of Air/Water channel	≥ 2.8 (2.8 to 5)	20 Light blue	—●	20	6
CJ-EDB-120-06	Endoscope and cystoscope	≥ 2.8 (2.8 to 5)	120 Blue	●—●	20 / 20	6
CJ-ESB-230-06	Endoscope adult	≥ 2.8 (2.8 to 5)	230 Blue	—●	20	6
CJ-EDB-230-06	Endoscope adult	≥ 2.8 (2.8 to 5)	230 Blue	●—●	20 / 20	6
CJ-EDB-300-06*	Enteroscope	≥ 2.8 (2.8 to 5)	300 Blue	●—●	20 / 20	6
Valve cleaning brush						
CJ-GDB-11-05	Valve cleaning brush	/	20 Sky blue	●=●	32 / 20	11 / 5
Triple cleaning brush including a sweeper brush						
NLS/ETB-230 06S05	Endoscope for adult	≥ 2.8 (2.8 to 5)	230 Blue	●—●—●	20 / 20 /125	6

* Special request