

# Neo Protoezim Plus 500

## Multienzymatic concentrated detergent and disinfectant

### COMPOSITION

Total 8.0 g of enzymes, cationic and anionic surfactants 51.0 g, 0.3 g DPTA, coformulants 13.0 g, purified water q.s. to 100.0 ml.

### EFFECTIVENESS

Detergent with bactericidal, virucidal, fungicidal action. Contact times:

10 min., at a dilution of 1: 500 (10 ml in 5 litres of water).

In ultrasonic tanks reduced contact time - 5 min.

### USE

Perfect solution for:

- Decontamination of soiled instruments immediately after use;
- Cleaning and disinfection of laparoscopy, arthroscopy, gastroscopes and other medical devices, critical and semi-critical;

### AVAILABLE IN

Plastic bottle 1000 ml.



**Cantel Medical (Italy) S.r.l. a socio unico**

Via Laurentina, 169

00071 Pomezia (RM)-Italia