# **FORCEPS**





SOFT IOL FOLDING FORCEPS ANGLED, TITANIUM



Lot 196





■ 0.08 mm deep Pierse type tips for grasping lens



Round handle

#### **OF 103**

HALL LENS INSERTION FORCEPS ANGLED, TITANIUM OF 103L (WITH LOCK)



Round handle



Overall length 100 mm, 3.9 inches

ICL CARTRIDGE LOADING FORCEPS, TITANIUM



Overall length 130 mm, 5.1 inches

- Tips  $0.45 \times 0.3 \text{ mm}$
- Angled 30 degree, shafts 17 mm long
- Designed to remove the ICL from the container and push it inside the cartridge

≥ 0F092.0

Round handle

#### OF 028

IOL IMPLANTATION FORCEPS ANGLED, TITANIUM



■ Shafts 6.6 mm long



# Overall length 105 mm, 4.1 inches

■ Flat handle

### OF 092

FAULKNER IOL IMPLANTATION FORCEPS ANGLED, TITANIUM



Tips 0.3 mm dia



Overall length 102 mm, 4 inches

- Flat handle
- Round shafts 11 mm long

# OF 092.01

GODFREY IOL IMPLANTATION FORCEPS ANGLED, TITANIUM



Tips  $0.8 \times 0.3$  mm wide Shafts 12.3 mm long



#### Overall length 104 mm, 4.1 inches

- Flat handle



CARTRIDGE LOADING FORCEPS STRAIGHT, TITANIUM



# Overall length 110 mm, 4.3 inches

- Flat handle Tips 1.2 mm wide
- Designed for inserting Acrylic IOL into A, B, C, D To be used with injectors OJ 004, OJ 005, OJ 008, OJ 009 cartridges





# **FORCEPS**







Round handle

Overall length 155 mm, 6.0 inches

#### OF 121

IOL FOLDING FORCEPS ANGLED TITANIUM



Lot 196

Overall length 110 mm, 4.3 inches ■ Flat handle



OF 153

SIEBEL IOL FOLDING FORCEPS ANGLED, TITANIUM Lot 196





■ Flat handle

Overall length 110 mm, 4.3 inches

OF 154

BURATTO IOL IMPLANTATION FORCEPS ANGLED, TITANIUM



Round handle



Overall length 105 mm, 4.1 inches

OF 149

IOÙ FOLDING FORCEPS ANGLED, TITANIUM



Round handle



Overall length 110 mm, 4.3 inches

**OF 155** 

IOL FOLDING FORCEPS ANGLED, TITANIUM



Lot 196



Overall length 110 mm, 4.3 inches

Round handle



IOL IMPLANTATION FORCEPS ANGLED, TITANIUM



■ Flat handle



Overall length 110 mm, 4.3 inches