



**CENTRE OF TESTING SERVICE
INTERNATIONAL**

OPERATE ACCORDING TO ISO/IEC 17025

TEST REPORT

TEST REPORT NUMBER : CGZ3170720-02355-L

CTS (Ningbo) Testing Service Technology Co., Ltd.
Fl.1 & 8 West, Bldg. B, No. 66, Qingyi Rd., Hi-Tech Zone, Ningbo, Zhejiang, China

TEST REPORT

Table of contents

Applicant	Dongguan Hengsheng Polybag Co.,Ltd
Address	Li Jia Fang Industrial Park, Shipai Town, Dongguan, Guangdong, 523336, China
	Contact: /
	TEL: +86- 0769-8659 3128; FAX: +86- 0769-8659 3138
Manufacturer/ factory	Dongguan Hengsheng Polybag Co.,Ltd
Buyer	---
Sample Description	Autoclave bags
Trade Mark	---
Model/Type reference	AB-001
Taking/Receiving Date	20 July 2017
Test Address	CTS (Ningbo) Testing Service Technology Co., Ltd. Fl.1 & 8 West, Bldg. B, No. 66, Qingyi Rd., Hi-Tech Zone, Ningbo, Zhejiang, China Tel: +86-574-87912121 (16 lines) Fax: +86-574-87907993
Test Date	21 July 2017
Test Item*	High Temperature test
Test Requested*	High temperature test : 135°C, 30 min.
Test Results	See follow result
Note	All information was provided by the applicant or detected at the sample. This report shall not be copied partly without the written approval of CTS.

Tested by:

Reviewed by:


Approved by:



Date:
25 July 2017



Date:
25 July 2017



Date:
25 July 2017

* - means the content is provided by applicant.



Copyright of this report is owned by Centre of Testing Service and may not be reproduced other than in full except with the written approval of the issuing Company.

CTS (Ningbo) Testing Service Technology Co., Ltd.

Fl.1 & 8 West, Bldg. B, No. 66, Qingyi Rd., Hi-Tech Zone, Ningbo, Zhejiang, China

Tel: +86-574-87912121 (16 lines)

Fax: +86-574-87907993

Complaint line: +86-574-87908003

E-mail: cts@cts-lab.com

See Reverse For Terms And Conditions of Service

TEST REPORT**1. Test results**

No.	Test Sample Name	Test Item	Test Results	Verdict
1.	Autoclave bags AB-001	High Temperature test: 135 °C , 30 min.	After test, the black font and pattern on autoclave bags surface has some fade, and the high temperature discoloration became dark brown. The inside of the autoclave bags is slightly glued together; no visible deformation, blistering, melting.	---
Remark: Meet the requirement of the client.				

Copyright of this report is owned by Centre of Testing Service and may not be reproduced other than in full except with the written approval of the issuing Company.

CTS (Ningbo) Testing Service Technology Co., Ltd.

Fl.1 & 8 West, Bldg. B, No. 66, Qingyi Rd., Hi-Tech Zone, Ningbo, Zhejiang, China

Tel: +86-574-87912121 (16 lines)

Fax: +86-574-87907993

Complaint line: +86-574-87908003

E-mail: cts@cts-lab.com

See Reverse For Terms And Conditions of Service

TEST REPORT

2. Photos



Figure 1 (Before test - 1)



Figure 2 (Before test - 2)

Copyright of this report is owned by Centre of Testing Service and may not be reproduced other than in full except with the written approval of the issuing Company.

CTS (Ningbo) Testing Service Technology Co., Ltd.

Fl.1 & 8 West, Bldg. B, No. 66, Qingyi Rd., Hi-Tech Zone, Ningbo, Zhejiang, China

Tel: +86-574-87912121 (16 lines)

Fax: +86-574-87907993

Complaint line: +86-574-87908003

E-mail: cts@cts-lab.com

See Reverse For Terms And Conditions of Service

TEST REPORT



Figure 3 (After test - 1)



Figure 4 (After test - 2)

-----end of report-----

Copyright of this report is owned by Centre of Testing Service and may not be reproduced other than in full except with the written approval of the issuing Company.

CTS (Ningbo) Testing Service Technology Co., Ltd.

Fl.1 & 8 West, Bldg. B, No. 66, Qingyi Rd., Hi-Tech Zone, Ningbo, Zhejiang, China

Tel: +86-574-87912121 (16 lines)

Fax: +86-574-87907993

Complaint line: +86-574-87908003

E-mail: cts@cts-lab.com

See Reverse For Terms And Conditions of Service