

# **PRODUCT** CATALOG 2019

### Cleanis SASU

22 avenue Aristide Briand 94110 ARCUEIL Tel.: 01 46 15 18 00 – Fax: 01 82 66 13 92

www.cleanis.com - cleanis@cleanis.fr

### **CLEANIS COMPANY**



CLEANIS, subsidiary of INNOTHERA Group, offers products and services designed to prevent the spread of pathogenic agents and improve hygiene practices.



**CEANIS** foundation

e are both ambitious and modest, working over the long term with patience, determination, and good sense to continuously improve the treatment of common illnesses around the world.


We rely first and foremost on our energy, our individual and collective commitment, founded on the assimilation of many cultures and disciplines, in a spirit of trust and mutual aid.

**M. Arnaud Gobet**CEO of INNOTHERA Group



Japan subsidiary foundation

### **CLEANIS IS A FRENCH COMPANY**





### **ARCUEIL**

The CLEANIS - INNOTHERA group headquarters operates in Arcueil, Ile-de-France. It now has 170 employees with a great diversity of professions.



#### **NOMEXY**

CLEANIS factory, in the Vosges region, uses several advanced technologies to deliver premium products. It now gathers 28 employees.

# A CUTTING-EDGE FACTORY FOR HIGH-QUALITY PRODUCTS

Cleanis factory, in Vosges region, was designed to deliver products that meet your expectations and comply to ISO 13485 requirements. Before their commercialisation wordlwide, Cleanis products are carefully inspected, in order to confirm their compliance.











## A CONSTANT IMPROVEMENT OF QUALITY PROCESS

#### **CLEANIS** is certified ISO 13485

Guaranteeing that the company follows strict rules to ensure a trustful process to the benefit and protection of the patient





Our suppliers are selected, audited and evaluated



A controlled environment for production



A strict
Quality Control



A responsive and attentive Medical Device Vigilance Service



 ${\sf Email:information.vigilance@innothera.com}$ 

### A DEDICATED COMMERCIAL SERVICE

- Expert in body wastes management in healthcare facilities
- Constant support in:
   Analysis of your needs,
   Improvement of the body waste management in healthcare facilities
   Solution proposals





Our International Sales Manager is available for all request or order :

Camille PEZRON +33 6 31 08 52 64 cpezron@cleanis.com

### Devices for body waste collection with super absorbent pad



CareBag® products contribute to the prevention of Hospital-Acquired Infections (HAIs). They are devices for body waste collection and containement, provided with a super absorbent pad. This allows to protect the patient environment from pathogens causing HAIs, and to confine body waste that includes cytotoxic derivates of chemotherapy. CareBag® products are then disposed of according to facility protocols (biohazardous waste or waste similar to household one).



### **Bedpan Liner\***

Designed to adapt to bedpan and commode pail.

Product reference		771 11 41
Classification		Class I medical device, compliant with directive 93/42/CEE
Primary packaging		x 20 (roll in dispenser pack)
Secondary packaging		box of 18 packs
Bag	Dimensions (h x l)	38 x 59 cm
	Dimensions (h x l)	24 x 13 cm
Pad	Biologic fluids absorption (0.9% saline solution)	450 ml min. Absorption tested with a methode adapted from ISO 17190-5
Expiry date		None



### Male Urinal with holder\*

Product reference		961 01 18
Classification		Class I medical device, compliant with directive 93/42/CEE
Primary packaging		x 12 (unitary presentation in dispenser pack)
Secondary packaging		box of 21 packs
Bag	Dimensions (h x l)	38 x 16 cm
Holder	Dimensions (h x l)	13 x 6,5 cm
Holder	Material	Carton
	Dimensions (h x l)	24 x 13 cm
Pad	Biologic fluids absorption (0.9% saline solution)	450 ml min. Absorption tested with a methode adapted from ISO 17190-5
Expiry date		None



### Devices for body waste collection with super absorbent pad



### Vomit Bag with holder\*

The holder allows to rigidify the Vomit Bag, thus improving its handiness for the patient.

Product reference		961 01 44
Classification		Class I medical device, compliant with directive 93/42/CEE
Primary packaging		x 12 (unitary presentation in dispenser pack)
Secondary packaging		box of 21 packs
Bag	Dimensions (h x l)	25 x 25 cm
Halala.	Dimensions (h x l)	20 x 4,3 cm
Holder	Material	Carton
	Dimensions (h x l)	24 x 13 cm
Pad	Biologic fluids absorption (0.9% saline solution)	450 ml min. Absorption tested with a methode adapted from ISO 17190-5
Expiry date		None



### **Toilet Bowl Liner\***

It is adapted to all type of toilets, with or without toilet lid.

Product reference		953 85 34
Classification		Class I medical device, compliant with directive 93/42/CEE
Primary packaging		x 20 (roll in dispenser pack)
Secondary packaging		box of 18 packs
Bag	Dimensions (h x l)	38 x 65 cm
	Dimensions (h x l)	24 x 13 cm
Pad	Biologic fluids absorption (0.9% saline solution)	450 ml min. Absorption tested with a methode adapted from ISO 17190-5
Expiry date		None

Training kit to raise awareness regarding infectious risks linked to the emptying of bodily waste in toilets





### Kit CareBag SplashUV™

1<sup>st</sup> training kit to raise awareness regarding infectious risks, the Kit CareBag SplashUV™ is an innovative method to track the risks associated to traditional practices – flushing, manual cleaning.

The Kit allows to simulate body waste: micro droplets, resulting from pouring and flushing the mix of tap water and fluorescent solution into the toilets, become visible under the UV beam of the flashlight.

Product reference	982 92 00
Primary packaging	Unitary box with instructions for use integrated
Dimensions	75,6 x 36 x 15,3 cm
Kit content	5 syringes filled with a fluorescent liquid
	3 pairs of glasses
	1 UV flashlight
Warning	Read the instructions for use carefully. The syringes must be kept away from sunlight and humidity, at room temperature.



### Refill for Kit CareBag SplashUV™

Syringes of the Kit SplashUV<sup>™</sup> can be bought separately.

Product reference	982 92 01
Primary packaging	Pack of 5 syringes
Warning	Read the instructions for use carefully. The syringes must be kept away from sunlight and humidity, at room temperature.



### Devices for body waste collection with super absorbent pad



### **Pre-Calibrated Scale**

- Allows to measure urinary output
- Provided with a hanging hook
- Easily cleaned and disinfected through your facility guidelines

Product reference	961 02 18
Primary packaging	Unitary box
Secondary packaging	Box of 10 scales
Material	Acrylic tube (PMMA)
	Hanging hook: stainless steel
	Weighing hook: stainless steel (not in the scope of legal metrology)
Measure unity and accuracy	ml (accuracy: +/- 10ml)
Classification	Compliant with directive 2014/32/EU
Additional information	Reusable product. Compatible with your disinfectant solution



### Dispenser

Metal bag dispenser, compatible with CareBag® products in roll.

Product reference	602 00 01
Primary packaging	Box of 10 dispensers
	Dimensions : 21,3 x 12,5 x 13,2 cm
Characteristics	Pre-made holes for an easy installation on the wall (screws not included), or hanged onto the cart of the caregiver if hooks or other hanging systems are already present.

### The super absorbents for the safe management of hazardous fluids



GelMax<sup>®</sup> Super Absorbents are effective to limits risky manipulations which can contaminate the environment. These products are cost-effective, since they reduce the costs of hazardous liquid waste disposal. By turning fluids into a gel, GelMax<sup>®</sup> products are time-saving, and provide a safer work environment by making cleaning easier.





### Pad\*

GelMax<sup>®</sup> Pad is used to improve the absorption capacity of CareBag<sup>®</sup> products. It is also used in different containers like basins, in order to make liquid disposal easier.

Product reference	958 21 04
Classification	Class I medical device, compliant with directive 93/42/CEE
Primary packaging	x 25
Secondary packaging	Box of 21 packs
Dimensions (h x l)	24 x 13 cm
Biologic fluids absorption (0.9% saline solution)	450 ml min. Absorption tested with a methode adapted from ISO 17190-5
Expiry date	None



### The super absorbents for the safe management of hazardous fluids







### Water Soluble Pouch\*

GelMax® Water Soluble Pouch is used to lower the disposal cost of hazardous fluids. Its hydrosoluble cover dissolves upon contact with liquids. Thanks to its little size, it can be used with little containers - used ostomy pouches, suction canister - and in laboratories for medical biology.

Product reference	957 84 10
Classification	Class I medical device, compliant with directive 93/42/CEE
Primary packaging	x 100
Secondary packaging	Box of 18 packs
Dimensions (h x l)	12 x 3 cm
Material	Sodium polyacrylate powder (7g) Water soluble cover in poly-vynil alcohol (PVA)
Biologic fluids absorption (0.9% saline solution)	450 ml min. Absorption tested with a methode adapted from ISO 17190-5
Causes impacting absorption capacity	Liquid temperature Ionic strength of the liquid
Expiry date	3 years

#### **Product**

### Average absorption volume of 1 GelMax® Pouch

0.9% saline solution	450 ml
Hydrochloric acid 3% <sup>1</sup>	70 ml
Peracetic acid 10% <sup>1</sup>	210 ml
Picric acid 1.3% <sup>1</sup>	210 ml
Bleach 2%	280 ml
Sodium hydroxide 10% <sup>1</sup>	170 ml
Ethanol 20%¹	1180 ml
Formalin 3%	1700 ml
Aniosyme <sup>1, 2</sup>	1 500 ml

<sup>1.</sup> According to results obtained by LNE (Laboratoire national de métrologie et d'essais), following a method adapted from ISO 17190:5 standard. Test report and accurate results available upon request.

<sup>2.</sup> Disinfectant distributed by Anios at 0.5%



### Pre-moistened gloves for bathing without rinsing

Aqua® gloves are a real help for fragile or bedridden patients, and caregivers. Developed for fragile skins, Aqua® gloves are pre-moistened with a no-rinsing, dermatologically tested lotion, with Aloe vera. Alcool free, Aqua® gloves are microwaveable, providing more comfort.

Aqua® gloves are frequently used in palliative care, geriatrics, rehabilitation but also after diaper use.



### **Total Hygiene Glove**

Total Hygiene glove is intended to body and face cleaning of fragile or bedridden patients, without rinsing. It suits all kind of patient's skin, including children older than 3 years old.

Product reference	462 82 09
Primary packaging	x 12 (pack)
Secondary packaging	Carton of 12 packs
Glove	17 x 23,5 cm
Lotion	Dermatologically tested lotion, Aloe Vera-based
Expiry date	3 years from manufacture date



### Pre-moistened gloves for bathing without rinsing





### **Sensitive Glove**

Sensitive glove is intended to body and face cleaning of fragile or bedridden patients with sensitive skin, without rinsing.

Without perfume, it suits sensitive patient's skin, including children older than 3 years old.

Product reference	992 91 86
Primary packaging	x 12 (pack)
Secondary packaging	Carton of 12 packs
Glove	17 x 23,5 cm
Lotion	Dermatologically tested lotion, Aloe Vera-based <b>Perfume free</b>
Expiry date	3 years from manufacture date




### Shampoo Glove

Shampoo glove is intended to hair cleaning of fragile or bedridden patients, without rinsing. It suits any scalp type and any hair type of patient's skin, including children older than 3 years old.

Product reference	992 92 00
Primary packaging	x 12 (pack)
Secondary packaging	Carton of 12 packs
Glove	17 x 23,5 cm
Lotion	Dermatologically tested lotion, Aloe Vera-based
Expiry date	3 years from manufacture date

### **OUR AIM: YOUR SATISFACTION**



90% of respondents affirmed that CLEANIS products **improve** hygiene in hospital services 91% of healthcare facilities declare that CLEANIS products improve healthcare quality 92% of positive feedback concerning the improvement in quality work life thanks to CLEANIS products

Customer satisfaction survey 2018, conducted in the whole French territory, adressed to healthcare professionals.

Collected feedback: 93.

### TRACK YOUR ORDER

### For your first order or request:

Tel : **01 46 15 18 00** 

01 82 66 13 92

Email: cleanis@cleanis.fr

For your usual orders, please contact the person in charge of your account directly or send an email to cleanis@cleanis.fr. We will contact you as soon as possible.

#### **ORDER STREAM**



Receipt of order request

#### **CLEANIS Customer Service**

- For each order, the Sales Administration department will send you an order confirmation
- After signature of this order confirmation or advance payment (following credit line with our insurance company), our warehouse prepares your order
- The person in charge of your account will send you the availability date or the delivery date of your order



Order delivery

### **CLEANIS factory (Vosges, France)**

- Preparation of the order
- Quality control
- Stock management
- Delivery according to agreed conditions

#### **DELIVERY**

We stay at your disposal to propose you the best solutions of delivery.

Following your country and your requested delivery date, our team can ship your goods by express mail/air/sea/truck. A quotation is included in all our order confirmations where the sales incoterm isn't EXW French warehouse.



Invoicing

#### Service Client CLEANIS

Mailing by email and the originals by DHL (if originals are required) of all documents (invoice, packing list, certificate of origin/of analysis/of conformity, Waybill, AWB, legalised documents...).

