



We, STERILMED MEDICAL ELEKTRIK ELEKTRONIK OTOMASYON INSAAT GIDA

SANAYI VE DIS TICARET LIMITED SIRKETI, as the manufacturer of the Sterilmed WD8, declare that in order to ensure the proper functioning and longevity of the equipment, it is preferable to use the products listed below during the warranty period:

- 1. STERILMED MEDICAL ELEKTRIK ELEKTRONIK OTOMASYON INSAAT GIDA SANAYI VE DIS TICARET LIMITED SIRKETI: Care Base & Care NP
- 2. DTH Sağlık Hizmetleri Tic. A.Ş. Kırıkkale Şb: Care Enzyme & Care NP
- 3. DETRO HEALTHCARE: Detro Matic A & Detro Matic NS

The use of the recommended products helps maintain optimal equipment performance and prevents potential malfunctions caused by non-compliant substances. Failure to use these approved washing solutions may result in the following risks:

- 1. Recirculation pump failure
- 2. Damage to heating elements
- 3. Deterioration of peristaltic pumps

Furthermore, STERILMED MEDICAL ELEKTRIK ELEKTRONIK OTOMASYON INSAAT GIDA SANAYI VE DIS TICARET LIMITED SIRKETI, does not guarantee the disinfection quality or the integrity of the instruments if unapproved washing solutions are used. The preset washing programs and pump calibration are specifically designed for the recommended solutions. Any deviation from these may compromise both the effectiveness of the disinfection process and the proper operation of the equipment.

For further information, please refer to the technical documentation or contact our technical support department.

General Manager

Burhan ZORL

ELEKTRIKÆE MED MEDDEAL

ELEKTRIKÆEEKTRONK OTØBASYON NSAAT
GIDA SANKTI VE DIS TICABET (MITTED STREET
BESJON Mallson SANKTI NE DIS TICABET (MITTED STREET
BESJON Mallson SANKTI NE DIS TICABET MITTED STREET
LE 1/90 312 37837 10 0 7 65 49 301 304 45 77001
KORMY DIS TABETAGET MITTED NO. 1781 104 4577 MOIS
TICABET SET DE 28029 9 100 Gestatirned John TI

STERİLMED MEDİCAL ELK. ELKT. OTOMASYON İNŞ. GIDA SAN. VE DIŞ. TİC. LTD. ŞTİ.

Başkent Organize Sanayi Bölgesi 18. Cad. No: 43 06909 Malıköy, Sincan ANKARA / TÜRKİYE

Tel: +90 312 375 81 00 Faks: +90 312 375 92 92