

Instruction for Use



[Product] Patient Examination Gloves (Nitrile Examination Gloves)

—Sizes: XS, S, M, L, XL

[REF] LMNEGBE01, LMNEGVT01, LMNEGBK01, LMNEGBP01, LMNEGCB01, LMNEGWE01, LMNEGPK01

[Intended Use] Patient Examination Gloves (Nitrile Examination Gloves) is a non - sterile disposable device intended for medical purpose that is worn on the examiner's hand to protect patient by preventing contamination between patient and examiner.

[Storage] It shall be stored in a dry and cool place with good ventilation.

[Shelf Life] 3 years












[Cautions]


Do not use if the product is expired. Do not use if the package is damaged. Do not reuse.








[Disposal]

Disposal restrictions depend upon the contamination incurred during use and are subject to national or local legislation.

[Explanation of product labels and packaging marks]:

1		Medical device
2		Manufacturer
3		Consult instructions for use
4		Do not re-use
5		Batch code
6		Date of manufacture
7		Use-by date
8		Caution
9		Non-sterile
10		Keep away from sunlight
11		Keep dry

	Patient Examination Gloves (Nitrile Examination Gloves)	
	Instruction for use	
	File No.: L-CE0007-08-02	Version: A/1

12		Do not use if package is damaged and consult instructions for use	
13		Catalogue Number	
14		Authorized representative in the European Community	
15		CE mark	
16		Latex free	
17		Fragile, handle with care	



Lyncmed Medical Technology (Beijing) Co., Ltd.
Room 1601, Building No. 2, Zhubang 2000 Business Building, Balizhuang Xili 99,
Chaoyang District, 100022 Beijing, P.R. China



MedNet EC-REP GmbH
Borkstrasse 10, 48163 Münster, Germany