

ESR

ESR Tube contain 3.8% sodium citrate for ESR test.ESR measurements refer to the Westergren method.

The tube are available with an open or closed sedimentation system.

Usage: ESR Blood Collection Tube are intended to be exclusively used for venous blood collection for determination of the Erythrocyte Sedimentation Rate.

Specimen: whole blood.

REF	Tube Size (mm)	Draw Volume (mL)	Tube Materia	Pack/Carton
C110016	13x75	1.6mL	PET	100PCS/1200PCS
C310016			GLASS	100PCS/1200PCS
C110020	13x75	2mL	PET	100PCS/1200PCS
C310020			GLASS	100PCS/1200PCS
C110024	13x75	2.4mL	PET	100PCS/1200PCS
C310024			GLASS	100PCS/1200PCS
C110030	13x75	3mL	PET	100PCS/1200PCS
C310030			GLASS	100PCS/1200PCS
C110040	13x75	4mL	PET	100PCS/1200PCS
C310040			GLASS	100PCS/1200PCS
C111012	8x120	1.28mL	PET	100PCS/1200PCS
C311012			GLASS	100PCS/1200PCS
C111016	8x120	1.6mL	PET	100PCS/1200PCS
C311016			GLASS	100PCS/1200PCS



No Additive

Plain Without Additives Blood Collection Tube are intended to be used in laboratory mainly for biochemistry, immunological & serological tests.

Usage: researches of biological fluids in clinical chemistry, serology, immunology.

Specimen: serum or other biological fluids.

Clotting time: 60 minutes.

Conditions of centrifugation: 1300 g for 10 minutes.Before centrifuging the tube with serum it is necessary to wait until blood clotting.

REF	Tube Size (mm)	Draw Volume (mL)	Tube Materia	Pack/Carton
C100010	13x75	1mL	PET	100PCS/1200PCS
C300010			GLASS	100PCS/1200PCS
C100020	13x75	2mL	PET	100PCS/1200PCS
C300020			GLASS	100PCS/1200PCS
C100030	13x75	3mL	PET	100PCS/1200PCS
C300030			GLASS	100PCS/1200PCS
C100040	13x75	4mL	PET	100PCS/1200PCS
C300040			GLASS	100PCS/1200PCS
C100050	13x100	5mL	PET	100PCS/1200PCS
C300050			GLASS	100PCS/1200PCS
C100060	13x100	6mL	PET	100PCS/1200PCS
C300060			GLASS	100PCS/1200PCS
C100070	13x100	7mL	PET	100PCS/1200PCS
C300070			GLASS	100PCS/1200PCS
C100080	16x100	8mL	PET	100PCS/1200PCS
C300080			GLASS	100PCS/1200PCS
C100090	16x100	9mL	PET	100PCS/1200PCS
C300090			GLASS	100PCS/1200PCS
C100100	16x100	10mL	PET	100PCS/1200PCS
C300100			GLASS	100PCS/1200PCS

