


ARTHROPLASTY SYSTEMS




normmed
A Piece of You

HECTOR


NEW
TOTAL HIP
SYSTEM



Titanium
Plasma+HA
Coating
Acetabular Cup




Titanium
Plasma+HA
Coating
Stem System



CoCr Material
Femoral
Heads

PREVENT
ROTATION
**BEST
DESIGN**


Polyaxial
Acetabular
Screws



High Cross Linked
Polyethylene (XLPE)
Material




Total
Hip System



Cementless Stem

Version Size

Standart 01,02,03,04,05,06,07,08,09,10
Lateral 01,02,03,04,05,06,07,08,09,10



Cemented Stem

Version Size

Standart 03,04,05,06,07,08,09,10
Lateral 03,04,05,06,07,08,09,10



CoCr Femoral Head

Diameter (mm) Size

28 -3,+0,+4,+8,+12
32 -3,+0,+4,+8,+12
36 -3,+0,+4,+8,+12



Bipolar Head

Diameter (mm)

44,46,48,50,52,
54,56,58,60



Acetabular Screw

Diameter (mm) Length(mm)

6.5 15,20,25,
30,35,40



Acetabular Cup

Diameter (mm)

46,48,50,52,54,
56,58,60,62,64



Acetabular Liner 0°

Inner Ø(mm) Outer Ø(mm)

28,32,36 46,48,50,52,54,
56,58,60,62,64



Acetabular Liner 20°

Inner Ø(mm) Outer Ø(mm)



28,32,36 46,48,50,52,54,
56,58,60,62,64



Centraliser

Size

01,05,09



NORMMED MEDICAL INDUSTRY AND TRADE INC.

İvedik OSB Mah. Havalandırmacılar Cad. No:193

Yenimahalle / ANKARA / TURKEY

info@normmed.com.tr

www.normmed.com.tr

**Total
Hip System**

NEW

TOTAL KNEE
SYSTEM

GORDION

CoCr Material
Femoral
Components

Glossy Surface

Sandblasting
Titanium Material
Tibial Components

SINGLE
RADIUS
HIGH
MOBILITY

Ultra High
Molecular Weight
Polyethylene
Tibial Inserts

Total
Knee System



Knee Prosthesis PS Cemented Femoral Component

Size	Side
01,02,03,04,05,06	R
01,02,03,04,05,06	L



Knee Prosthesis CR Cemented Femoral Component

Size	Side
01,02,03,04,05,06	R
01,02,03,04,05,06	L




Fix Cemented Tibial Component

Size
01,02,03,04,05,06



Standart Patellar Component

Diameter(mm)
29,32,35,38



Fix CR Tibial Insert

Size	Height (mm)
01	9,11,13,15
02	9,11,13,15
03	9,11,13,15
04	9,11,13,15
05	9,11,13,15
06	9,11,13,15



Fix PS Tibial Insert

Size	Height (mm)
01	9,11,13,15
02	9,11,13,15
03	9,11,13,15
04	9,11,13,15
05	9,11,13,15
06	9,11,13,15

