Disposable skin staplers



■ W W W . J S R I P . C O M

Disposable endoscopic specimen retrieval bag



■ W W W . J S R I P . C O M

PRODUCT FEATURE

- Easy operation, flexible and convenient. Nail out and suture smoothly
 and accurately, safe, high-efficiency and quickly;
- The wound of skin suturing is small. Avoid re-damage of epidermis tissue, and more helpful for would healing;
- 3. Beautiful and smooth suturing the wound is without "centipede feet", Healing scar is slight and neat;
- 4. Shorten suturing time, reduce pain of patients, no painful feeling when remove the staples;
- 5. Special staple remover, reduce injury of epidermal tissue and comfort the pain.

SCOPE OF APPLICATION

Be applicable for closing skin wounds.

SPECIFICATIONS AND TECHNICAL PARAMETERS

Specifications	Formed size of staple (mm)	No. of staples
RYPF-15W/R	7.0*4.5/6.0*4.0	15
RYPF-25W/R	7.0*4.5/6.0*4.0	25
RYPF-35W/R	7.0*4.5/6.0*4.0	35





PRODUCT FEATURE

- 1. New ergonomic design makes a more comfortable hand-feel, controllability and stability;
- 2. To prevent tumor seeding and inflammation by tissue fluid penetration;
- 3. The bag can be expanded automatically when push the handle, it's easy to collect specimen;
- 4. Endo bag is designed in accordance to the minimally invasive theory, it can go through small trocar (10mm);
- 5. The bag is transparent and soft, with makes a good visually.

SCOPE OF APPLICATION

Endo bag is applicable for minimally invasive surgeries, to remove specimens from body.

SPECIFICATIONS AND TECHNICAL PARAMETERS

Model	Volume of Bag (ml)
RYQW40	≥40
RYQW60	≥75
RYQW80	≥200
RYQW100	≥450
RYQW130	≥700
RYQF100	≥800

