

## Amplatz Dilators single

**Hydrophilic coated**

Used for progressive dilation of the nephrostomy tract prior to percutaneous stone removal. Radiopaque dilators.

### Single

Art. No.	Ø (Fr/CH)	Length (cm)
<b>RD 18</b>	18	30
<b>RD 20</b>	20	30
<b>RD 22</b>	22	30
<b>RD 24</b>	24	30
<b>RD 26</b>	26	30
<b>RD 28</b>	28	30
<b>RD 30</b>	30	30

## Amplatz Sheaths single

**Hydrophilic coated**

Used during renal dilatation to provide atraumatic working tracks after removal of the dilator.




Art. No.	Dilators Ø (Fr/CH)	Length (cm)
<b>AM 10</b>	10	30
<b>AM 24</b>	24	17
<b>AM 26</b>	26	17
<b>AM 28</b>	28	17
<b>AM 30</b>	30	17

## Initial Puncture Needles

Used for percutaneous Nephrostomy and all other percutaneous applications where a 0.035" guidewire has to be used.

### 2-Part Puncture Needle, Trocar tip



Art. No.	Needle Ø (GA)	Length (cm)
<b>IPN-18-20-2T</b>	18	20

### 3-Part Puncture Needle, Bevel tip



<b>IPN-18-20-3B</b>	18	20
---------------------	----	----

### Chiba Needles



<b>CH-22-20</b>	22	20
<b>CH-22-28</b>	22	28