

Gemini - Semi Automatic Microtome
for histological samples



Technical features

S5	Sectioning thickness setting range	S6	Trimming thickness setting range
	Range of the slice thickness: 0,25 - 100 µm		Range of the slice thickness: 1 - 600 µm
	0,25 to 2,5 µm in 0,25 µm increments		1 to 10 µm in 1 µm increments
	2,5 to 5 µm in 0,5 µm increments		10 to 20 µm in 2 µm increments
	5 to 10 µm in 1 µm increments		20 to 50 µm in 5 µm increments
	10 to 30 µm in 2 µm increments		50 to 150 µm in 10 µm increments
	30 to 60 µm in 5 µm increments		150 to 600 µm in 50 µm increments
S8	Movable range for the blade holder basal	S14	Maximum specimen size
	Forward/Backward: ±65 mm		Standard clamp (HxWxD): 55x50x30 mm
	Left / Right: ±40 mm		Super mega cassette clamp (HxWxD): 68x48x15 mm
S4	Vertical stroke	S14	Specimen speed
	72 mm ±1 mm		Variable speed: 0,5 – 420 mm/s ± 10%
S12	Specimen feed	S14	Size and weight
	30 mm ±1 mm		Width: 475 mm
S12	Retraction	S14	Depth (including tray): 580 mm
	Manual mode: 5-100 µm in increments of 5 µm		Depth (excluding tray): 510 mm
S12	Automatic mode: according to speed, can be deactivated	S14	Height: 340 mm
	Weight: 40 kg		Weight: 40 kg
S12	Spieccimen orientation	S14	Certifications
	Horizontal: 8 , vertical: 8 , rotation 360°		CE-IVD

For up-to-date featurress, please refer to the instrument data sheet.
Technical specifications are subject to change without notice.

Distribuitied by:



Histo-Line Laboratories
Histo-Line Laboratories s.r.l.
Viale G. di Vittorio, 30 | 20048 Pantigliate (MI) ITALY
Phone +39 0255230061 | Fax +39 0255213764
histoline@histoline.com | www.histoline.com



Gemini - Semi Automatic Microtome
for histological samples

Applications of routine histopathology

MRS3500



English

Gemini - Semi Automatic Microtome

for histological samples

The system advantages



Convenient storage area
Recessed staging area for convenient placement of blades, brushes and other ancillary tools.

S1

For optimal cutting comfort.

S15

Integrated waste trays to save you time and energy



Flexible
Long stroke length allows high quality sections of Macro / SuperMega cassettes



ERGONOMICS

Multi-function mode
Complete function of automatic sectioning, slice counting, microns counting, fast trimming. **Two different trimming memory positions** can be set by user. This allow a fast sample position.

S16



The precision to cut
Microtome for cut quality sections with consistent thickness.
Produce high quality paraffin sections for H&E, IHC and molecular techniques.



Stability
Silent retraction reduces noise during prolonged use.
Maximum rigidity and stability are provided by the optimized microtome base plate.

S1



Lateral knife holder adjustment
Precise lateral knife holder adjustment, with three predefined stop positions.
Blade removal system and blade holder with cutting protection for operator's safety

S1



Control panel
All functions are controlled, space-saving control panel with display and buttons allow you to make changes in use.

The perfect combination of efficiency and user comfort

Simple and smart microtome for routine histopathology applications and research field



S10



Adapted to your needs
Choose between conventional, full hand-wheel rotation, manual sectioning or "rocking mode", where the hand-wheel is turned back and forth over a short distance.



Safety during sectioning
The hand wheel can be locked at any position, so as to ensure the safety of the operation.
Safe quick braking system; precise electronic mechanic design

S14

Precision specimen orientation
Accurately align the block face to the blade edge, using the specimen orientation system with zero-position centering.
The blade ejector and magnetic brush enable the user to safely remove the blade from the blade holder without touching it.

S16

Blade holder with adjustable cutting angle to the use of any type of disposable blades.



Ergonomics and safety
A hand-wheel locking device for one-thumb operation locks the specimen head in the uppermost position, which is ideal for changing the specimen and/or blade.

Technical support and after sale

Histo-Line offers a complete technical support for you.
We find together the best solution for you and your instrument in certain and quick times.

Technical Support Products specialist@histoline.com
Technical Support Equipments techservice@histoline.com



Wide range of accessories and consumables on request

Histo-Line complete range of solutions for Microtomy



MRS3750
Aries
Automatic Microtome



MR3000
Libra
Manual Microtome