

Angiographic Syringes

Medical Imaging / www.scw-medicath.com



| Product No. | | Volume(ml) | Description | Models | |
|-------------|-------------|------------|--------------------------------------|---------------|---------------|
| Single | Dual | | Used for Libel-Flarsheim | Single | Dual |
| 42.16.20000 | 42.16.20001 | 200 | Used for CT: CT 9000 & CT 9000 ADV | CT-200-B1 | CT-200/200-B1 |
| 42.16.20002 | | 150 | Used for DSA: Angiomat 6000 | DSA-150-B1 | |
| 42.16.20003 | | 150 | Used for CT & DSA: Angiomat Illumena | CT/DSA-150-B1 | |
| 42.16.20005 | 42.16.20006 | 60 | Used for MRI: Optistar | MRI-60-B1 | MRI-60/60-B1 |

| Product No. | | Volume(ml) | Description | Models | |
|-------------|-------------|------------|----------------------|-----------|---------------|
| Single | Dual | | Used for EZEM | Single | Dual |
| 42.16.50000 | 42.16.50001 | 200 | Used for CT: Empower | CT-200-EM | CT-200/200-EM |

Medical Imaging / www.scw-medicath.com



| Product No. | | Volume(ml) | Description | Models | |
|-------------|-------------|------------|---------------------------------|------------|---------------|
| Single | Dual | | Used for Nemoto Kyorindo | Single | Dual |
| 42.16.30000 | 42.16.30001 | 100 | Used for CT: A-25 | CT-100-NE | CT-100/100-NE |
| 42.16.30002 | 42.16.30003 | 200 | Used for CT: A-25 A-60 | CT-200-NE | CT-200/200-NE |
| 42.16.30005 | | 120 | Used for DSA: Nemoto 120S | DSA-120-NE | |
| 42.16.30006 | 42.16.30007 | 60 | Nemoto Dual Shot. Sonic Shot 50 | MRI-60-NE | MRI-60/60-NE |
| | 42.16.30008 | 100/200 | Used for CT:A-25 | | CT-100/200-NE |
| | 42.16.30012 | 60/200 | Used for CT: A-25,A-60 | | CT-60/200-NE |

| Product No. | | Volume(ml) | Description | Models | |
|-------------|-------------|------------|------------------------|------------|----------------|
| Single | Dual | | Used for Metron | Single | Dual |
| 42.16.40000 | 42.16.40001 | 200 | Used for CT: ACCUTRON | CT-200-MED | CT-200/200-MED |
| 42.16.40002 | 42.16.40003 | 65 | Used for MRI: ACCUTRON | MRI-65-MED | MRI-65/65-MED |
| | 42.16.40004 | 65/200 | Used for MRI: ACCUTRON | | MRI-65/200-MED |