



ECLIPSE[®] BCP

Composite interference screw



Distributed by



AMPLITUDE[®]

Features


- Canulated composite interference screw
- Progressive insert
- Range :
 - Diameters 7/8/9 mm Lenght 20/25 et 30mm
 - Diameters 10mm Lenght 25 et 30mm
 - Diameters 11/12 mm lenght 35mm



Material


PLADLLA 75% / BCP 25%

- PLDLLA, polylactic acide L & D
- BCP, 20% hydroxyapatite
80% phosphate tricalcique




Study

- Comparative study** by CT-Scan at 2 years follow up show PLDLLA interference screw with BCP can reduce the extent of resorption while improving ossification and reducing tunnel widening. can reduce the extent of resorption while improving ossification, reducing tunnel widening,



Screwing instrumentation : the choice

- Ratcheting trinkle driver handle : Amplitude
 - Shank screwdriver
 - Shank starter



- Screw Guidewire 1.1 mm diameter
 - Rigid : Amplitude
 - Nitinol : Biomatlante



- Monobloc: Biomatlante
 - screwdriver
 - starter

References

ECLIPSE® BCP Composite interference screw

Manufactured by Biomatlante  0123

Composite Interference Screw ECLIPSE® BCP Ø 7 Lg 20 mm	1-0402771
Composite Interference Screw ECLIPSE® BCP Ø 7 Lg 25 mm	1-0402772
Composite Interference Screw ECLIPSE® BCP Ø 7 Lg 30 mm	1-0402773
Composite Interference Screw ECLIPSE® BCP Ø 8 Lg 20 mm	1-0402781
Composite Interference Screw ECLIPSE® BCP Ø 8 Lg 25 mm	1-0402781
Composite Interference Screw ECLIPSE® BCP Ø 8 Lg 30 mm	1-0402783
Composite Interference Screw ECLIPSE® BCP Ø 9 Lg 20 mm	1-0402791
Composite Interference Screw ECLIPSE® BCP Ø 9 Lg 25 mm	1-0402792
Composite Interference Screw ECLIPSE® BCP Ø 9 Lg 30 mm	1-0402793
Composite Interference Screw ECLIPSE® BCP Ø 10 Lg 25 mm	1-0402702
Composite Interference Screw ECLIPSE® BCP Ø 10 Lg 30 mm	1-0402703
Composite Interference Screw ECLIPSE® BCP Ø 11 Lg 35 mm	1-0401004
Composite Interference Screw ECLIPSE® BCP Ø 12 Lg 35 mm	1-0401024

Reference : FT.GB.013/3.0

Instrumentation Set

Manufactured by AMPLITUDE  2797

Modular Instrumentation Set In-Out	2-0299940
Modular Instrumentation Set Out-In	2-0299957

Creation date : March 2019 – Composite interference screw BCP Legal manufacturer : Biomatlante • 5 rue Édouard Belin • Z.A. Les IV Nations • 44360 • Vigneux-de-Bretagne • France – Destination : medical device for Knee ligament surgery
– instrumentation set legal manufacturers : Amplitude , 11 cours Jacques Offenbach, 26000 Valence, France and Biomatlante • 5 rue Édouard Belin • Z.A. Les IV Nations • 44360 • Vigneux-de-Bretagne • France

Customer Service-France :

Porte du Grand Lyon,
01700 Neyron, France
Tel. : +33 (0)4 37 85 19 19
Fax : +33 (0)4 37 85 19 18

E-mail : amplitude@amplitude-ortho.com

Customer Service-Export :

11, cours Jacques Offenbach. Zone Mozart 2,
26000 Valence, France
Tel. : +33 (0)4 75 41 87 41
Fax : +33 (0)4 75 41 87 42

Internet : www.amplitude-ortho.com