

ABDOMINAL WALL REINFORCEMENT FIXATION IN LAPAROTOMY AND LAPAROSCOPY



Non-absorbable
fixation device

Product specifications

Surgical approach	Laparoscopy and laparotomy
Single use	Yes
Material	Helical titanium tack
Diameter	5 mm
Sterilization	Ethylene oxide
Number of tacks	ProTack™: 30 Stat Tack™: 15

References and descriptions

ProTack™ 5 mm

Fixation device with non-absorbable tacks



Reference	Description	Qty/Box
174006	ProTack™ Fixation device containing 30 titanium tacks	6

Stat Tack™

Fixation device



Reference	Description	Qty/Box
OMS-TTSS	Stat Tack™ Fixation device containing 15 titanium tacks	5

PLEASE REFER TO THE PRODUCT INSTRUCTIONS
FOR USE BEFORE USING THE PRODUCT.

FOR HEALTHCARE PROFESSIONALS ONLY.

© Medtronic, France. September 2017. CARL: 1303.
All rights reserved. Photo credit: Medtronic.
Printed by RR Donnelley. www.rrdgs.com.
Produced outside the system

www.medtronic.fr

Medtronic France S.A.S
27 Quai Alphonse Le Gallo | CS30001
92513 Boulogne-Billancourt Cedex
Tel. +33(0)1 55 38 17 00
Fax: +33(0)1 55 38 18 00

RCS Nanterre 722 008 232

Medtronic